Elite

POSITIONING BRODA

Tilt & Positioning Wheelchairs



Shown with 6" Multi-Directional Locking Rear Casters, and long armrests options

Patient, Facility, & Caregiver Benefits

- Facilitates positioning for comfort and respiratory function
- Tilt and recline provides ease of use for frequent re-positioning
- Proper seating alignment decreases postural deviations



Tilt-in-Space functionality allows easy repositioning of users



Removable armrests make transfers safer for caregivers and patients



Posterior Tilt effectively opens the diaphragm improving digestion, oxygenation, and organ function



Comfort Tension Seating® increases support and sitting tolerance





Industry-leading Positioning, Versatility, and Functionality

UPHOLSTERY OPTIONS









FRAME OPTIONS






SPECIFICATIONS

Degrees of Tilt	Up to 40° Posterior
Weight Capacity	Up to 350 lbs. (159 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45 cm, 50cm, 55cm, 60cm, 65cm)
Seat Depths	15.5", 17",18.5", 20" (37.5cm, 43cm, 46cm, 50cm)
Seat Heights	19", 20.5",22", 23.5" (48cm, 51cm, 55cm, 59cm)
Wheels	5" Total Lock Front Casters (13cm) 14" or 20" Rear Mag Wheels (36cm or 50cm) 14" Mid Mag Wheels (36cm) 6" Multi-directional Locking Casters (15cm)
Lifetime Frame Warranty	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Lower side pads, upper side pads, seat and back pad with cover neckrest, upper laterals and, sole/calf pads

OPTIONAL ACCESSORIES











Clear Tray

Padded Thigh Strap

Additional accessories available: foot supports, trays, and positioning straps