

# Encore

## Pedal Wheelchair

PEDAL  
**BRODA**®



**TIME**  
**B**  
New Product  
Award 2019

2019  
**TODAY'S CAREGIVER**  
**FRIENDLY**  
**AWARD**  
CAREGIVER.COM


**mobility**  
product award 2019



*New frame design, including adjustable wings and tall back option*



*Front pivot seat tilt maintains a proper position for self-propulsion*




*Optional dynamic rocking feature calms and creates contentment*

*Shown with short back, Matrix® Elan headrest and lateral support options*

## Patient, Facility, & Caregiver Benefits


- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning



**Comfort Tension Seating®**  
increases support as well as sitting tolerance

## *The Encore of Comfort, Functionality, and Style*

### UPHOLSTERY OPTIONS



### FRAME OPTIONS



### SPECIFICATIONS

<b>Tilt</b>	Up to 33° (with Dynamic Rocking)
<b>Weight Capacity</b>	Up to 350 lbs. (114 kg)
<b>Seat Widths</b>	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm) 16", 18", 20", 22" with Dynamic Rocking (40cm, 45cm, 50cm, 55cm)
<b>Seat Depths</b>	17", 18.5", 20" (43cm, 46cm, 50cm) 15.5", 17", 18.5", 20" with Dynamic Rocking (39cm, 43cm, 46cm, 50cm)
<b>Back Lengths</b>	22", 32", 33.5", 35" (55 cm, 81cm, 85cm, 89cm)
<b>Seat Surface Heights</b>	Adjustable 12.5" - 17.5" (31cm-44cm) Adjustable 14" - 19" with Dynamic Rocking (35-48cm)
<b>Wheels</b>	5" Front Casters: Total Lock (12.5cm) 20" or 22" Front Mag (50cm or 55cm) 14" Rear Mag (36cm) 5" Rear Casters (12.5cm)
<b>Lifetime Warranty Frame</b>	10 years on 16 Gauge tubular steel frame
<b>Pressure-Reducing Seat &amp; Backrest</b>	Integrated with Comfort Tension Seating®
<b>Standard Padding</b>	Tall back: Neckrest, shoulder bolsters, lower side pads, seat and back pad with cover Short back: lower side pads, seat and back pad with cover

Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps