

Encore

TRANSPORT
BRODA®

Pedal Wheelchair WC-19 Transport Option



**Crash-tested and WC-19
compliant for vehicle transport**



**Front pivot seat tilt keeps users
engaged with surroundings**



**Brackets enable safe patient
transport between locations**






**Tilts for in-transit contentment
and safety**

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning

Comfort & Safety on the Go




UPHOLSTERY OPTIONS



SPECIFICATIONS

Tilt	Up to 25°
Weight Capacity	Up to 300 lbs. (136 kg)
Seat Widths	16", 18", 20", 22", 24", (40cm, 45cm, 50cm, 55cm, 61cm)
Seat Depths	17", 18.5", 20" (43cm, 46cm, 50cm)
Back Heights	32", 33.5", 35" (81cm, 85cm, 89cm)
Seat Surface Heights	Adjustable 12.5" - 17.5" (31cm-44cm)
Wheels	20" or 22" Front Mag (50cm or 55cm) 5" Rear Casters (12.5cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Tall back: Neckrest, shoulder bolsters, lower side pads, seat and back pad with cover

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps