

Centric

Positioning Wheelchair

POSITIONING
BRODA[®]



Posterior tilt aids in pressure redistribution




Tilt with semi-recline accommodates hip flexion limitations, allowing patients to be out of bed



Optional contracture sling accommodates various levels of knee contractures

Patient, Facility, & Caregiver Benefits

- Tilt and recline provides ease of use for frequent repositioning
- Design provides safety and stability helping to reduce risk of injury
- Proper seating alignment decreases postural deviations reducing the use of restraints



Comfort Tension Seating[®] improves postural support and sitting tolerance

UPHOLSTERY OPTIONS



FRAME OPTIONS



OPTIONAL ACCESSORIES



Padded Thigh Belt



Tray



O² Holder



Half Tray

Essential Positioning And Peace of Mind

SPECIFICATIONS

Degrees of Tilt	Up to 25°
Degrees of Recline	Up to 45°
Weight Capacity	250 lbs. (114 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 63cm, 66cm)
Seat Depths	15.5", 17", 18.5", 20" (39cm, 43cm, 46cm, 50cm)
Seat Surface Heights	Adjustable 19.5", 20.5", 21.5" (49cm, 51cm, 54cm)
Wheels	5" (12.5cm) Total Lock Front or Rear Caster Wheels 20", 22", 24" (50cm, 55cm, 60cm) Front/Rear Mount Mag Wheels 5" (12.5cm) Directional Rear Caster Wheels
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating [®]
Standard Padding	Head support cushion, lateral shoulder supports, lower side cushions, upper side cushions, calf cushion, sole cushion, seat and back cushions with cover

Additional accessories available: foot supports, specialty padding packages, and positioning straps