

Matrx® Backs

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height: <input type="radio"/> 16" <input type="radio"/> 18" (excludes PB & PB Deep) <input type="radio"/> 20"	
Back Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Back Width: <input type="radio"/> 22"
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$866	<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$1,151
<input type="radio"/> Elite - Recommended for more active users \$766	<input type="radio"/> Elite - Recommended for more active users \$1,051
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs \$918	<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs \$1,095
<input type="radio"/> PB - Recommended for less active users seeking comfort \$734	
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support \$892	
Lateral Support Accessories: <input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L \$116	Privacy Flap Accessories: <input type="radio"/> 16" Width Small Privacy Flap \$53 <input type="radio"/> 18" Width Medium Privacy Flap \$53 <input type="radio"/> 20" Width Large Privacy Flap \$53 <input type="radio"/> 22" Width Extra Large Privacy Flap \$53
<input type="radio"/> Elite Back Swing away Lateral Hardware \$243	
<input type="radio"/> Swing away Lateral Pad (each) 4.25"Hx6.5"L \$122	

Matrx® Headrest

Head Support Mounting: <input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support \$285
Head Support Options: <input type="radio"/> Elan Head Support <input type="radio"/> 10" <input type="radio"/> 14" \$310
<input type="radio"/> Elan 4 Point Head Support <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10" \$448
<input type="radio"/> Elan Headband with 4 Point Cushion <input type="radio"/> Small <input type="radio"/> Large \$610

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorials Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorials services needed to deliver are requested when an order is submitted to Broda.

Lift Gate Service: Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

Residential Delivery Service: Delivery to a residential location

Notification Before Delivery: Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

Inside Delivery Service: Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

Limited Access Service: Delivery locations that are not open to the walk-in public, during normal business hours, Mon-Fri, 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.

587V3: SEAT HEIGHT & WHEEL COMBINATIONS								
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB
14"	20	5	SMALL	1	1	TOP	3	2
	22	5	SMALL	1	2	TOP	2	1
15.5"	20	6	MEDIUM	2	1	TOP	2	1
	22	6	MEDIUM	1	1	TOP	2	1
	24	6	SMALL	1	2	TOP	1	1
	24	5	MEDIUM	1	2	TOP	1	1
16"	24	6	MEDIUM	1	1	TOP	1	1
	20	5	SMALL	3	3	TOP	1	1
	20	5	MEDIUM	3	3	TOP	1	1
	20	6	MEDIUM	3	2	TOP	1	1
	22	5	SMALL	2	3	TOP	1	1
	22	6	SMALL	2	2	TOP	1	1
16.5"	22	5	MEDIUM	2	3	TOP	1	1
	22	6	MEDIUM	2	2	TOP	1	1
	24	6	SMALL	2	3	BOTTOM	3	1
	24	5	MEDIUM	2	4	BOTTOM	3	1
17"	24	6	MEDIUM	2	3	BOTTOM	3	1
	20	5	MEDIUM	4	5	BOTTOM	3	1
	20	6	MEDIUM	4	4	BOTTOM	3	1
	22	5	MEDIUM	3	5	BOTTOM	3	1
17.5"	22	6	MEDIUM	3	4	BOTTOM	3	1
	20	8	SMALL	4	3	BOTTOM	3	1
	22	8	SMALL	3	3	BOTTOM	3	1
	24	8	SMALL	3	3	BOTTOM	2	1
18"	24	6	MEDIUM	3	5	BOTTOM	2	1
	24	8	MEDIUM	3	3	BOTTOM	2	1
	20	8	MEDIUM	5	4	BOTTOM	2	1
	22	8	MEDIUM	4	4	BOTTOM	2	1
19"	24	8	MEDIUM	4	5	BOTTOM	1	1

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.
Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG-404 587V3 Rev. 1