

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
<ul style="list-style-type: none"> - Tubular Steel, Powder Coated Frame - 45 Degrees, Infinitely Adjustable Seat Tilt - Fixed Height Arm Supports - Adjustable Seat Depth 16" - 20" - Adjustable Seat Surface Height 14" - 19" - Choice of Wheel, Caster and Fork Sizes - Push Lock Parking Brakes 	<ul style="list-style-type: none"> - Drop In, Swing Away Lower Leg Supports - Composite Foot Support Plates - Cushioned Side Panels - Adjustable Backrest Angle - 250 lbs. Maximum Occupant Weight - Seat, back and head support cushions not included

Comfort Rehab Wheelchair (587) - MSRP \$2,724.75


Frame Color Options

Frame Color: Black Blue Bronze Nickel Orange Red Sand Yellow

Seat and Back Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	N/C \$204.75 \$204.75 N/C	Seat Depth:
			<input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20"
			N/C N/C N/C N/C

Seat Surface Height: <input type="radio"/> 14" <input type="radio"/> 14.5" <input type="radio"/> 15" <input type="radio"/> 15.5" <input type="radio"/> 16" <input type="radio"/> 16.5" <input type="radio"/> 17" <input type="radio"/> 17.5" <input type="radio"/> 18" <input type="radio"/> 19"	Back Cane Height: <input type="radio"/> Ultra-Low Double Bend (16") Fixed <input type="radio"/> Ultra Low (16") Fold Down <input type="radio"/> Low (18") Fold Down <input type="radio"/> Medium (21") Fold Down <input type="radio"/> Tall (24") Fold Down
N/C N/C N/C N/C N/C N/C N/C N/C N/C N/C	Standard \$231 \$231 \$231 \$231 \$231

Refer to Seat Height & Wheel Combinations Matrix for assistance.



Front Mount Options	Wheel Options	Rear Mount Options
<input type="radio"/> 5" Fixed Caster Wheels - Pair Standard <input type="radio"/> 6" Fixed Caster Wheels - Pair N/C <input type="radio"/> 8" Fixed Caster Wheels - Pair N/C	<input type="radio"/> 20" Mag Wheels - Pair Standard <input type="radio"/> 22" Mag Wheels - Pair N/C <input type="radio"/> 24" Mag Wheels - Pair N/C	

Caster Forks: Required <input type="radio"/> Small - Recommended for Small Mag Wheel (20" or 22") and Lower Seat Surface Heights N/C <input type="radio"/> Medium - Recommended for Large Mag Wheel (24") Higher Seat Surface Heights N/C	Mag Wheel Position: Factory Set (Adjustable in the future) <input type="radio"/> A (Most Forward) - Smaller Turning Radius N/C <input type="radio"/> B N/C <input type="radio"/> C N/C <input type="radio"/> D* (Most Rearward) - Larger Turning Radius N/C <small>*Do not recommend with 20" Mag Wheels</small>
--	--

Parking Brake Extension Handles - Pair N/C Anti-Tip Devices - Required Safety Features - Factory Installed N/C

Height Adjustable, Removable Arm Supports N/C

Lower Leg Support Assembly and Lower Leg Support Options

Lower Leg Support Assembly Style: <i>must choose one</i>	Lower Leg Support Accessories:
<input type="radio"/> Drop In Straight - Pair N/C <input type="radio"/> Drop In Flared - Pair N/C <input type="radio"/> Pin and Saddle Straight - Pair N/C <input type="radio"/> Pin and Saddle Flared - Pair N/C <input type="radio"/> Drop In IKON Style - Pair N/C <input type="radio"/> F-SW75ADJ-587V3 Swing Away 75° Lower Leg Supports - Pair N/C <input type="radio"/> F-SW60ADJ-587V3 Swing Away 60° Lower Leg Supports - Pair N/C <input type="radio"/> F-SW-IKN N/C <input type="radio"/> F-SW-IKN-ADJ Swing Away Lower Leg Supports Black w/ Angle Adjustable Flip-up Foot Supports N/C <input type="radio"/> 2" Extended Footrest Tubes - <i>Only compatible on the F-SW-IKN and F-SW-IKN-ADJ</i> \$63.00	<input type="radio"/> Posterior Foot Supports - Pair \$21.00 <input type="radio"/> Lower Leg Support Panel (for calf support) \$63.00 <input type="radio"/> Lower Leg Support Foot Board \$115.50 <input type="radio"/> F-SW-ELV Swing Away, Elevating Lower Leg Supports Black - Pair \$341.25 <input type="radio"/> F-SW-ELV-IKN Swing Away, Elevating Lower Leg Supports Silver - Pair \$425.25 <input type="radio"/> F-SW-ANGLE Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair \$283.51 <input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support w/ Angle Adjustable Footplates Black - Pair

Common Accessories

Pelvic Strap <input type="radio"/> Pelvic Positioning Strap Grommet Mount <input type="radio"/> Cir Buckle <input type="radio"/> Side Buckle \$54.60	Trays: <input type="radio"/> ABS Tray with Strap (16", 18" 20" 22" - Matches Corresponding Seat Width) <input type="radio"/> Sewn Strap <input type="radio"/> Grey <input type="radio"/> Black \$319.20 <input type="radio"/> Padded ABS Tray Cover (ABS Tray only) \$162.75 <input type="radio"/> Clear Tray with Strap (16", 18", 20" - Matches Corresponding Seat Width) \$329.70
Other: <input type="radio"/> IV Pole and Assembly \$102.90 <input type="radio"/> Oxygen Holder and Assembly \$173.25 <input type="radio"/> Accessory Bag (Black) Each \$43.05 <input type="radio"/> Accessory Bag Personalizing 2 Initials \$27.30 <input type="radio"/> Cup Holder \$48.30	<input type="radio"/> Sewn Strap <input type="radio"/> Front Mounted OR <input type="radio"/> Rear Mounted \$173.25 <input type="radio"/> Brush Font OR <input type="radio"/> Block Font \$27.30

Matrix® Seating Series

Matrix® Seat Cushions

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" (Only available w/ PS Cushion)	Seat Cushion Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"
<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown \$745 <input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$624 <input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown \$643 <input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$569	<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown \$1,025 <input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$681 <input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown \$643 <input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$975

Libra Cushion Accessories: <input type="radio"/> Libra Fluid Supplement \$120 <input type="radio"/> Libra 1" Hemi Wedge \$101 <input type="radio"/> Libra 1" Anterior Seat Wedge \$202 <input type="radio"/> Libra Tapered Wedge Left \$85 <input type="radio"/> Libra Tapered Wedge Right \$85 <input type="radio"/> Libra Abductor \$85	Vi, PSP & PS Cushion Accessories: <input type="radio"/> Pelvic Obliquity Kit \$85 Cushion Rigidizer: <input type="radio"/> 16" Depth <input type="radio"/> 18" Depth <input type="radio"/> 20" / 22" Depth <table border="1" style="width:100%"> <tr> <th colspan="2">All Cushions</th> <th>16" Width</th> <th>18" Width</th> <th>20" Width</th> <th>22" Width</th> </tr> <tr> <td></td> <td>16" Depth</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> </tr> <tr> <td></td> <td>18" Depth</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> </tr> <tr> <td></td> <td>20" Depth</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> </tr> <tr> <td></td> <td>22" Depth</td> <td>N/A</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> <td><input type="radio"/> \$115</td> </tr> </table>	All Cushions		16" Width	18" Width	20" Width	22" Width		16" Depth	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115		18" Depth	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115		20" Depth	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115		22" Depth	N/A	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
All Cushions		16" Width	18" Width	20" Width	22" Width																										
	16" Depth	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115																										
	18" Depth	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115																										
	20" Depth	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115																										
	22" Depth	N/A	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115																										



Matrix® Backs

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height: <input type="radio"/> 16" <input type="radio"/> 18" <i>(excludes PB & PB Deep)</i> <input type="radio"/> 20"	
Back Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	
<input type="radio"/> E2 - Recommended for users seeking optimal immersion and airflow	\$965
<input type="radio"/> E2 Deep - Recommended for users seeking optimal immersion w/ contoured support needs	\$1,125
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications	\$866
<input type="radio"/> Elite - Recommended for more active users	\$766
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs	\$918
<input type="radio"/> PB - Recommended for less active users seeking comfort	\$734
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support	\$892
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications	\$1,151
<input type="radio"/> Elite - Recommended for more active users	\$1,051
<input type="radio"/> Elite Deep - Recommended for active users with contoured support	\$1,095
Matrix Back	<input type="radio"/> Air Force
Military	<input type="radio"/> Army
Service	<input type="radio"/> Marine
Decals:	<input type="radio"/> Navy
Privacy Flap	<input type="radio"/> 16" Width Small Privacy Flap
Accessories:	<input type="radio"/> 18" Width Medium Privacy Flap
	<input type="radio"/> 20" Width Large Privacy Flap
	<input type="radio"/> 22" Width Extra Large Privacy Flap
Lateral Support	<input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L
Accessories:	<input type="radio"/> Elite Back Swing Away Lateral Hardware
	<input type="radio"/> Swing Away Lateral Pad (Each) 4.25"Hx6.5"L
E2 Back Accessories:	<input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair)
	<input type="radio"/> E2 Chest Harness Interface
	<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support
	<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32"L
	<input type="radio"/> Quad Release Tab for PostureFlex
Back Slipcovers:	<input type="radio"/> E2 Infection Control Slipcover
	<input type="radio"/> E2 Deep Infection Control Slipcover
	<input type="radio"/> Infection Control Slipcover for Matrix Cushions

Matrix® Headrest

Head Support Mounting:	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support	\$285
Choose One	<input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support	\$275
Head Support Options:	<input type="radio"/> Head Support	\$310
	<input type="radio"/> 4 Point Head Support	\$448
	<input type="radio"/> Headband with 4 Point Cushion	\$610

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
 - Residential Delivery Service:** Delivery to a residential location
 - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
 - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
 - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.*

Knee-to-Heel Range (Comfort Rehab ONLY)			
Front Rigging Model	Description	Minimum	Maximum
Range Adjusts in 3/4" increments			
F-SW 75-587V3	75" Drop-in, black	14.25"	17.25"
F-SW 75ADJ-587V3	75" Drop-in w/ angle adjustable footplates, black	14.25"	17.25"
F-SW 60ADJ-587V3	60" Drop-in w/ angle adjustable footplates, black	15"	18"
F-SW-ELV-587V3	Drop-in, elevating, black	15"	18"
F-SW-ELV-ANGLE-587V3	Drop-in, elevating w/ angle adjustable footplates, black	15"	18"

587V3: SEAT HEIGHT & WHEEL COMBINATIONS									
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB	
14"	20	5	SMALL	1	1	TOP	3	2	
	22	5	SMALL	2	2	TOP	2	1	
15.5"	20	6	MEDIUM	2	1	TOP	2	1	
	22	6	MEDIUM	1	1	TOP	2	1	
	24	6	SMALL	1	2	TOP	1	1	
	24	5	MEDIUM	1	2	TOP	1	1	
16"	24	6	MEDIUM	1	1	TOP	1	1	
	20	5	SMALL	3	3	TOP	1	1	
	20	5	MEDIUM	3	3	TOP	1	1	
	20	6	MEDIUM	3	2	TOP	1	1	
	22	5	SMALL	2	3	TOP	1	1	
	22	6	SMALL	2	2	TOP	1	1	
16.5"	22	5	MEDIUM	2	3	TOP	1	1	
	22	6	MEDIUM	2	2	TOP	1	1	
	24	5	SMALL	1	3	TOP	1	1	
	20	6	SMALL	3	3	TOP	1	1	
17"	24	6	SMALL	2	3	BOTTOM	3	1	
	24	5	MEDIUM	2	4	BOTTOM	3	1	
	24	6	MEDIUM	2	3	BOTTOM	3	1	
	20	5	MEDIUM	4	5	BOTTOM	3	1	
17.5"	20	6	MEDIUM	4	4	BOTTOM	3	1	
	22	5	MEDIUM	3	5	BOTTOM	3	1	
	22	6	MEDIUM	3	4	BOTTOM	3	1	
	20	8	SMALL	4	3	BOTTOM	3	1	
18"	22	8	SMALL	3	3	BOTTOM	3	1	
	24	8	SMALL	3	3	BOTTOM	2	1	
	24	6	MEDIUM	3	5	BOTTOM	2	1	
	24	8	MEDIUM	3	3	BOTTOM	2	1	
19"	20	8	MEDIUM	5	4	BOTTOM	2	1	
	22	8	MEDIUM	4	4	BOTTOM	2	1	
19"	24	8	MEDIUM	4	5	BOTTOM	1	1	

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.
Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels
Form08CONFIG- 404 587V3 Rev. 1