

Comfort

Rehab Wheelchair

REHAB
BRODA



Shown with optional head support, Matrix Elite back support, Matrix Vi seat cushion, and swing-away lower leg supports

Patient, Facility, & Caregiver Benefits

- Pressure redistribution that aids in maintaining skin integrity
- Durable frame to withstand heavy workload healthcare environments
- Adaptable for wide range of patients and care needs



Front pivot seat tilt maintains user comfort with smooth, effortless tilt



Posterior tilt optimizes pressure redistribution



Solid seat pan offers support and is compatible with a variety of seat cushion options



Range of simple adjustments accommodates a variety of needs

Matrix® seating options available to ship preconfigured on-chair.

matrix® SEATING SERIES

Elan Headrest



Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.

Elite TR Back



Shell – Contoured aluminum with cutouts for lifting and breathability. Available in standard (3”) OR deep (6”) contours.

Foam – HR (high resilient) foam enclosed in cover.

Libra Cushion




Foam – HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.

Simplified Durability for Rehabilitation Care


SPECIFICATIONS

Degrees of Tilt	Up to 45°
Weight Capacity	Up to 250 lbs. (114 kg)
Angle Adjustable Back Canes	5°, 10°, 15°, 20°
Seat Widths	16”, 18”, 20”, 22” (40cm, 45cm, 50cm, 55cm)
Seat Depths	Adjustable: 16”, 17”, 18”, 19”, 20 (40cm, 43cm, 45cm, 48cm, 50cm)
Seat to Surface Heights	Adjustable: 14”, 15”, 15.5”, 16”, 16.5”, 17”, 17.5”, 18”, 19” (35cm, 38cm, 39cm, 42cm, 43cm, 44cm, 45cm, 48cm)
Wheels	5”, 6” or 8” Front Caster Wheels with fork (13cm, 15cm, 20cm) 20”, 22” or 24” Rear Mag Wheels (50cm, 55cm, 60cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Government Coding (US & CA)	HCPCS: E1161 ADP (ON): WMA5B-025

FRAME OPTIONS



BLACK



BRONZE



TUNGSTEN
SAND




SATIN
NICKEL



CARMINE
RED



DEEP
ORANGE



ZINC
YELLOW



OCEAN
BLUE

Additional accessories available: lower leg supports, trays, and positioning belts