

<b>PO:</b> _____	<b>Date:</b> _____	<b>Customer #:</b> _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
<ul style="list-style-type: none"> <li>• Solid Seat Pan</li> <li>• Tubular Steel, Powder Coated Frame</li> <li>• Seat cushion, back &amp; head support not included</li> <li>• Independent Height Adjustable Arm Supports</li> <li>• Rocking Mechanism (<i>allows up to 33° of tilt</i>)</li> <li>• Fold Down and Angle Adjustable Back Canes</li> <li>• 33° Infinitely Adjustable Seat Tilt</li> </ul>	<ul style="list-style-type: none"> <li>• 250lb Maximum Occupant Weight</li> <li>• Adjustable Seat Surface Height</li> <li>• 20" Front Mount Mag Wheels</li> <li>• 2 x 5" Rear Total Lock Caster Wheels</li> <li>• Lower Leg Support Assembly</li> <li>• Removable Arm Supports</li> </ul>

**Encore Rehab Wheelchair® - Model: 500RW (E1161) MSRP \$3,499.65**

Seat and Back Support Options			
<b>Seat Width:</b>	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" (E2201) <input type="radio"/> 22" (E2201)	N/C \$420 \$420	<b>Seat Depth:</b> <input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20" (E2203)
			<input type="radio"/> Solid Seat Pan (E2231) \$157.50
			<input type="radio"/> Omit Seat Pan N/C
<b>Seat Surface</b>	<input type="radio"/> 14" <input type="radio"/> 15" <input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19"	N/C Standard N/C	<b>Back Cane</b> <input type="radio"/> Ultra-Low Double Bend (16") Fixed Standard
<b>Height:</b>			<b>Height:</b> <input type="radio"/> Low (18") Fold Down N/C
			<input type="radio"/> Medium (21") Fold Down N/C
			<input type="radio"/> Tall (24") Fold Down N/C

Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair Standard	<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels - Pair \$205.80
<input type="radio"/> 22" Mag Wheels - Pair N/C	<input type="radio"/> Front Anti-Tip Devices (E0971) Required Safety Fee \$105.00	<input type="radio"/> Parking Brake Extension Handles - Pair (E0961) \$69.30
<input type="radio"/> Quick Release Axles (Pair) \$126.00		

Arm Support Options			
<b>Arm Support</b> <input type="radio"/> Arm Supports - Pair Standard	<b>Arm Support</b> <input type="radio"/> Height Adjustable Arm Supports (E0973) \$149.10	<b>Style:</b> <input type="radio"/> *Height and Depth Adjustable Arm Supp \$149.10	
<b>Cushion</b> <input type="radio"/> Arm Support Trough - Each \$157.38	<b>Options:</b> <input type="radio"/> Cushioned Arm Supports* -Pair \$86.10	*Not compatible with 22" mag wheels in the 14", 15", & 16" seat heights	
*Not compatible with a tray		<b>Arm Support Mount:</b> <input checked="" type="radio"/> Seat Mounted Arm Supports Standard	

Frame Color Options			
<b>Frame Color:</b> <input type="radio"/> Black <input type="radio"/> Orange <input type="radio"/> Carmine Red	<input type="radio"/> Tungsten Sand <input type="radio"/> Ocean Blue	<input type="radio"/> Bronze <input type="radio"/> Satin Nickel	<input type="radio"/> Yellow Zinc

Drop In Style	Must Choose Lower Leg Support Assembly Style	Pin and Saddle Style
<input type="radio"/> SM-HANGER-48RV4-RKR-IKN (pair) N/C	<input type="radio"/> SM-HANGER-48RV4-RKR-PIN (pair) Less \$54	N/C
Lower Leg Support Options		
<input type="radio"/> F-SW-IKN Swing Away Lower Leg Support Silver - Pair \$102.90	<input type="radio"/> F-SW Swing Away Lower Leg Support Black - Pair \$222.60	
<input type="radio"/> F-SW-IKN ADJ Swing Away Lower Leg Support Black - Pair \$211.05	<input type="radio"/> F-SW-ELV Elevating Swing Lower Leg Support Black - Pair (E0990) \$279.03	
w/ Angle Adjustable Flip-up Foot Supports (K0040) \$105	<input type="radio"/> F-SW-ELV-ANGLE Elevating Swing Away Lower Leg Supports (pair) (E0995) \$220.50	
<input type="radio"/> F-SW-ELV-IKN Elevating Swing Away Lower Leg Support Silver (E0990) \$309.75	w/ Angle Adjustable Flip-up Foot Supports(K0040) \$105.00	
w/ Angle Adjustable Flip-up Foot Supports(K0040) \$105	<input type="radio"/> F-SW-ANGLE Swin Away 70 Degree Lower Leg Support \$222.60	
<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support Panel (E0995) \$86.10	w/ Angle Adjustable Footplates (K00040) \$105.00	
	<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support Panel (E0995) \$86.10	

\*Slight ground interference may occur in the lowest seat surface height. Slightly tilt chair or elevate lower leg support as needed

2" Extended Footrest Tubes - *Only compatible on the F-SW-IKN and F-SW-IKN-ADJ* \$67.41

Knee-to-Heel Range			
Front Rigging Mod	Description	Minimum	Maximum
Range Adjusts in 1" increments			
F-SW-IKN	70° Drop-in, silver	13"	17" (19" with optional extender)
F-SW-IKN-ADJ	70° Drop-in w/ angle adjustable footplates, black	13"	17" (19" with optional extender)
Range Adjusts in 3/4" increments			
F-SW	70° Pin & Saddle, black	11"	14"
F-SW-ANGLE	70° Pin & Saddle w angle adjustable footplates, black	11"	14.75"
F-SW-ELV	Pin & Saddle, elevating, black	15.25"	18.25"
F-SW-ELV-ANGLE	Pin & Saddle, elevating w/ angle adjustable footplates, black	15.25"	18.25"
Infinite Range Adjustments			
F-SW-ELV-IKN	Drop-in, elevating w/ angle adjustable footplates, silver	13"	20"

Common Accessories				
<b>Pelvic Belts:</b>	<input type="radio"/> S-BELT	Pelvic Belt (E0978)	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle	\$56.70
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle (E0978)	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle	\$56.70
	<input type="radio"/> S-BELT-587V3	Pelvic Belt Grommet Mount (E0978)	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle	\$56.70
<b>Trays:</b>	<input type="radio"/> TRAY-ABS	Tray ABS with Strap (E0950)	<input type="radio"/> Sewn Strap <input type="radio"/> Grey <input type="radio"/> Black	\$337.05
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width) (E0950)	<input type="radio"/> Black Only	\$518.44
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap (E0950)	<input type="radio"/> Sewn Tray Strap	\$348.60
	<input type="radio"/> TRAY-HALF	Half Tray (E0950) (side specific)(with upper side pad)	<input type="radio"/> Grey <input type="radio"/> Black	\$463.05
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover		\$171.15
<b>Other:</b>	<input type="radio"/> IV-POLE (K0105)	IV Pole and Assembly - <input type="radio"/> Front or <input type="radio"/> Rear Mount		\$108.15
	<input type="radio"/> O2-HOLDER (E2208)	Oxygen Holder and Assembly <input type="radio"/> Front or <input type="radio"/> Rear Mount		\$181.65
	<input type="radio"/> STABIL-FR	Rear Stabilizer Frame		\$97.65
<b>Positioning:</b>	<input type="radio"/> SR	Side Reduction Cushions - Pair (E0956)		\$108.15

**Matrx® Seating Series**



**Matrx® Seat Cushions**

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

<b>Seat Cushion Depth:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	<input type="radio"/> 22" (Only available w/ PS Cushion)
----------------------------	---------------------------	---------------------------	---------------------------	--

<b>Seat Cushion Width:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"
<input type="radio"/> Libra (E2624) - Recommended for highest risk of or existing skin l	\$485		
<input type="radio"/> Vi (E2607) - Recommended for moderate to high risk of skin bre:	\$425		
<input type="radio"/> PSP (E2607) - Recommended for moderate to high risk of skin bi	\$395		
<input type="radio"/> PS (E2605) - Recommended for low to moderate risk of skin bre:	\$345		

<b>Seat Cushion Width:</b>	<input type="radio"/> 22"
<input type="radio"/> Libra (E2625) - Recommended for highest risk of or existing skin break	\$595
<input type="radio"/> Vi (E2608) - Recommended for moderate to high risk of skin breakdown:	\$475
<input type="radio"/> PS (E2606) - Recommended for low to moderate risk of skin breakdown:	\$545

<b>Libra Cushion</b>	<input type="radio"/> Libra Fluid Supplement	\$75
<b>Accessories:</b>	<input type="radio"/> Libra 1" Hemi Wedge	\$65
	<input type="radio"/> Libra 1" Anterior Seat Wedge	\$125
	<input type="radio"/> Libra Tapered Wedge Left	\$50
	<input type="radio"/> Libra Tapered Wedge Right	\$50
	<input type="radio"/> Libra Abductor	\$60

<b>Vi, PSP &amp; PS Cushion</b>	
<b>Accessories:</b>	<input type="radio"/> Pelvic Obliquity Kit \$50

<b>Cushion Rigidizer:</b>	<b>E0992</b>	16" Depth	18" Depth	20" / 22" Depth
<b>All Cushions</b>	16" Width	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
	18" Width	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
	20" Width	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
	22" Width	N/A	<input type="radio"/> \$115	<input type="radio"/> \$115

**Matrix® Backs**

*For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions*

<b>Back Height:</b>	<input type="radio"/> 16" <input type="radio"/> 18"(excludes PB & PB Deep) <input type="radio"/> 20"
<b>Back Width:</b>	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"
<input type="radio"/> E2 (E2615) - Recommended for users seeking optimal immersion and airflow	\$645
<input type="radio"/> E2 Deep (E2620) - Recommended for users seeking optimal immersion w/ contoured support needs	\$775
<input type="radio"/> Elite TR (E2620) - Recommended for tilt in space & recline applications	\$650
<input type="radio"/> Elite (E2620) - Recommended for more active users	\$625
<input type="radio"/> Elite Deep (E2620) - Recommended for active users with contoured support needs	\$760
<input type="radio"/> PB (E2615) - Recommended for less active users seeking comfort	\$625
<input type="radio"/> PB Deep (E2615) - Recommended for less active users seeking comfort & support	\$645
<b>E2 Back</b>	<input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair) (E0956) \$290
<b>Accessories:</b>	<input type="radio"/> E2 Chest Harness Interface \$75
	<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support \$55
	<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32"L \$245
	<input type="radio"/> Quad Release Tab for PostureFlex \$25
<b>Back</b>	<input type="radio"/> E2 Infection Control Slipcover \$65
<b>Slipcovers:</b>	<input type="radio"/> E2 Deep Infection Control Slipcover \$65
	<input type="radio"/> Infection Control Slipcover for Matrix Cushions \$65
<b>Back Width:</b>	<input type="radio"/> 22"
<input type="radio"/> Elite TR (E2621) - Recommended for tilt in space & recline application	\$845
<input type="radio"/> Elite (E2621) - Recommended for more active users	\$795
<input type="radio"/> Elite Deep (E2621) - Recommended for active users with contoured sup	\$930
<b>Matrix Back</b>	<input type="radio"/> Air Force N/C
<b>Military Service</b>	<input type="radio"/> Army N/C
<b>Decals:</b>	<input type="radio"/> Marine N/C
	<input type="radio"/> Navy N/C
<b>Privacy Flap</b>	<input type="radio"/> 16" Width Small Privacy Flap \$40
<b>Accessories:</b>	<input type="radio"/> 18" Width Medium Privacy Flap \$40
	<input type="radio"/> 20" Width Large Privacy Flap \$40
	<input type="radio"/> 22" Width Extra Large Privacy Flap \$40
	<b>Lateral Support</b>
	<input type="radio"/> Elite Back Swing Away Lateral Hardware (E0128) \$22
<b>Accessories:</b>	<input type="radio"/> Swing Away Lateral Pad (Each) 4.25"Hx6.5"L (E0956) \$1

**Matrix® Headrest**

<b>Head Support Mounting:</b>	<input type="radio"/> Elan Head Support Hardware (E0128) - Required with Selection of Head Support \$220
<b>Choose One</b>	<input type="radio"/> Loxx Head Support Hardware (E0128) - Required with Selection of Head Support \$215
<b>Support Options:</b>	<input type="radio"/> Head Support Cushion (E0955) <input type="radio"/> 10" <input type="radio"/> 14" \$235
	<input type="radio"/> Head Support 4 Point Cushion (E0955) <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10" \$325
	<input type="radio"/> Headband with 4 Point Cushion (E0955) <input type="radio"/> Small <input type="radio"/> Large \$450

**Additional Options:**

**Delivery Instructions or Special Requests:**

Total: \$

**Accessorial Services:** Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

**any accessorial services needed to deliver are requested when an order is submitted to Broda.**

**Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

**Residential Delivery Service:** Delivery to a residential location

**Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

**Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building,

upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

**Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

**Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**

