



Shown with Matrix Libra cushion, Matrix Elite TR back, Matrix Elan Headrest, and swing-away footrests options. Frame in Carmine Red

Patient, Facility, and Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning



Back canes fold down for easy storage options



Dynamic rocking function soothes and stimulates users



Front pivot seat tilt maintains a proper position for self-propulsion



Open seating system allows for almost any seating configuration

matrix®

SEATING SERIES

The Encore of Comfort, Functionality, and Seating Customization



Elan Headrest

Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.



Elite TR Back

Shell - Contoured aluminum with cutouts for lifting and breathability. Available in standard (3") OR deep (6") contours.

Foam - HR (high resilient) foam enclosed in cover.




Libra Cushion

Foam - HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.


SPECIFICATIONS

Degrees of Tilt	Up to 33° (with Dynamic Rocking)
Weight Capacity	Up to 250 lbs. (114 kg)
Seat Widths	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
Seat Depths	16", 17", 18", 19", 20" (40cm, 43cm, 45cm, 48cm, 50cm)
Back Cane Heights	16", 18", 21", 24" (40cm, 45cm, 53cm, 60cm) Angle adjustable to 5°, 10°, 15°, 20°
Adjustable Floor-to-Seat Heights	14", 15", 16", 17", 18", 19" (35cm, 38cm, 40cm, 43cm, 45cm, 48cm)
Wheels	20" or 22" Front Mag (50cm or 55cm) 5" Rear Casters: Total Lock (12.5cm) 5" Rear Aluminum: Non-Locking (12.5cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame

FRAME OPTIONS



BLACK




BRONZE


TUNGSTEN
SAND

SATIN
NICKEL

CARMINE
RED

DEEP
ORANGE

ZINC
YELLOW

OCEAN
BLUE

Additional accessories available: foot supports, trays, and positioning straps