

# BRODA<sup>®</sup>

# Encore

## Rehab Wheelchair



mobility  
product award 2021

HCPCS  
E1161



*Back canes fold down for easy storage options*



*Dynamic rocking function soothes and stimulates users*



*Front pivot seat tilt maintains a proper position for self-propulsion*



*Open seating system allows for almost any seating configuration*

*Shown with Matrix Libra cushion, Matrix E2 back, and Matrix Loxx Head support*

## **Patient, Facility, and Caregiver Benefits**

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning

**brodaseating.com | 1-844-552-7632**

### matrix<sup>®</sup>

SEATING SERIES



#### Loxx Headrest

Highly adjustable hardware with tool-less adjustment levers for easy fitting and position changes.



#### E2 Back

**Shell** – Contoured aluminum with cutouts for lifting and breathability. Available in standard (3”) OR deep (6”) contours.

**Foam** – HR (high resilient) foam enclosed in cover.



#### Libra Cushion



**Foam** – HR (high resilient) molded foam made with UltraFresh antimicrobial for odor protection. Contoured pelvic well, waterfall front edge.

MOBILITY WITH COMFORT,  
FUNCTIONALITY, &  
CUSTOMIZED SEATING

### SPECIFICATIONS

<b>Degrees of Tilt</b>	Up to 33° (with Dynamic Rocking)
<b>Weight Capacity</b>	Up to 250 lbs. (114 kg)
<b>Seat Widths</b>	16”, 18”, 20”, 22” (40cm, 45cm, 50cm, 55cm)
<b>Seat Depths</b>	16”, 17”, 18”, 19”, 20” (40cm, 43cm, 45cm, 48cm, 50cm)
<b>Back Cane Heights</b>	16”, 18”, 21”, 24” (40cm, 45cm, 53cm, 60cm) Angle adjustable to 5°, 10°, 15°, 20°
<b>Adjustable Floor-to-Seat Heights</b>	14”, 15”, 16”, 17”, 18”, 19” (35cm, 38cm, 40cm, 43cm, 45cm, 48cm)
<b>Wheels</b>	20” or 22” Front Mag (50cm or 55cm) 5” Rear Casters: Total Lock (12.5cm) 5” Rear Aluminum: Non-Locking (12.5cm)
<b>Lifetime Warranty Frame</b>	10 years on 16 Gauge tubular steel frame

### FRAME OPTIONS



Additional accessories available: foot supports, trays, and positioning straps

[brodaseating.com](http://brodaseating.com) | 1-844-552-7632