

PO: _____

Date: _____

Customer #: _____

Customer Name: _____
 Address: _____
 Phone: _____
 Contact: _____
 Email: _____

Ship To Name: _____
 Address: _____
 Phone: _____
 Contact: _____
 Tag: _____



Standard Features	
• 250lb Weight Capacity	• 2 x 5" Total Lock Rear Casters
• Caregiver Activated Rocking	• 2 x 20" Front Mounted Mag Wheels
• Comfort Tension Seating®	• Seat Mounted Footrest Hangers (standard with 18.5" or 20" seat depth)
• 14" - 19" Adjustable Seat to Floor Height	• Base Mounted Footrest Hangers (standard with 15.5" or 17" seat depth)
• 32" Back Height	• Swing Away Footrests
• Up to 33° of Seat Tilt	

Latitude Rehab Wheelchair - Model: 500 MR MSRP \$3,164.95

Chair Color Options

Strap Color: Grey Blue Black N/C Frame Color: Grey Black N/C

Upholstery Color Options (Select one)

Traditional Upholstery Colors: Burgundy Charcoal Sapphire Blue

*Traditional Upholstery Colors come with seat and back pad in matching vinyl color.

New Upholstery Colors:

Vinyl:

Indigo*** Espresso*** Graphite***
 Black Smoke House Chocolate

Polyurethane:

Portobello Geyser Fossil Forge

Camo Nylon: (with performance finish)

Digital Camo Green Retro Camo Jungle \$150.43

Upholstery Colors come with seat and back pad in black vinyl. *Indigo, Espresso, and Graphite available with seat and back pad in matching vinyl color

Seat Options

Seat Width: <input type="radio"/> 16" or <input type="radio"/> 18" N/C	Seat Depth: <input type="radio"/> 15.5" <input type="radio"/> 20" \$109.94
<input type="radio"/> 20" N/C	<input type="radio"/> 17" N/C
<input type="radio"/> 22" \$457.10	<input type="radio"/> 18.5" Standard

Seat Height: 14" 15" 16" 17" 18" 19" N/C

Seat Option: Comfort Tension Seating® (CTS) Standard Solid Seat Pan (Not compatible with 15.5" seat depth) \$157.50

See CTS Seat Options below

See Matrx Seat Cushion Options below

Comfort Tension Seating® (CTS) Seat Padding Options

Broda Seat Padding Includes: Vinyl Cover Terry Cover (Grey or White) N/C

Full Vinyl Wrap Seat \$54.98

Additional Positioning Padding (APP) Seat Cushion 16" 18" 20" Seat Depth 17" 19" 21" \$225.66

Seat Width: 22" 17" 19" 21" \$238.38

Solid Seat Pan: Matrx® Seat Cushions

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions



Seat Cushion Depth:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" (Only available w/ PS Cushion)	
Seat Cushion Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Cushion Width: <input type="radio"/> 22"
<input type="radio"/> Libra-Recommend for highest risk of or existing skin breakdown	\$745	<input type="radio"/> Libra-Recommend for highest risk of or existing skin breakdown \$1,025
<input type="radio"/> Vi-Recommend for moderate to high risk of skin breakdown	\$624	<input type="radio"/> Vi-Recommend for moderate to high risk of skin breakdown \$681
<input type="radio"/> PSP-Recommend for moderate to high risk of skin breakdown	\$645	
<input type="radio"/> PS - Recommend for low to moderate risk of skin breakdown	\$569	<input type="radio"/> PS -Recommend for low to moderate risk of skin breakdown\$975

Libra Cushion Accessories: <input type="radio"/> Libra Fluid Supplement \$120 <input type="radio"/> Libra 1" Hemi Wedge \$101 <input type="radio"/> Libra 1" Anterior Seat Wedge \$202 <input type="radio"/> Libra Tapered Wedge Left \$85 <input type="radio"/> Libra Tapered Wedge Right \$85 <input type="radio"/> Libra Abductor \$85	Vi, PSP & PS Cushion Accessories: <input type="radio"/> Pelvic Obliquity Kit \$85
	Cushion Rigidizer: 16" D 18" D 20" / 22" D For All Cushions
	16" Width <input type="radio"/> \$115 <input type="radio"/> \$115 <input type="radio"/> \$115
	18" Width <input type="radio"/> \$115 <input type="radio"/> \$115 <input type="radio"/> \$115
	20" Width <input type="radio"/> \$115 <input type="radio"/> \$115 <input type="radio"/> \$115
	22" Width N/A <input type="radio"/> \$115 <input type="radio"/> \$115

<input type="radio"/> Infection Control Slipcover for Matrx Cushion	\$95
Back Options	
Back Height: <input type="radio"/> 32" Standard	Back Seating <input checked="" type="radio"/> Comfort Tension Seating® (CTS) Standard
<input type="radio"/> 33.5" \$109.94	
<input type="radio"/> 35" \$109.94	

Back Padding Options	
<input checked="" type="radio"/> Broda Padding Includes: Headrest, Shoulder Bolsters and Wings, Lower Side Pads, Back Pad with Cover	N/C
<input type="radio"/> Vinyl Back <input type="radio"/> Terry Back (Grey or White)	N/C
<input type="radio"/> Full Vinyl Wrap Back	\$54.98
<input type="radio"/> Additional Positioning Padding (APP) Back Cushion	Back Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" \$168.96
	<input type="radio"/> 22" \$217.55

Common Padding Upgrades	
Headrests: <input type="radio"/> NR-CV Headrest Visco Foam Contoured Pad <input type="radio"/> Triptex <input type="radio"/> Omit Standard	\$115.72
<input type="radio"/> NS Neck Support Soft Visco Cushion	\$115.72
Positioning: <input type="radio"/> SR Side Reduction Pads (pair)	\$109.94
<input type="radio"/> SAS Side Arm Supports (pair)	\$135.45
<input type="radio"/> CC-LTS Comfort Company Lateral Supports	\$516.12
<input type="radio"/> LTS Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)	\$636.47

Armrest Options (Select One)	
Armrest Options: <input type="radio"/> Height Adjustable Non-removable	Standard
Armtop Style: <input type="radio"/> Soft Armrests	Standard
<input type="radio"/> Padded Armrests* (pair)	\$86.79
<input type="radio"/> Arm Trough (each)	\$150.43

*Not compatible with a tray

Front Mount Configuration	Wheel Options	Rear Mount Configuration
<input type="radio"/> 20" Front Mount Mag Wheels Standard	<input type="radio"/> 5" Total Lock Rear Casters Standard	
<input type="radio"/> 22" Front Mount Mag Wheels N/C		
<input type="radio"/> 5" Total Lock Front Casters (when ordered with rear mags)	<input type="radio"/> 20" Rear Mags N/C	
<input type="radio"/> 5" Front Aluminum Casters \$219.87	<input type="radio"/> 22" Rear Mags* N/C	
<input type="radio"/> Wheel Lock Extension Handles (pair) \$74.06		

Anti-Tippers (*factory required safety feature*) N/C

Footrest Options

Footrest (allows for footrests to tilt with chair)(configuration specific - please contact Customer Support for full details)

Hangers:	<input type="checkbox"/> SM-HANGER	Flared Seat Mount - Short		\$127.29
	<input type="checkbox"/> SM-HANGER	Flared Seat Mount - Long		\$127.29
	<input type="checkbox"/> SM-HANGER	Straight Seat Mount - Short		\$127.29
	<input type="checkbox"/> SM-HANGER	Straight Seat Mount - Long		\$127.29
Footrests:	<input type="checkbox"/> F-SW-ELV	Elevating Swing Away Footrests (<i>pair</i>)		\$283.51
	<input type="checkbox"/> F-SW	Swing Away Footrests (<i>pair</i>)		\$185.16
	<input type="checkbox"/> LEG-FOOT-48/48R	Leg and Foot Board		\$260.37
	<input type="checkbox"/> LEGPANEL-48	Largest Panel Pad Black		\$86.79

Common Accessories

Seat Belts:	<input type="checkbox"/> S-BELT	Pelvic Positioning Strap	<input type="checkbox"/> Center Buckle	<input type="checkbox"/> Side Buckle	\$57.87
	<input type="checkbox"/> S-BELT-A	Pelvic Positioning Strap Auto Buckle	<input type="checkbox"/> Center Buckle	<input type="checkbox"/> Side Buckle	\$57.87
	<input type="checkbox"/> S-BELT-PD	Pelvic Positioning Strap Padded			\$83.32
	<input type="checkbox"/> S-BELT-PD-A	Pelvic Positioning Strap Padded Auto Buckle			\$83.32
Thigh Belts:	<input type="checkbox"/> T-BELT	Padded Thigh Belt			\$133.08
(not compatible	<input type="checkbox"/> T-BELT-A	Padded Thigh Belt Auto Buckle			\$138.86
with Broda APP	<input type="checkbox"/> T-BELT-NEO	Padded Thigh Belt Neoprene			\$133.08
padding option)	<input type="checkbox"/> T-BELT-NEO-A	Padded Thigh Belt Neoprene Auto Buckle			\$138.86
Trays:	<input type="checkbox"/> TRAY-ABS	Tray ABS with Strap	<input type="checkbox"/> Sewn Tray Strap	<input type="checkbox"/> Grey <input type="checkbox"/> Black	\$341.38
	<input type="checkbox"/> TRAY-CLEAR	Tray Clear with Strap	<input type="checkbox"/> Sewn Tray Strap		\$352.95
	<input type="checkbox"/> TRAY-HALF	Half Tray (<i>side specific</i>)(includes upper side pad)			\$468.67
	<input type="checkbox"/> TRAY-COVER	Tray Cover (<i>padded</i>)(for TRAY-ABS)			\$173.59
Stabilizers:	<input type="checkbox"/> STABIL-FR	<input type="checkbox"/> Front or <input type="checkbox"/> Rear Mounted Stabilizers (<i>select one</i>)(<i>pair</i>)			\$98.36
	<input type="checkbox"/> STABIL-LAT	Lateral Frame Stabilizers			\$173.59
Other:	<input type="checkbox"/> IV-POLE	IV Pole & Holder (<i>includes mounting bracket</i>) - <input type="checkbox"/> Front or <input type="checkbox"/> Rear Mount			\$109.94
	<input type="checkbox"/> O2-HOLDER	O2 Holder (<i>includes mounting bracket</i>) - <input type="checkbox"/> Front or <input type="checkbox"/> Rear Mount			\$185.16
	<input type="checkbox"/> ACCESSORY-BAG	Accessory Bag (Black) Each			\$46.28
	<input type="checkbox"/> BAG PERSONALIZIN(2 Initials)	<input type="checkbox"/> Brush Fon OR <input type="checkbox"/> Block Font			\$28.93
	<input type="checkbox"/> CUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)			\$52.07

Delivery Instruction or

Special Requests:

Total: \$



Latitude Rehab Wheelchair 500MR

AU Order Form
January 2022



Latitude Rehab Wheelchair 500MR

AU Order Form
January 2022