



*Dynamic Rocking* feature helps calm agitation and absorb rigid movements




*Posterior tilt* improves posture, redistributes pressure and aids in fall prevention



*Front pivot seat tilt* enhances visual orientation, and encourages self-mobility with proper foot on floor positioning




### **Patient, Facility, & Caregiver Benefits**

- Tilt enables ease of use for frequent repositioning
- Designed for safety and stability helping to reduce risk of injury
- Provides stimulation and contentment - excellent for memory care



*Comfort Tension Seating<sup>®</sup>* improves postural support and sitting tolerance

## UPHOLSTERY OPTIONS




## Enhanced Mobility and Freedom for Memory Care

## SPECIFICATIONS

Degrees of Tilt	up to 33° (with dynamic rocking)
Weight Capacity	250 lbs. (113 kg)
Seat Widths	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
Seat Depths	15.5", 17", 18.5" (39cm, 43cm, 46cm)
Floor-to-Seat Heights	14", 15", 16", 17", 18", 19", (35cm, 38cm, 40cm, 43cm, 45cm, 48cm)
Wheels	5" (12cm) Aluminum Front Casters 5" (12cm) Total Lock Front Casters 20" or 22" (35cm or 55cm) Front Mags 5" (15cm) Total Lock Rear Casters 20" or 22" (35cm or 55cm) Rear Mags
Lifetime Frame Warranty	10 years on 16 Gauge tubular steel frame
Standard Padding	Upper & Lower Side Pads, Back Pad with Cover

## FRAME OPTIONS



## CUSHION OPTIONS

**matrix**®  
SEATING SERIES



**Libra:** Contoured foam prevents sliding and provides support.



**Vi:** Foam base provides stability and prevents bottoming out.

Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps