

BRODA[®]

Latitude

Rehab Wheelchair



Dynamic rocking feature helps calm agitation and absorb rigid movements



Posterior tilt improves posture and redistributes pressure



Front pivot seat tilt enhances visual orientation, and encourages self-mobility with proper foot on floor positioning



Comfort Tension Seating[®] improves postural support and sitting tolerance


Patient, Facility, & Caregiver Benefits

- Tilt enables ease of use for frequent repositioning
- Designed for safety and stability helping to reduce risk of injury
- Provides stimulation and contentment - excellent for memory care

brodaseating.com | 1-844-552-7632

ENHANCED MOBILITY &
FREEDOM FOR MEMORY CARE

UPHOLSTERY OPTIONS



Additional accessories available: foot supports, trays, and positioning belts

SPECIFICATIONS

Degrees of Tilt	up to 33° (with dynamic rocking)
Weight Capacity	250 lbs. (113 kg)
Seat Widths	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
Seat Depths	15.5", 17", 18.5" (39cm, 43cm, 46cm)
Seat Surface Heights	14", 15", 16", 17", 18", 19", (35cm, 38cm, 40cm, 43cm, 45cm, 48cm)
Wheels	5" (12cm) Aluminum Front Caster Wheels 5" (12cm) Total Lock Front Caster Wheels 20" or 22" (35cm or 55cm) Front Mag Wheels 5" (15cm) Total Lock Rear Caster Wheels 20" or 22" (35cm or 55cm) Rear Mag Wheels
Lifetime Frame Warranty	10 years on 16 Gauge tubular steel frame
Standard Padding	Upper & Lower Side Pads, Back Pad with Cover

FRAME OPTIONS



CUSHION OPTIONS

Available with solid seat pan option.

matrix®
SEATING SERIES



Libra: Contoured foam prevents sliding and provides support.



Vi: Foam base provides stability and prevents bottoming out.



PS: Lightweight cushion creates superior comfort for users.