

# BRODA<sup>®</sup>

## Product Lookbook



# BRODA<sup>®</sup>

brodaseating.com  
1.844.552.7632

## ***TABLE OF CONTENTS***

3	<i>TYPES OF CHAIRS</i>
6	<i>SEATING &amp; POSITIONING ESSENTIALS</i>
7	<i>POSITIONING</i>
12	<i>PEDAL</i>
15	<i>REHAB (E1161)</i>
20	<i>NON-EMERGENCY MEDICAL TRANSPORT</i>
22	<i>HYGIENE</i>
25	<i>BEHAVIORAL CARE (GLIDERS)</i>
27	<i>PRODUCT PACKAGES &amp; ACCESSORIES</i>
30	<i>CHAIR SPECIFICATIONS</i>
33	<i>UPHOLSTERY &amp; FRAME COLORS</i>
34	<i>CONFIGURE YOUR WHEELCHAIR</i>

***Essential  
Comfort***



# PRODUCT CATEGORIES

Broda's full line of tilt positioning and pedal wheelchairs have been improving lives since 1981. Our chairs provide long-term and complex care with the highest level of comfort, resulting in the highest quality of life.

## POSITIONING



Tilting and reclining functionality provide the highest level of comfort and care

## PEDAL



Wheelchairs that enhance an individual's ability to safely move and be independent

## REHAB



Durable, versatile wheelchairs give support in rehab and rehabilitation therapy

## HYGIENE



Shower and commode chairs facilitate a dignified, multi-positional hygiene experience

## NON-EMERGENCY MEDICAL TRANSPORT (NEMT)



All the benefits of a standard Broda wheelchair, certified for vehicular travel


## BEHAVIORAL CARE (GLIDERS)



Stationary chairs with safe gliding that promote calming and contentment

# Make Sure

## PATIENT BENEFITS



## THE BRODA DIFFERENCE

### Comfort for Care



Comfort is essential, especially for those who spend much of their time in a wheelchair. Make comfortable seating part of long-term care with Broda's **ready out-of-the-box** wheelchairs.


### Proven for Providers

Broda's wheelchairs are meticulously designed to alleviate stress and worry for caregivers and healthcare providers. With easy adjustments and customizable options, Broda products offer **best-in-class** functionality.



# It's a Broda.

## CAREGIVER & FACILITY BENEFITS



## Durable by Design



With four decades of experience (est. 1981) and an **industry-leading**, limited lifetime frame warranty on 16-gauge steel, Broda provides safe, comfortable seating that stands the test of time.

Broda's wheelchairs are trusted by patients, families, and caregivers to improve each patient's physical and emotional wellbeing, leading to **better quality of life**.

## Time-tested for Trust



# SEATING & POSITIONING ESSENTIALS

Features and functionality for most of Broda's positioning, pedal, and rehab wheelchairs

## SEAT TILT

### Keeps individual's hips at a fixed angle

- Facilitates pressure redistribution
- Allows anterior and posterior positioning for comfort
- Angle-adjustable seat accommodates lifts and transfers



## BACK RECLINE

### Opens angle of individual's hips

- Accommodates hip flexion limitations
- Decreases postural deviations, such as lateral lean and slumping
- Trendelenburg-capable for maximum comfort

## ADJUSTABLE & REMOVABLE ARM SUPPORTS

### Provides support for individual's arms

- Improves postural stability with support of upper extremities
- Creates easy patient access for frequent repositioning
- Gives additional support for conditions such as edema



## ADJUSTABLE & REMOVABLE LEG SUPPORTS


### Creates support for individual's legs and/or feet

- Improves overall comfort with support of lower extremities
- Provides lower extremity management and improved circulation
- Facilitates safe front standing transfers

## ADJUSTABLE & REMOVABLE LATERAL SHOULDER SUPPORTS

### Helps maintain individual's postural support

- Promotes a midline posture and minimizes postural deviations
- Easily removable for safe transfers
- Enhances users' comfort and safety



# POSITIONING



The **Synthesis Positioning Wheelchair** has been a blessing for my father, his caregivers and my family. The Synthesis Positioning Wheelchair **improves his quality of life** and I can't imagine him not being in it over the past few years. We are fortunate that the Veterans Association will fund Broda's chairs for veterans like him who need them. - Kent P.

Family Member

# Synthesis

## Positioning Wheelchair

**A FULL RANGE OF POSITIONING FEATURES  
WITH CONTEMPORARY DESIGN**



*Shown with front aluminum caster wheel option*

### **PRODUCT FEATURES**

DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 90° & Trendelenburg Capable
SEAT WIDTHS	16", 18", 20", 22", 24", 26"
SEAT DEPTHS	17", 18.5", 20"
WEIGHT CAPACITY	Up to 350 lbs.



*Tilt functionality provides additional pressure offloading*



*Front pivot seat tilt maintains a forward line of sight, enhancing engagement*



*Adjustable lateral shoulder supports fit a variety of users*



*Easy side transfers for safe patient handling*

# Elite

## Positioning Wheelchair & Tilt Wheelchair

### ENHANCED POSITIONING AND FUNCTIONALITY



Shown with mid-mounted mag wheel option



Tilt and recline functions enhance pressure redistribution and postural support



Accommodates **safe patient handling** and changing of incontinence products



Removable arm supports make transfers safer for caregivers and patients



WC-19 Transport Package available for safe patient vehicle transport

### PRODUCT FEATURES

ONLY AVAILABLE IN CANADA

DEGREES OF TILT  
DEGREES OF RECLINE  
SEAT WIDTHS  
SEAT DEPTHS  
WEIGHT CAPACITY

**Positioning:** Up to 40° | **Tilt:** Up to 40°  
**Positioning:** Up to 90° | **Tilt:** N/A  
16", 18", 20", 22", 24", 26", 28"  
15.5", 17", 18.5", 20"  
Up to 350 lbs.

# Midline

## Positioning Wheelchair

### ENHANCED POSITIONING WITH REFINED FEATURES




*Full tilt and recline provides superior repositioning*



*Adjustable seat surface height adapts the Midline to support users at multiple care levels*



*Tilt functionality enhances comfort*



*Comfort Tension Seating® increases support as well as sitting tolerance*

### PRODUCT FEATURES

DEGREES OF TILT  
DEGREES OF RECLINE  
SEAT WIDTHS  
SEAT DEPTHS  
WEIGHT CAPACITY

Up to 38°  
Up to 90°  
16", 18", 20", 22", 24", 26", 28"  
15.5", 17", 18.5", 20"  
Up to 350 lbs.

# Vanguard

## Bariatric Wheelchair

### BARIATRIC POSITIONING CARE REDEFINED



Shown with upgraded full vinyl wrap option



Anterior and posterior tilt aids in safe patient handling



Swing away, removable arm supports facilitate side-transfers



Manual foot pump allows safe, no-lift tilt positioning



Concave back for gluteal shelf improves user support & dignity

### PRODUCT FEATURES

DEGREES OF TILT  
DEGREES OF RECLINE  
SEAT WIDTHS  
SEAT DEPTHS  
WEIGHT CAPACITY

Up to 30° Posterior / Up to 10° Anterior  
Up to 45°  
24", 26", 28", 30"  
Adjustable 21" or 23"  
Up to 600 lbs.

# PEDAL



Prior to the **Encore wheelchair**, we had tried numerous positioning devices for the resident, which were ineffective. She had a tendency to lean forward and to one side. Once using the Encore, the resident's posture immediately improved. She, herself, was very thankful because she was **"so comfortable"**. The staff couldn't believe the result. Besides this resident sitting in proper alignment, she could still self-propel with her feet. - Linda P.  
Restorative Coordinator



# Encore

## Pedal Wheelchair

**MOBILITY WITH COMFORT,  
FUNCTIONALITY, & STYLE**



**TIME  
B**

**New Product  
Award 2019**

**2019  
TODAY'S  
CAREGIVER  
FRIENDLY  
AWARD  
CAREGIVER.COM**

**mobility**   
product award 2019

*Shown with short back, adjustable head support, and lateral support options*



*Dynamic frame design includes adjustable shoulder supports and tall back option*



*Front pivot seat tilt optimizes hand and foot propulsion*



*Optional dynamic rocking calms and creates contentment*



*Optional WC-19 transportation system for secure vehicle travel*

### PRODUCT FEATURES

DEGREES OF TILT	Up to 25°, Up to 33° with Dynamic Rocking
SEAT WIDTHS	16", 18", 20", 22", 24", 26", 28" 16", 18", 20", 22" with Dynamic Rocking
SEAT DEPTHS	17", 18.5", 20" 15.5", 17", 18.5", 20" with Dynamic Rocking
WEIGHT CAPACITY	350 lbs., 250 lbs. with Dynamic Rocking
SEAT HEIGHTS	12.5", 13.5", 14.5", 15.5", 16.5", 17.5" 14", 15", 16", 17", 18", 19" with Dynamic Rocking

FOR MORE INFORMATION, VISIT [BRODASEATING.COM](http://BRODASEATING.COM)

# Latitude Sashay

**ENHANCED MOBILITY  
AND FREEDOM**

## SASHAY PEDAL WHEELCHAIR



*Shown with swing away foot supports*

## LATITUDE PEDAL WHEELCHAIR



*Shown with swing away, elevating lower leg supports and padded arm supports*

### SASHAY FEATURES

DEGREES OF TILT  
SEAT WIDTHS  
SEAT DEPTHS  
WEIGHT CAPACITY  
SEAT HEIGHTS

Up to 18°  
16", 18", 20", 22", 24", 26", 28"  
15.5", 17", 18.5", 20"  
Up to 350 lbs.  
12.5", 13.5", 14.5", 15.5", 16.5", 17.5"

### LATITUDE FEATURES

Up to 33° with Dynamic Rocking  
16", 18", 20", 22"  
15.5", 17", 18.5", 20"  
Up to 250 lbs.  
14", 15", 16", 17", 18", 19"



*Sashay low seat surface height facilitates foot propulsion*



*Latitude front mount mag wheels enable hand propulsion*



*Latitude gentle rocking motion soothes agitation and creates contentment*

# REHAB

(E1161)



Broda's wheelchairs make it **easier for caregivers** to provide optimal care with less stress on the client and the caregiver. For example, easier transfers, easier re-positioning, less maintenance, and easier mobility. - Carole P. MS/OTR.

# Encore

## Rehab Wheelchair

**MOBILITY WITH COMFORT,  
FUNCTIONALITY, &  
CUSTOMIZED SEATING**

**mobility**   
product award 2020



*Shown with Matrx Libra cushion, Matrx E2 back support, Matrx Loxx head support, and swing-away lower leg supports*

### PRODUCT FEATURES

DEGREES OF TILT	Up to 33° with Dynamic Rocking
SEAT WIDTHS	16", 18", 20", 22"
SEAT DEPTHS	16", 17", 18", 19", 20"
SEAT HEIGHTS	Adjustable 14", 15", 16", 17", 18", 19"
WEIGHT CAPACITY	Up to 250 lbs.



*Back canes fold down for easy storage options*



*Dynamic rocking function soothes and stimulates users*



*Quick-release mag wheels can be removed for easy storage*



*Open seating system allows for almost any seating configuration*

# Synthesis

## Rehab Wheelchair



**Tilt functionality** offloads pressure and promotes upright posture



**Front pivot seat tilt** maintains a forward line of sight, enhancing engagement



**Adjustable upper lateral supports** fit a variety of users



**Solid seat pan** provides support for a variety of cushions

### PRODUCT FEATURES

DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 90°
SEAT WIDTHS	16", 18", 20", 22", 24", 26"
SEAT DEPTHS	17", 18.5", 20"
SEAT HEIGHTS	19.5", 20.5", 22", 23.5"
WEIGHT CAPACITY	Up to 350 lbs.

# Comfort

## Rehab Wheelchair

**SIMPLIFIED DURABILITY FOR REHABILITATION CARE**



*Shown with Matrix Libra cushion, Matrix Elite TR back support, Matrix Elan Head Support, and swing-away lower leg supports*

### PRODUCT FEATURES

DEGREES OF TILT	Up to 45°
SEAT WIDTHS	16", 18", 20", 22"
SEAT DEPTHS	Adjustable 16", 17", 18", 19", 20"
SEAT HEIGHTS	Adjustable 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19"
WEIGHT CAPACITY	Up to 250 lbs.



*Front pivot seat tilt maintains user comfort with smooth, effortless tilt*



*Posterior tilt optimizes pressure redistribution*



*Solid seat pan offers support and is compatible with a variety of seat cushion options*



*Range of simple adjustments accommodates a variety of needs*

# Elite

---

# Latitude

**ENHANCED MOBILITY  
AND FREEDOM**

## ELITE REHAB WHEELCHAIR



*Shown with optional solid seat pan, Matrix Libra cushion, locking twin casters, and pelvic belt.*

## LATITUDE REHAB WHEELCHAIR



*Shown with optional solid seat pan and Matrix PS cushion.*

### ELITE FEATURES

DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 15°
SEAT WIDTHS	16", 18", 20", 22", 24", 26", 28"
SEAT DEPTHS	17", 18.5", 20"
SEAT HEIGHTS	20.5"
WEIGHT CAPACITY	Up to 350 lbs.

### LATITUDE FEATURES

DEGREES OF TILT	Up to 33° with Dynamic Rocking
DEGREES OF RECLINE	N/A
SEAT WIDTHS	16", 18", 20", 22"
SEAT DEPTHS	15.5", 17", 18.5", 20"
SEAT HEIGHTS	Adjustable 14" - 19"
WEIGHT CAPACITY	Up to 250 lbs.



*Elite's recline helps decrease postural deviations*



*Latitude's Dynamic Rocking functionality provides contentment*



*Solid seat pan option accommodates a variety of Matrix cushions for additional support*

# NON-EMERGENCY MEDICAL TRANSPORT



For facilities, we are **using Broda to reduce their cost of patient transportation**. For the consumers it's just comfort. Broda provides 100% comfort; you can't even compare the difference between Broda wheelchairs to other wheelchairs.

- Weslava P.  
Graceful Living Medical Transportation



# Synthesis

## Transport Wheelchair

COMFORT AND SAFETY  
ON THE GO



*Crash tested and WC-19 compliant  
for vehicle transport*



*Easy to load in a variety of  
vehicles*



*Brackets enable safe patient  
transport between locations*



*Tilting functionality allows for  
adjustable comfort in transit*

### PRODUCT FEATURES

DEGREES OF TILT  
DEGREES OF RECLINE  
SEAT WIDTHS  
SEAT DEPTHS  
WEIGHT CAPACITY

Up to 40°  
Up to 90°  
16", 18", 20", 22", 24"  
17", 18.5", 20"  
Up to 350 lbs.

# HYGIENE

A Broda's Revive Shower Commode chair is shown in a bathroom. The chair is white with blue slatted backrest and seat, and black padded armrests and seat. It has four wheels and is positioned in front of a shower stall with a tiled wall and a glass block wall. A showerhead is visible on the wall to the left.

Broda's **Revive Shower Commode** meets our patient's needs. He indicates he feels safe and he is able to participate in his own care each morning. He is **safe at all times** while in the chair, and the caregivers are safe at all times while providing his care. All staff report being pleased with the commode and find it much easier to assist him. - Pam  
Director of Nursing



# Revive

## Shower Commode

**DURABLE BATHING SYSTEM  
FOR SAFE PATIENT CARE**



**Anterior and posterior tilt aid in fall prevention and safe patient handling**



**Multi-positional functionality with tool-free adjustability for ease of use**



**Multiple configurations available to accommodate unique needs**



**Rust-resistant frame and components, with a 10-year warranty**

### **PRODUCT FEATURES**

DEGREES OF TILT	Up to 12° Posterior, Up to 10° Anterior
DEGREES OF RECLINE	Up to 67°
SEAT WIDTH	18", 22"
SEAT DEPTH	18"
WEIGHT CAPACITY	Up to 400 lbs.

# Revive

## Bariatric Shower Commode

**BARIATRIC HYGIENE CARE IN  
A CLASS ALL ITS OWN**



**Optional concave back**  
accommodates and supports  
gluteal shelf



**Mid Mag Wheel** provides  
increased maneuverability



**Anterior tilt** aids in patient  
transfers and ease of patient care



**Rust-resistant frame and components**  
are durable and long-lasting

### PRODUCT FEATURES

DEGREES OF TILT	Up to 14° Posterior, Up to 8° Anterior
DEGREES OF RECLINE	Up to 67°
SEAT WIDTH	27", 33" with width adjustable armrests
SEAT DEPTH	21", 23"
WEIGHT CAPACITY	Up to 800 lbs.

# BEHAVIORAL (GLIDERS)



Our Tranquille gliders are always disappearing out of our common areas and end up in our patients' rooms. We don't know how they get there, but they love them and **prefer sitting in them** more than anything else. One patient with a back injury told us this was the only chair he could sit comfortably in. He wanted to take it home with him after his rehabilitation was complete.

- Josiah D.  
Facility Administrator

# Tranquille (Adult)

## Aspire (Pediatric)

### TRANQUILLE GLIDER



### ASPIRE GLIDER




Gentle *gliding motion*



*Stimulates active patients and reduces the need to wander*



*Auto-locking feature increases safety and prevents falls*



*Comfort Tension Seating® for patients who remain seated throughout the day*

#### TRANQUILLE

#### ASPIRE

SEAT WIDTHS	18", 20"	12", 14", 16"
SEAT DEPTHS	17"	10", 12", 14"
BACK LENGTHS	26", 32"	18", 20", 22"
WEIGHT CAPACITY	Up to 250 lbs.	Up to 200 lbs.

# Padding Packages

## ADDITIONAL POSITIONING PADDING (APP)



- Seat Cushion
- Back and Head Support
- Foot and Calf Support
- Lateral Shoulder Support Cushions

APP shown on Synthesis Positioning Wheelchair

### RECOMMENDED FOR:

Use with positioning chairs to give patients extra comfort.

## SLIM POSTURE PADDING (SPP)



- Seat Cushion
- Back Support
- Extended Back Support\*  
\*for tall back height options

SPP shown on Encore Pedal Wheelchair

### RECOMMENDED FOR:

Use with pedal chairs for extra comfort with a lower profile.

## HUNTINGTON'S SPECIAL PADDING (HSP)



- Long Lateral Shoulder Supports
- Thigh Positioning Strap
- Long Footbox

HSP shown on Synthesis Positioning Wheelchair

### RECOMMENDED FOR:

Users with neurological and movement disorders.

## ACUTE CARE CONFIGURATION (ACC)



- Full Vinyl Covers
- ABS Footrest with Calf and Sole Pads
- Covers Comfort Tension Seating® straps for easy cleaning

ACC shown on Synthesis Positioning Wheelchair

### RECOMMENDED FOR:

Healthcare environments where infection control and durability are a high priority.

# Rehab Seating

## **matrx**<sup>®</sup> SEATING SERIES

For convenience, all Matrx<sup>®</sup> options ship on-chair for Broda Rehab Wheelchairs

**Libra Cushion**



**Vi Cushion**



**PSP Cushion**



**PS Cushion**



**Elite E2 Back**



**Elite Backs (TR/Deep)**



**PB Back (Deep)**



**Loxx Headrest**



**Elan Headrest**



The Matrx<sup>®</sup> Seating series pairs perfectly with Broda's Rehab wheelchairs. Additional cushion, back, and headrest accessories available.

# Accessories

Call us for additional accessory options.  
**1-844-552-7632**

## FRAME ACCESSORIES

- O<sup>2</sup> Holder
- IV Pole
- Rear, Front, or Lateral Stabilizers
- Lateral Torso Supports



O<sup>2</sup> Holder



Lateral Support

## PADDING AND ARM SUPPORTS

- Long Armrests
- Side Arm Support
- Side Reduction Pads
- Arm Trough



Arm Trough



Side Arm Support

## TRAYS

- Slide-On Full Tray
- Half Tray
- Padded Tray Cover



Half Tray



Full Tray

## LEG & FOOT SUPPORTS

- ABS Legrest
- Foot Box
- Swing-Away Footrest
- Legrest Panel
- Leg Support



ABS Legrest



Foot Box

## HEAD & NECK SUPPORTS

- Contoured Neckrest Pad
- Neck Support
- Neckrest with Two Round Bolsters



Neckrest with Bolsters



Contoured Neckrest

## POSITIONING STRAPS

- Thigh Positioning Strap
- Pelvic Positioning Strap



Padded Thigh Strap



Pelvic Positioning Strap

# Chair Specifications

Available Only in Canada



Synthesis Positioning Wheelchair  
Synthesis Transport Wheelchair



Elite Positioning Wheelchair  
Elite Tilt Wheelchair



Midline Positioning Wheelchair



Centric Positioning Wheelchair



Vanguard® Bariatric Wheelchair

Weight Capacity	Up to 350 lbs.	Up to 350 lbs.	Up to 350 lbs.	Up to 250 lbs.	Up to 600 lbs.
Degrees of Tilt	Up to 40°	Up to 40°	Up to 38°	Up to 25°	Up to 30° Posterior Up to 10° Anterior
Degrees of Recline	Up to 90° - Trendelenburg capable	<b>Positioning:</b> Up to 90° <b>Tilt:</b> N/A	Up to 90° - Trendelenburg capable	Up to 45°	Up to 45°
Seat Widths Available	16", 18", 20", 22", 24", 26"	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22", 24", 26", 28"	24", 26", 28", 30"
Seat Depths	<b>Positioning:</b> 17", 18.5", 20", 22" <b>Transport:</b> 17", 18.5", 20"	15.5", 17", 18.5", 20", 22"	15.5", 17", 18.5", 20"	15.5", 17", 18.5", 20"	Adjustable 21" or 23"
Back Heights	32", 33.5", 35"	32", 33.5", 35"	32", 33.5", 35"	32", 33.5", 35"	Adjustable 32", 34", 36", 38", 40"
Seat Surface Heights	19.5", 20.5", 22", 23.5"	19.5", 20.5", 22", 23.5"	Adjustable 17", 18", 19", 20", 21", 22"	Adjustable 19.5", 20.5", 21.5"	Adjustable 20", 21", 22", 23", 24"
Mag Wheels	14" Rear Mounted 14" Mid Mounted 20" Rear Mounted	14" Mid Mounted 14" Rear Mounted 20" Rear Mounted	20", 22", 24" Front Mounted 20", 22" Rear Mounted	20", 22", 24" Front Mounted 20", 22", 24" Rear Mounted	N/A
Recommended Seating	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating®
Recommended Back Support	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating®
Recommended Head Support	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support

# Chair Specifications



*Encore Pedal Wheelchair*



*Sashay Pedal Wheelchair*



*Latitude Pedal Wheelchair*



*Revive Shower Commode  
Revive Bariatric Shower Commode*



*Tranquille Glider  
Aspire Pediatric Glider*

Weight Capacity	Up to 350 lbs. Up to 250 lbs. (with rocking)	Up to 350 lbs.	Up to 250 lbs.	Up to 400 lbs. <b>Bari:</b> up to 800 lbs.	Up to 250 lb. <b>Aspire:</b> up to 200 lb.
Degrees of Tilt	Up to 25° Up to 33° (with dynamic rocking)	Up to 18°	Up to 33°	Up to 12° Posterior Up to 10° Anterior <b>Bari:</b> Up to 14° Posterior Up to 8° Anterior	N/A
Degrees of Recline	N/A	N/A	N/A	Up to 67°	N/A
Seat Widths	16", 18", 20", 22", 24", 26" 28" 16", 18", 20", 22" (with dynamic rocking)	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22"	18", 22" <b>Bari:</b> 27"-33"	18", 20" <b>Aspire:</b> 12", 14", 16"
Seat Depths	17", 18.5", 20"	15.5", 17", 18.5" or 20"	15.5", 17", 18.5" or 20"	18" <b>Bari:</b> 21", 23"	17" <b>Aspire:</b> 10", 12", 14"
Back Heights	22", 32", 33.5", 35"	32", 33.5", 35"	32", 33.5", 35"	26"	26", 32" <b>Aspire:</b> 18", 20", 22"
Seat Surface Heights	Adjustable 12.5" - 17.5" Adjustable 14" - 19" (with dynamic rocking)	Adjustable 12.5", 13.5", 14.5", 15.5, 16.5", 17.5"	Adjustable 14", 15", 16", 17, 18" 19"	24.5"	16", 18"
Mag Wheels	20" or 22" Front Mounted	20", 22" Front Mounted 20", 22" Rear Mounted	20", 22" Front Mounted 20", 22" Rear Mounted	N/A	N/A
Recommended Seating	Comfort Tension Seating* Slim Positioning Padding (SPP) Cushion	Comfort Tension Seating* Slim Positioning Padding (SPP) Cushion	Comfort Tension Seating* Slim Positioning Padding (SPP) Cushion	N/A	N/A
Recommended Back Support	Comfort Tension Seating* Slim Posture Padding (SPP) Back	Comfort Tension Seating* Slim Posture Padding (SPP) Back	Comfort Tension Seating* Slim Posture Padding (SPP) Back	N/A	N/A.
Recommended Head Support	Matrix® Elan Headrest Broda® Neck Support (tall back)	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support Broda® Head Support	Broda® Neck Support

# Chair Specifications



*Synthesis Rehab Wheelchair*



*Encore Rehab Wheelchair*



*Comfort Rehab Wheelchair*



*Elite Rehab Wheelchair*



*Latitude Rehab Wheelchair*

Weight Capacity	Up to 350 lbs.	Up to 250 lbs.	Up to 250 lbs.	Up to 350 lbs.	Up to 250 lbs.
Degrees of Tilt	Up to 40°	Up to 33° (with dynamic rocking)	Up to 45°	Up to 40°	Up to 33° (with dynamic rocking)
Degrees of Recline	Up to 90°	N/A	N/A	15°	N/A
Seat Widths	16", 18", 20", 22", 24", 26"	16", 18", 20", 22"	16", 18", 20", 22"	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22"
Seat Depths	17", 18.5", 20"	16", 17", 18", 19", 20"	16", 17", 18", 19", 20"	17", 18.5", 20", 22"	15.5", 17", 18.5", 20"
Back Cane Angles	N/A	5°, 10°, 15°, 20°	5°, 10°, 15°, 20°	N/A	N/A
Back Heights	32", 33.5", 35"	16", 18", 21", 24"	16", 18", 21", 24"	32"	32", 33.5", 35"
Seat Surface Heights	19.5", 20.5", 22", 23.5"	Adjustable 14", 15", 16", 17", 18", 19"	Adjustable 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19"	20.5"	Adjustable 14", 15", 16", 17", 18", 19"
Mag Wheels	14", 20" Rear Mounted 14" Mid Mounted	20", 22" Front Mounted	20", 22", 24" Rear Mounted	14", 20", 22" Rear Mounted	20", 22" Front or Rear Mounted
Recommended Seating	Matrx® Libra Cushion	Matrx® Libra Cushion	Matrx® Libra Cushion	Matrx® Libra Cushion	Matrx® Libra Cushion
Recommended Back Support	Matrx® E2 Back Support	Matrx® E2 Back Support	Matrx® E2 Back Support	Comfort Tension Seating® Slim Posture Padding (SPP) Back	Comfort Tension Seating® Slim Posture Padding (SPP) Back
Recommended Head Support	Broda Neck Support	Matrx® Loxx Headrest	Matrx® Loxx Headrest	Broda Neck Support	Broda Neck Support

## UPHOLSTERY OPTIONS



## FRAME FINISH OPTIONS




*Some upholstery and frame configurations not available with every chair.*

Call us for  
configuration questions.  
**1-844-552-7632**

# Configure Your Wheelchair

## Step 1: Select Chair Type



## Step 2a: Select Positioning Wheelchair (Caregiver-Propelled)

Features, Functionality, and Customization Level

Bariatric

Vanguard Bariatric Wheelchair



Best  Basic



Synthesis Positioning Wheelchair  
Synthesis Transport Wheelchair



Elite Positioning Wheelchair  
Elite Tilt Wheelchair



Midline Positioning Wheelchair



Centric Positioning Wheelchair

350 lbs. weight capacity Up to 40° Tilt Up to 90° Recline	350 lbs. weight capacity Up to 40° Tilt Up to 90° Recline	350 lbs. weight capacity Up to 38° Tilt Up to 90° Recline	250 lbs. weight capacity Up to 25° Tilt Up to 45° Recline
---	---	---	---

### Optional Add-ons:

- Additional Positioning Padding (APP)
- Huntington's Special Padding (HSP)
- WC-19 Transport Package
- Acute Care Configuration (ACC)
- Slim Posture Padding (SPP)

## Step 2b: Select Pedal Wheelchair (Self-Propelled [Foot or Hand])

Features, Functionality, and Customization Level

Best  Basic



Encore Pedal Wheelchair with  
Dynamic Rocking feature



Encore Pedal Wheelchair



Latitude Pedal Wheelchair with  
Dynamic Rocking Feature



Sashay Pedal Wheelchair

250 lbs. weight capacity Up to 33° Tilt	350 lbs. weight capacity Up to 25° Tilt	250 lbs. weight capacity Up to 33° Tilt	350 lbs. weight capacity Up to 18° Tilt
--	--	--	--

### Optional Add-ons:

- Slim Posture Padding (SPP)
- Additional Positioning Padding (APP)
- WC-19 Transport Package

## Step 3: Select Accessories, Colors, & Customizations

See page 29 for common accessories.

See page 33 for frame and upholstery color options.

Call us for ordering options.  
**1-844-552-7632**

# Configure Your Rehab Chair

## Step 2c: Select Chair Frame

**BRODA®**

Independent

Dependent



*Encore Rehab Wheelchair*



*Comfort Rehab Wheelchair*



*Latitude Rehab Wheelchair*



*Elite Rehab Wheelchair*



*Rehab Wheelchair*

250 lbs. weight capacity Up to 33° Tilt Dynamic Rocking HCPCS E1161	250 lbs. weight capacity Up to 45° Tilt HCPCS E1161	250 lbs. weight capacity Up to 33° Tilt Dynamic Rocking HCPCS E1161	350 lbs. weight capacity Up to 40° Tilt Up to 15° Tilt HCPCS E1161	350 lbs. weight capacity Up to 40° Tilt Up to 90° Tilt HCPCS E1161
--	---	--	---	---

For convenience, all Matrix® options ship on chair.

## Step 2c.ii: Select Seating

**matrix®**  
SEATING SERIES

More Pressure Reduction

Less Pressure Reduction



*Libra*

HCPCS E2624/E2625



*Vi*

HCPCS E2607/E2608



*PSP*

HCPCS E2607



*PS*

HCPCS E2605/E2606

Encore Rehab or Comfort Rehab

Latitude Rehab or Elite Rehab

## Step 2c.iii: Select Back

**matrix®**  
SEATING SERIES

More Supportive

Less Supportive



*Elite E2™*

HCPCS E2620/E2621



*Elite Backs (TR/Deep)*

HCPCS E2620/E2621



*PB & PB Deep*

HCPCS E2615

**BRODA®**



*Comfort Tension Seating®  
Additional Positioning Padding (APP)  
Slim Posture Padding (SPP)*

HCPCS E2601/E2611

Encore Rehab or Comfort Rehab

## Step 2c.iv: Select Headrest

**matrix®**  
SEATING SERIES

More Supportive

Less Supportive



*Loxx Headrest*

HCPCS E0955/E1028



*Elan Headrest*

HCPCS E0955/E1028

## Step 3: Select Additional Accessories

Call us for more options!  
**1-844-552-7632**

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.



# BRODA<sup>®</sup>

*Essential Comfort*

[brodaseating.com](http://brodaseating.com)

**Broda USA**

15455 Conway Rd, Suite 200  
Chesterfield, MO 63017, USA

**Broda Canada**

560 Bingemans Centre Dr  
Kitchener, ON N2B 3X9, Canada