



*Anterior and posterior tilt aid in safe patient handling*



*Tilting and reclining permit safe and thorough bathing*



*Length-adjustable lower leg support and depth-adjustable head support provide enhanced comfort*



*Rust-resistant frame and components with a 10-year warranty*

### **SPECIFICATIONS**

<i>Degrees of Tilt</i>	<i>Up to 12° posterior tilt Up to 10° anterior tilt</i>
<i>Degrees of Recline</i>	<i>Up to 67°</i>
<i>Seat Widths</i>	<i>18", 22" (45mm, 55mm)</i>
<i>Seat Depths</i>	<i>18" (45mm)</i>
<i>Weight Capacity</i>	<i>400 lbs. (181kg)</i>
<i>Lifetime Warranty Frame</i>	<i>10 years on 16-gauge tubular stainless steel white frame</i>

# BRODA<sup>®</sup>

# Revive

## Bariatric Shower Commode



*Concave back option accommodates and supports gluteal shelf*



*Dynamic tilting and reclining permits safe and thorough bathing*



*Anterior and posterior tilt allows for easy user positioning*



*Rust-resistant frame and components are durable and long-lasting*



### **SPECIFICATIONS**

<i>Degrees of Tilt</i>	<i>Up to 14° posterior tilt Up to 8° anterior tilt</i>
<i>Degrees of Recline</i>	<i>Up to 65°</i>
<i>Seat Widths</i>	<i>27", 33" (68mm, 83mm)</i>
<i>Seat Depths</i>	<i>21", 23" (53mm, 58mm)</i>
<i>Weight Capacity</i>	<i>800 lbs. (362kg)</i>
<i>Lifetime Warranty Frame</i>	<i>10 years on 16-gauge tubular stainless steel white frame</i>

**brodaseating.com | 1-844-552-7632**

V2.0 12/2023