

# Sashay

*Pedal Wheelchair*

PEDAL  
**BRODA**®



*Shown with optional front and back caster wheels*



*Front pivot seat tilt maintains user comfort with smooth, effortless tilt*




*Adjustable arm-rests aid in user comfort and positioning*



*Adjustable seat surface height allows for proper foot propulsion*

## **Patient, Facility, & Caregiver Benefits**





- Tilt enables ease of use for frequent repositioning
- Design for safety and stability helping to reduce risk of injury
- Provides stimulation and contentment - excellent for memory care



**Comfort Tension Seating®**  
*improves postural support and sitting tolerance*

## Enhanced Mobility and Freedom


### UPHOLSTERY OPTIONS



### SPECIFICATIONS

Degrees of Tilt	up to 18°
Weight Capacity	350 lbs. (158 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 63cm, 66cm, )
Seat Depths	15.5", 17", 18.5", 20" (39cm, 43cm, 46cm, 50cm)
Seat Surface Heights	Adjustable: 12.5", 13.5", 14.5", 15.5", 16.5", 17.5" (31cm, 34cm, 36cm, 39cm, 41cm, 44cm)
Wheels	5" (12cm) Aluminum Non-locking Front Caster Wheels 5" (12cm) Fixed Front Caster Wheels 20" or 22" (50cm or 55cm) Front or Rear Mag Wheels 5" (12cm) Total Lock Rear Caster Wheels
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Back Support	Integrated with Comfort Tension Seating®
Standard Padding	Head Support, Adjustable Lateral Shoulder Supports, Lower Side Cushions, Seat & Back Cushions with Covers

### FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning belts