

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
<ul style="list-style-type: none"> • 2 x 24" Rear Mag Wheels • 2 x 5" Total Lock Front Stainless Caster Wheels • Adjustable Seat Surface Height • Adjustable Seat Depth • Adjustable Rear Mag Position (14"-21") • Adjustable Seat Angle (0°-6°) • Tension Adjustable Back Straps 	<ul style="list-style-type: none"> • Swing Away Footrests w/ Heel Loops (14"-18" knee-to-heel range) • 10° Back Angle (adjustable to 5° and 15°) • Swing Up Arm Supports • Stainless Steel Adjustable Anti-tippers • 300lb Maximum Occupant Weight • Stainless Steel, Powder Coated Frame

Sentinel Shower Commode - Model: CS313 MSRP \$2200.00			
Seat and Back Support Options			
Seat Width:	<input type="radio"/> 18" Standard <input type="radio"/> 20" \$400 <input type="radio"/> 22" \$500	Seat Depth:	<input type="radio"/> 17" Standard <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 19" \$100
Seat Height:	<input type="radio"/> 19.5" Standard <input type="radio"/> 23.5" N/C	Back Height:	<input type="radio"/> 18" Standard <input type="radio"/> 16" <input type="radio"/> 20" N/C
<input type="radio"/> Parking Brake Extension Handles		\$78.51	

Front Options	Wheel Options	Rear Options
<input checked="" type="radio"/> 5" Front Total Lock Stainless Caster Wheels Standard	<input type="radio"/> 24" Mag Wheels- Pair Standard 2 rear caster options below: (\$400) <input type="radio"/> 5" Total Lock Stainless Caster Wheels (2) N/C <input type="radio"/> 5" Directional Lock Stainless Steel Caster Wheel (1) & 5" Total Lock Stainless Caster Wheel (1)	

Seat Options	Commode Options
<input checked="" type="radio"/> Two in One Commode/Shower Seat Standard	<input type="radio"/> Commode pan with Basket \$62

Common Accessories	
Lateral Supports:	<input type="radio"/> Molded Laterals - Pair (pre-cut slots in back included) \$98.13
Pelvic Belt:	<input type="radio"/> Pelvic Belt (18", 20" seat widths= small, 22"= medium) \$63.79 <input type="radio"/> Center Buckle <input type="radio"/> Side Buckle

Additional Options:

Delivery Instructions or Special Requests:

Total: \$ _____

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.