

BRODA[®]

Sentinel

Shower Commode




Shown with standard mag wheel configuration and optional commode pan and basket



Features a range of simple adjustments for seat depth, wheel position, seat height, and more.



2-in-1 solid seat and commode is highly-versatile and easily mountable in all four directions



Configurable with either user-propelled mag wheels or cost-effective casters

Patient, Facility, & Caregiver Benefits

- 10-year warranty on a highly-durable, rust-resistant, stainless steel frame
- Safe, dignified hygiene with anti-microbial, tension-adjustable upholstery
- Robust, stable frame with stainless steel anti-tip devices for a secure bathing and hygiene experience



Fold-up arm supports and swing-away leg supports allow for easy bathing and transfers

brodaseating.com | 1-844-552-7632

SAFE, COMFORTABLE, INDEPENDENT BATHING

OPTIONS & ACCESSORIES



Molded Lateral Supports



Commode pan with Basket



Pelvic Positioning Strap

SPECIFICATIONS

Seat Width	18", 20", 22" (46cm, 51cm, 56cm)
Seat Depth	16" - 19" (Adjustable) (40cm - 48cm)
Seat Height	19.5" & 23.5" (50cm, & 60cm)
Adjustable Seat Angle	0° - 6° options
Back Height	16", 18", 20" (Adjustable) (41cm, 46cm, 51cm)
Angle-Adjustable Back	5°, 10°, 15° (Adjustable)
Wheel Options	24" (61cm) Rear Mount Mag Wheels 5" (12cm) Total Lock Front Caster Wheels 5" (12cm) Total Lock Rear Caster Wheels 5" (12cm) Directional lock rear
Weight Capacity	300 lbs. (136 kg)
Lifetime Warranty Frame	10-year on 16 gauge tubular stainless steel frame

COMPANION WHEELCHAIR FOR DAILY USE



Encore Rehab Wheelchair

The Encore Rehab Wheelchair offers a universal seating system, customizable for each user's seating and positioning needs. With 33 degrees of infinitely-adjustable tilt and an optional dynamic rocking feature, the Encore Rehab is perfect for a wide variety of care environments.