

PO: _____ **Date:** _____ **Customer #:** _____

Customer Name: _____
Address: _____
Phone: _____
Contact: _____
Email: _____

Ship To Name: _____
Address: _____
Phone: _____
Contact: _____
Tag: _____



Standard Features
<ul style="list-style-type: none"> • Comfort Tension Seating® • Removable Height and Width Adjustable Upper Laterals (Wings) • Up to 40° Infinitely Adjustable Seat Tilt • Up to 90° Infinitely Adjustable Back Recline • Trendelenburg Capable • Height Adjustable, Swing Away Removable Arms • Elevating and Articulating Legrest • 2 x 5" Aluminum Non-Locking Front Casters • 2 x 20" Rear Mag Wheels • 350lb Weight Capacity

Synthesis Transport Wheelchair - Model: V4 WC19 MSRP \$5,308.54

(WC19 Compliant Transit Tie Down Vehicle Transport Package - Factory Installed Only)

Seat and Back Options

Seat Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C	Seat Depth: <input type="radio"/> 17" Standard
<input type="radio"/> 22" \$457.10	<input type="radio"/> 18.5" N/C
<input type="radio"/> 24" \$1,035.70	<input type="radio"/> 20" \$109.94
Seat Height: <input checked="" type="radio"/> 20.5" Standard	Back Height: <input type="radio"/> 32" Standard
	<input type="radio"/> 33.5" \$109.94
	<input type="radio"/> 35" \$109.94

Front Options	Wheel Options	Rear Options
<input type="radio"/> 5" Aluminum Non-Locking Front Casters Standard <i>(available with rear mags only)</i>	<input type="radio"/> 20" Rear Mags Standard	<input type="radio"/> Omit Hand Rims on 20" Rear Mags N/C
<input type="radio"/> 5" Front Locking Twin Wheel Casters N/C <i>(required if rear casters or mid mags selected)</i>	<input type="radio"/> 14" Rear Mags <input type="radio"/> 14" Mid Mags N/C	<input type="radio"/> 6" Multi-Directional Locking Rear Casters #####
	<input type="radio"/> Wheel Lock Extension Handles (pair) \$74.06	

Elevating Footrest Options (Select One)

<input type="radio"/> Adjustable Elevating Split Footrest <i>(not available in 16" seat width)</i> Standard
Choose Number of Straps (footrest length range) <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Adjustable Elevating Single Footrest N/C
Choose Number of Straps (footrest length range) <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")

Armrest Options

Armrest Options: <input type="radio"/> Soft Armrest Standard	Armrest Style: <input checked="" type="radio"/> Height Adjustable Arms Standard
<input type="radio"/> Long Armrest \$44.08	Armrest Mount: <input checked="" type="radio"/> Seat Mounted Arms Standard
<input type="radio"/> Padded Armrests*(pair)(black) \$82.66	
<i>*Not compatible with a tray</i>	

Base Frame Finish Color (Select one)

<input type="radio"/> Carmine Red	<input type="radio"/> Ocean Blue	<input type="radio"/> Midnight Black	N/C
<input type="radio"/> Satin Nickel	<input type="radio"/> Tungsten Sand	<input type="radio"/> Bronze	N/C

Upholstery Color Options: (Select one)

* Standard padding configuration: Headrest, exterior of shoulder bolsters, upper and lower side pads, sole pads come in Accent Color. Seat, back and calf pads, and interior of shoulder bolsters come standard in Savoy Black Vinyl.

** Indigo, Espresso & Graphite are available on seat, back & calf pad & most accessories. See below under padding for selection.

*** In honor of our Veterans Service, Camo options are now available.



Standard Padding*: Shown with Graphite Upholstery & Carmine Red Frame Color *Shown with Optional Magwheel

Accent Padding**: Shown with Indigo Upholstery & Black Frame Color

Upholstery Accent Color*: (Select One)

- Vinyl:** Interlace Indigo** Interlace Graphite**
 Kaleidoscope Smokehouse Kaleidoscope Chocolate
 Savoy Black Interlace Espresso**

- Polyurethane:** Portobello Geyer Fossil Forge
Camo Nylon:* (with performance finish)**
 Digital Camo Green Retro Camo Jungle \$150.43

Padding - Choose Option 1, 2, OR 3

- Option 1 - BRODA Padding** N/C

Standard BRODA Padding:

Includes: **In Accent Color:** Headrest, Exterior Shoulder Bolsters, Upper & Lower Side Pads & Sole Pads.*

In Black Vinyl: Seat & Back Pads with Covers, Calf Pad & Interior of Shoulder Bolsters.

Seat, Back & Calf	<input type="radio"/> Black Vinyl	Standard
Pad Cover Options:	<input type="radio"/> Accent Color*	N/C
	*** <i>(For Indigo, Espresso or Graphite only)</i>	
	<input type="radio"/> Full vinyl Wrap Seat & Back Pads <i>(Select Black Vinyl OR Accent Color Vinyl Above)</i>	\$109.94
	<input type="radio"/> Calf Pad Color <i>(Select Black Vinyl OR Accent Color Vinyl Above)</i>	N/C
	<input type="radio"/> Grey Terry <input type="radio"/> White Terry	N/C
HSP Padding:	<input type="radio"/> Includes Shoulder Bolsters Long, Thigh Belt & Footbox Long	\$595.96
	T-Belt Color: <input type="radio"/> Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	Footbox Long Color: <input type="radio"/> Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	<input type="radio"/> Optional Elevating Fixed Footrest	N/C

- Option 2 - Broda APP (Additional Positioning Padding)** (Not compatible with 35" back height)

Broda APP Padding Includes: Seat Cushion, Back and Head Support, Shoulder Bolsters, Foot & Calf Support, Upper & Lower Side Pads

Seat Width:

- 16" 18" 20"
 22" 24"

Seat Depth:

- 17" 19" 21" \$604.07
 17" 19" 21" \$726.73

Option 3 - Broda SPP (Slim Posture Padding)

Broda SPP Padding Options: Seat, Back, and Back Extension Pads

Seat Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"		Seat Depth:	<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$199.98
	<input type="radio"/> 22"	<input type="radio"/> 24"				<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$220.21
Back:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$159.54	Back Extension:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$117.97
13" Height	<input type="radio"/> 22"	<input type="radio"/> 24"		\$180.88	12" Height	<input type="radio"/> 22"	<input type="radio"/> 24"		\$128.08
Other									
	<input type="radio"/> ACCESSORY-BAG Accessory Bag (Black) Each								\$46.28
	<input type="radio"/> BAG PERSONALIZING (2 Initials)		<input type="radio"/> Brush Font	OR	<input type="radio"/> Block Font			\$28.93	
	<input type="radio"/> CUP-HOLDER Cup Holder								\$52.07

Additional Options:

Delivery Instructions or Special Requests:

Total: \$