

PO: _____

Date: _____

Customer #: _____

Customer Name: _____
 Address: _____
 Phone: _____
 Contact: _____
 Email: _____

Ship To Name: _____
 Address: _____
 Phone: _____
 Contact: _____
 Tag: _____



Standard Features	
• Comfort Tension Seating®	• Lumbar Support
• 16 Gauge Powder Coated Tubular Steel	• Auto-Locking Mechanism
• Height Adjustable Armrests	• 2 Rear Transport Wheels
• Soft Vinyl Covered Seat & Back Pads	• 250lb Weight Capacity
• Caregiver lock	• Headrest on Adult High - Back

Tranquille Glider - Model: 100 MSRP \$1,206.80

Seat and Back Options

Adult Small, Low - Back - Model: 100-10AL

18" Seat Width 17" Seat Depth
 26" Back Height 16" Seat Height

Adult High - Back - Model: 100-20AL

20" Seat Width 17" Seat Depth
 32" Back Height 18" Seat Height

Adult Low - Back - Model: 100-15AL

20" Seat Width 17" Seat Depth
 26" Back Height 18" Seat Height

Chair Color

Padding Color: Black Charcoal Smokehouse Burgundy Graphite Sapphire Blue Meteor Mocha Acid Green Indigo Geyser Carbon Espresso Chocolate Portobello

Frame Color: Satin Nickel Bronze Carmine Red Tungsten Sand Grey Ocean Blue

Armrest Options

Armrest Options: Soft Armrests Standard Arm Trough (each) \$143.27
 Padded Armrests* (pair) \$82.66
 *Not compatible with a tray

Armrest Style: Height Adjustable Arms Standard

Common Accessories

Padding: NR-C - Headrest Contoured Pad \$110.21
 NR-CV - Headrest Visco Foam Contoured Pad \$110.21
 NS - Neck Support Soft Visco Cushion \$110.21
 SR - Side Reductions Pads (pair) \$104.70
 CC-LTS - Comfort Company Laterals \$491.54

Additional Options:

Delivery Instructions or Special Requests: _____

Total: \$