

# Tranquille

*Glider*

SENIOR LIVING  
**BRODA®**



*Gentle gliding motion*




*Stimulates active patients and reduces the need to wander*



*Auto-locking feature increases safety and prevents falls*

## **Patient, Facility, & Caregiver Benefits**

- Provides necessary stimulation for memory care patients
- Increases resident safety and contentment
- Promotes proper spinal posture and pelvic positioning



*Comfort Tension Seating® increases support and sitting tolerance*

## UPHOLSTERY OPTIONS



## SPECIFICATIONS

<b>Weight Capacity</b>	250 lbs.* (114 kg)
<b>Wheels</b>	Rear transport wheel, wheels for non-patient transport**
<b>Lifetime Warranty Frame</b>	10 years on 16 Gauge tubular steel frame
<b>Common Accessories</b>	Padded armrests, full or half tray, adjustable armrests, contoured headrest
<b>Standard Padding</b>	Seat and back pad vinyl cover

\*Adjustments available to accommodate larger patients.  
\*\*The Tranquille Auto-Locking Glider is not a transport chair.

Model	Seat width	Seat depth	Seat height	Back height
Small low-back 100-10AL	18" (45cm)	17" (43cm)	16" (40cm)	26" (66cm)
Low-back 100-15AL	20" (50cm)	17" (43cm)	18" (45cm)	26" (66cm)
High-back 100-20AL	20" (50cm)	17" (43cm)	18" (45cm)	32" (81cm)

## FRAME OPTIONS

