

<b>PO:</b> _____	<b>Date:</b> _____	<b>Customer #:</b> _____
<b>Customer Name:</b> _____	<b>Ship To Name:</b> _____	
<b>Address:</b> _____	<b>Address:</b> _____	
<b>Phone:</b> _____	<b>Phone:</b> _____	
<b>Contact:</b> _____	<b>Contact:</b> _____	
<b>Email:</b> _____	<b>Tag:</b> _____	



Standard Features	
<ul style="list-style-type: none"> <li>- Tubular Steel, Powder Coated Frame</li> <li>- 45 Degrees, Infinitely Adjustable Seat Tilt</li> <li>- Fixed Height Arm Supports</li> <li>- Adjustable Seat Depth 16" - 20"</li> <li>- Adjustable Seat Surface Height 14" - 19"</li> <li>- Choice of Wheel, Caster and Fork Sizes</li> <li>- Push Lock Parking Brakes</li> </ul>	<ul style="list-style-type: none"> <li>- Drop In, Swing Away Lower Leg Supports</li> <li>- Composite Foot Support Plates</li> <li>- Cushioned Side Panels</li> <li>- Adjustable Backrest Angle</li> <li>- 250 lbs. Maximum Occupant Weight</li> <li>- Seat, back and head support cushions not included</li> </ul>

**Comfort Rehab Wheelchair (587) - MSRP \$2,888.24**

**Frame Color Options**

**Frame Color:**     Black     Blue     Bronze     Nickel     Orange     Red     Sand     Yellow

Seat and Back Options			
<b>Seat Width:</b>	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	<b>Seat Depth:</b>	<input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20"
	N/C		N/C
	\$206.70		N/C
	\$206.70		N/C
	N/C		N/C

**Seat Surface Height:**     14"     14.5"     15"     15.5"     16"     16.5"     17"     17.5"     18"     19"

N/C  
N/C  
N/C

**Back Cane Height:**     Ultra-Low Double Bend (16") Fixed    Standard

Ultra Low (16") Fold Down    \$233.20

Low (18") Fold Down    \$233.20

Medium (21") Fold Down    \$233.20

Tall (24") Fold Down    \$233.20

Refer to Seat Height & Wheel Combinations Matrix for assistance.

**Back Cane Angle:**     5°     10°     15°     20°

*Front Mount Options*

5" Fixed Caster Wheels - Pair    Standard

6" Fixed Caster Wheels - Pair    N/C

8" Fixed Caster Wheels - Pair    N/C

*Rear Mount Options*

20" Mag Wheels - Pair    Standard

22" Mag Wheels - Pair    N/C

24" Mag Wheels - Pair    N/C

**Caster Forks: Required**

Small - Recommended for Small Mag Wheel (20" or 22") and Lower Seat Surface Heights    N/C

Medium - Recommended for Large Mag Wheel (24") Higher Seat Surface Heights    N/C

**Mag Wheel Position:** Factory Set (Adjustable in the future)

A (Most Forward) - Smaller Turning Radius    N/C

B    N/C

C    N/C

D\* (Most Rearward) - Larger Turning Radius    N/C

\*Do not recommended with 20" Mag Wheels

Refer to Seat Height & Wheel Combinations Matrix for assistance.

Parking Brake Extension Handles - Pair    N/C     **Anti-Tip Devices** - Required Safety Features - Factory Installed    N/C

**Arm Support Options**

Height Adjustable, Removable Arm Supports    N/C

**Lower Leg Support Assembly and Lower Leg Support Options**

**Lower Leg Support Assembly Style:**    *must choose one*

Drop In Straight - Pair    N/C     Drop In Flared - Pair    N/C

Pin and Saddle Straight - Pair    N/C     Pin and Saddle Flared - Pair    N/C

Drop In IKON Style - Pair    N/C

F-SW75ADJ-587V3    Swing Away 75° Lower Leg Supports - Pair    N/C

F-SW60ADJ-587V3    Swing Away 60° Lower Leg Supports - Pair    N/C

F-SW-IKN    N/C

F-SW-IKN-ADJ    Swing Away Lower Leg Supports Black w/ Angle Adjustable Flip-up Foot Supports    N/C

2" Extended Footrest Tubes - *Only compatible on the F-SW-IKN and F-SW-IKN-ADJ*    \$71.45

**Lower Leg Support Accessories:**

Posterior Foot Supports - Pair    \$21.20

Lower Leg Support Panel (for calf support)    \$63.60

Lower Leg Support Foot Board    \$116.60

F-SW-ELV    Swing Away, Elevating Lower Leg Supports Black - Pair    \$344.50

F-SW-ELV-IKN    Swing Away, Elevating Lower Leg Supports Silver - Pair    \$429.30

F-SW-ANGLE    Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair    \$429.30

F-SW-ANGLE    Swing Away, 70 Degree Lower Leg Support w/ Angle Adjustable Footplates Black - Pair

**Common Accessories**

<b>Pelvic Strap</b> <input type="radio"/> Pelvic Positioning Strap Grommet Mount	<input type="radio"/> Cir Buckle <input type="radio"/> Side Buckle    \$55.12
<b>Trays:</b>	<input type="radio"/> Sewn Strap <input type="radio"/> Grey <input type="radio"/> Black    \$361.86
<input type="radio"/> ABS Tray with Strap (16", 18" 20" 22" - Matches Corresponding Seat Width)	<input type="radio"/> Padded ABS Tray Cover ( ABS Tray only)    \$184.01
<input type="radio"/> Pinned ABS Tray with Strap (16", 18", 20" - Matches Corresponding Seat Width)	<input type="radio"/> Sewn Strap    \$374.13
<b>Other:</b>	<input type="radio"/> Front Mounted    OR <input type="radio"/> Rear Mounted    \$116.54
<input type="radio"/> IV Pole and Assembly	<input type="radio"/> Oxygen Holder and Assembly    \$196.27
<input type="radio"/> Oxygen Holder and Assembly	<input type="radio"/> Accessory Bag (Black) Each    \$49.06
<input type="radio"/> Accessory Bag Personalizing 2 Initials	<input type="radio"/> Accessory Bag Personalizing 2 Initials    \$30.67
<input type="radio"/> O2-HOLDER-BAG Oxygen Holder, Soft Carry Bag	<input type="radio"/> Brush Font    OR <input type="radio"/> Block Font    \$104.94
<input type="radio"/> Cup Holder	<input type="radio"/> Cup Holder    \$55.19

**Matrix® Seating Series**



**Matrix® Seat Cushions**

*For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions*

**Seat Cushion Depth:**     16"     18"     20"     22" (Only available w/ PS Cushion)

**Seat Cushion Width:**     16"     18"     20"

**Seat Cushion Width:**     22"

**Libra** - Recommended for highest risk of or existing skin breakdown    \$745

**Vi** - Recommended for moderate to high risk of skin breakdown    \$624

**PSP** - Recommended for moderate to high risk of skin breakdown    \$645

**PS** - Recommended for low to moderate risk of skin breakdown    \$569

**Libra** - Recommended for highest risk of or existing skin breakdown    \$1,025

**Vi** - Recommended for moderate to high risk of skin breakdown    \$681

**PS** - Recommended for low to moderate risk of skin breakdown    \$975

**Libra Cushion Accessories:**

Libra Fluid Supplement    \$120

Libra 1" Hemi Wedge    \$101

Libra 1" Anterior Seat Wedge    \$202

Libra Tapered Wedge Left    \$85

Libra Tapered Wedge Right    \$85

Libra Abductor    \$85

**Vi, PSP & PS Cushion Accessories:**     Pelvic Obliquity Kit    \$85

**Cushion Rigidizer:**

All Cushions	Cushion Rigidizer:			
	16" Width	16" Depth	18" Depth	20" / 22" Depth
16" Width	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
18" Width	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
20" Width	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115

**Matrix® Backs**

*For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions*

<b>Back Height:</b> <input type="radio"/> 16" <input type="radio"/> 18" (excludes PB & PB Deep) <input type="radio"/> 20"			
<b>Back Width:</b> <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"		<b>Back Width:</b> <input type="radio"/> 22"	
<input type="radio"/> E2 - Recommended for users seeking optimal immersion and airflow \$965	<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$1,151		
<input type="radio"/> E2 Deep - Recommended for users seeking optimal immersion w/ contoured support needs \$1,125	<input type="radio"/> Elite - Recommended for more active users \$1,051		
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$866	<input type="radio"/> Elite Deep - Recommended for active users with contoured support \$1,095		
<input type="radio"/> Elite - Recommended for more active users \$766	<b>Matrix Back</b> <input type="radio"/> Air Force N/C		
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs \$918	<input type="radio"/> Military <input type="radio"/> Army N/C		
<input type="radio"/> PB - Recommended for less active users seeking comfort \$734	<input type="radio"/> Service <input type="radio"/> Marine N/C		
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support \$892	<input type="radio"/> Decals: <input type="radio"/> Navy N/C		
<b>E2 Back Accessories:</b>	<b>Privacy Flap Accessories:</b>		
<input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair) \$385	<input type="radio"/> 16" Width Small Privacy Flap \$53		
<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support \$75	<input type="radio"/> 18" Width Medium Privacy Flap \$53		
<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32"L \$295	<input type="radio"/> 20" Width Large Privacy Flap \$53		
<input type="radio"/> Quad Release Tab for PostureFlex \$35	<input type="radio"/> 22" Width Extra Large Privacy Flap \$53		
<b>Back Slipcovers:</b>	<b>Lateral Support Accessories:</b>		
<input type="radio"/> E2 Infection Control Slipcover \$95	<input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L \$116		
<input type="radio"/> E2 Deep Infection Control Slipcover \$95	<input type="radio"/> Elite Back Swing Away Lateral Hardware \$243		
<input type="radio"/> Infection Control Slipcover for Matrix Cushions \$95	<input type="radio"/> Swing Away Lateral Pad (Each) 4.25"Hx6.5"L \$122		

**Matrix® Headrest**

<b>Head Support Mounting:</b>	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support \$285		
<b>Choose One</b>	<input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support \$275		
<b>Head Support Options:</b>	<input type="radio"/> Head Support <input type="radio"/> 10" <input type="radio"/> 14" \$310		
	<input type="radio"/> 4 Point Head Support <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10" \$448		
	<input type="radio"/> Headband with 4 Point Cushion <input type="radio"/> Small <input type="radio"/> Large \$610		

**Additional Options:**

**Delivery Instructions or Special Requests:**

Total: \$

Accessory Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessory services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
  - Residential Delivery Service:** Delivery to a residential location
  - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
  - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
  - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri: 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.*

Knee-to-Heel Range (Comfort Rehab ONLY)			
Front Rigging Model	Description	Minimum	Maximum
Range Adjusts in 3/4" increments			
F-SW 75-587V3	75° Drop-in, black	14.25"	17.25"
F-SW 75ADJ-587V3	75° Drop-in w/ angle adjustable footplates, black	14.25"	17.25"
F-SW 60ADJ-587V3	60° Drop-in w/ angle adjustable footplates, black	15"	18"
F-SW-ELV-587V3	Drop-in, elevating, black	15"	18"
F-SW-ELV-ANGLE-587V3	Drop-in, elevating w/ angle adjustable footplates, black	15"	18"

587V3: SEAT HEIGHT & WHEEL COMBINATIONS									
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB	
14"	20	5	SMALL	1	1	TOP	3	2	
	20	5	SMALL	2	2	TOP	2	1	
15"	22	5	SMALL	1	2	TOP	2	1	
	20	6	MEDIUM	2	1	TOP	2	1	
15.5"	22	6	MEDIUM	1	1	TOP	2	1	
	24	6	SMALL	1	2	TOP	1	1	
	24	5	MEDIUM	1	2	TOP	1	1	
	24	6	MEDIUM	1	1	TOP	1	1	
16"	20	5	SMALL	3	3	TOP	1	1	
	20	5	MEDIUM	3	3	TOP	1	1	
	20	6	MEDIUM	3	2	TOP	1	1	
	22	5	SMALL	2	3	TOP	1	1	
	22	6	SMALL	2	2	TOP	1	1	
	22	5	MEDIUM	2	3	TOP	1	1	
16.5"	22	6	MEDIUM	2	2	TOP	1	1	
	24	5	SMALL	1	3	TOP	1	1	
	20	6	SMALL	3	3	TOP	1	1	
	24	6	SMALL	2	3	BOTTOM	3	1	
17"	24	5	MEDIUM	2	4	BOTTOM	3	1	
	24	6	MEDIUM	2	3	BOTTOM	3	1	
	20	5	MEDIUM	4	5	BOTTOM	3	1	
	20	6	MEDIUM	4	4	BOTTOM	3	1	
17.5"	22	5	MEDIUM	3	5	BOTTOM	3	1	
	22	6	MEDIUM	3	4	BOTTOM	3	1	
	20	8	SMALL	4	3	BOTTOM	3	1	
	22	8	SMALL	3	3	BOTTOM	3	1	
18"	24	8	SMALL	3	3	BOTTOM	2	1	
	24	6	MEDIUM	3	5	BOTTOM	2	1	
	24	8	MEDIUM	3	3	BOTTOM	2	1	
	20	8	MEDIUM	5	4	BOTTOM	2	1	
19"	22	8	MEDIUM	4	4	BOTTOM	2	1	
	24	8	MEDIUM	4	5	BOTTOM	1	1	

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.  
Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels  
Form08CONFIG-404 587V3 Rev. 1