

<b>PO:</b> _____	<b>Date:</b> _____	<b>Customer #:</b> _____
<b>Customer Name:</b> _____	<b>Ship To Name:</b> _____	
<b>Address:</b> _____	<b>Address:</b> _____	
<b>Phone:</b> _____	<b>Phone:</b> _____	
<b>Contact:</b> _____	<b>Contact:</b> _____	
<b>Email:</b> _____	<b>Tag:</b> _____	



Standard Features	
<ul style="list-style-type: none"> <li>· Tubular Steel, Powder Coated Frame</li> <li>· 45 Degrees, Infinitely Adjustable Seat Tilt</li> <li>· Fixed Height Arm Supports</li> <li>· Adjustable Seat Depth 16" - 20"</li> <li>· Adjustable Seat Surface Height 14" - 19"</li> <li>· Choice of Wheel, Caster and Fork Sizes</li> <li>· Push Lock Parking Brakes</li> </ul>	<ul style="list-style-type: none"> <li>· Drop In, Swing Away Lower Leg Supports</li> <li>· Composite Foot Support Plates</li> <li>· Cushioned Side Panels</li> <li>· Adjustable Backrest Angle</li> <li>· 250 lbs. Maximum Occupant Weight</li> <li>· Seat, back and head support cushions not included</li> </ul>


**Comfort Rehab Wheelchair (587) - (E1161) MSRP \$3,045**

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Frame Color Options	
<input type="radio"/> Black <input type="radio"/> Blue <input type="radio"/> Bronze <input type="radio"/> Nickel	<input type="radio"/> Orange <input type="radio"/> Red <input type="radio"/> Sand <input type="radio"/> Yellow

Seat and Back Options	
<b>Seat Width:</b> <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" (E2201)    \$217.40 <input type="radio"/> 22" (E2201)    \$217.40 <input type="radio"/> Solid Seat Pan (E2231) Factory Installed    \$166.95 <input type="radio"/> Omit Seat Pan    N/C	<b>Seat Depth:</b> <input type="radio"/> 16" <input type="radio"/> 17"    N/C <input type="radio"/> 18" <input type="radio"/> 19"    N/C <input type="radio"/> 20" (E2203)    N/C

<b>Seat Surface Height:</b> <input type="radio"/> 14" <input type="radio"/> 14.5" <input type="radio"/> 15" <input type="radio"/> 15.5"    N/C <input type="radio"/> 16" <input type="radio"/> 16.5" <input type="radio"/> 17" <input type="radio"/> 17.5"    N/C <input type="radio"/> 18" <input type="radio"/> 19"    N/C	<b>Back Cane Height:</b> <input type="radio"/> Ultra-Low Double Bend (16") Fixed    Standard <input type="radio"/> Ultra Low (16") Fold Down    \$217.40 <input type="radio"/> Low (18") Fold Down    \$217.40 <input type="radio"/> Medium (21") Fold Down    \$217.40 <input type="radio"/> Tall (24") Fold Down    \$217.40
<p style="text-align: center;"><b>Refer to Seat Surface Height and Wheel Combinations Matrix for assistance.</b></p>	<b>Back Cane Angle:</b> <input type="radio"/> 5° <input type="radio"/> 10° <input type="radio"/> 15° <input type="radio"/> 20°



Front Mount Options	Wheel Options	Rear Mount Options
<input type="radio"/> 5" Fixed Caster Wheels - Pair    Standard <input type="radio"/> 6" Fixed Caster Wheels - Pair    N/C <input type="radio"/> 8" Fixed Caster Wheels - Pair    N/C	<input type="radio"/> 20" Mag Wheels - Pair    Standard <input type="radio"/> 22" Mag Wheels - Pair    N/C <input type="radio"/> 24" Mag Wheels - Pair    N/C	

<b>Caster Wheel Forks: Required</b> <input type="radio"/> Small - Recommended for Small Mag Wheel (20" or 22") and Lower Seat Surface Heights    N/C <input type="radio"/> Medium - Recommended for Large Mag Wheel (24") Higher Seat Surface Heights    N/C	<b>Mag Wheel Position: Factory Set (Adjustable in the future)</b> <input type="radio"/> A (Most Forward) - Smaller Turning Radius    N/C <input type="radio"/> B    N/C <input type="radio"/> C    N/C <input type="radio"/> D* (Most Rearward) - Larger Turning Radius    N/C <small>*Do not recommended with 20" Mag Wheels</small>
--	--

<input type="radio"/> Parking Brake Extension Handles (pair) (E0961)    \$73.46	<input checked="" type="radio"/> Anti-Tip Devices (E0971) - Required Safety Features - Factory Installed    \$150.26
---	--

Arm Support Options
<input type="radio"/> Fixed Height Arm Supports    Standard <input type="radio"/> Height Adjustable, Removable Arm Supports (E0973)    \$278.25

Lower Leg Support Assembly and Lower Leg Support Options	
<b>Lower Leg Support Assembly Style: must choose one</b> <input type="radio"/> Drop In Straight - Pair    N/C <input type="radio"/> Drop In Flared - Pair    N/C <input type="radio"/> Pin and Saddle Straight - Pair    N/C <input type="radio"/> Pin and Saddle Flared - Pair    N/C <input type="radio"/> Drop In IKON Style - Pair    N/C <input type="radio"/> F-SW75ADJ-587V3 Swing Away 75° Lower Leg Supports - Pair    NC <input type="radio"/> F-SW60ADJ-587V3 Swing Away 60° Lower Leg Supports - Pair    NC <input type="radio"/> F-SW-IKN Swing Away Lower Leg Supports Silver - Pair    NC <input type="radio"/> F-SW-IKN-ADJ Swing Away Lower Leg Supports Black    NC <input type="radio"/> Steel w/ Angle Adjustable Flip-up Foot Supports - Pair (K004)    \$111.30 <input type="radio"/> 2" Extended Footrest Tubes - Only compatible on the F-SW-IKN and F-SW-IKN-ADJ    \$66.78	<b>Lower Leg Support Accessories:</b> <input type="radio"/> Posterior Foot Supports - Pair E0951    \$22.26 <input type="radio"/> Lower Leg Support Panel (for calf support) E0995    \$73.46 <input type="radio"/> Lower Leg Support Foot Board    \$132.45 <input type="radio"/> F-SW-ELV Swing Away, Elevating Lower Leg Supports Black - Pair (E0990)    \$296.06 <input type="radio"/> F-SW-ELV-IKN Swing Away, Elevating Lower Leg Supports Silver - Pair (E0990)    \$296.06 <input type="radio"/> Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair (K004)    \$779.10 <input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support W/ Angle Adjustable Footplates, Black - Pair (K0040)    \$111.30

Common Accessories	
<b>Pelvic Belt</b> <input type="radio"/> Pelvic Positioning Strap Grommet Mount (E0978)	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle    \$63.79
<b>Trays:</b> <input type="radio"/> Tray ABS with Strap (E0950) (16", 18" 20" 22" - Matches Corresponding Seat Width)    \$382.71 <input type="radio"/> Tray ABS Cushioned Cover    \$193.82 <input type="radio"/> Clear Tray with Strap (E0950) (16", 18", 20" - Matches Corresponding Seat Width)    \$394.98	<input type="radio"/> Sewn Strap <input type="radio"/> Grey <input type="radio"/> Black <input type="radio"/> Sewn Strap    \$394.98

<b>Other:</b> <input type="radio"/> IV Pole and Assembly - (K0105)    \$122.66 <input type="radio"/> Oxygen Holder and Assembly (E2208)    \$206.03 <input type="radio"/> Accessory Bag (Black) Each    \$49.70 <input type="radio"/> Accessory Bag Personalizing 2 Initials    \$30.68 <input type="radio"/> O2-HOLDER-BAG Oxygen Holder, Soft Carry Bag    \$97.52 <input type="radio"/> Cup Holder    \$55.20	<input type="radio"/> Front Mounted    OR <input type="radio"/> Rear Mounted <input type="radio"/> Brush Font    OR <input type="radio"/> Block Font
--	---

**Matrix® Seating Series**



**Matrix® Seat Cushions**

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth:  16"  18"  20"  22" (Only available w/ PS Cushion)

Seat Cushion Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"
<input type="radio"/> Libra (E2624) - Recommended for highest risk of or existing skin breakdown			\$485
<input type="radio"/> Vi (E2607) - Recommended for moderate to high risk of skin breakdown			\$425
<input type="radio"/> PSP (E2607) - Recommended for moderate to high risk of skin breakdown			\$395
<input type="radio"/> PS (E2605) - Recommended for low to moderate risk of skin breakdown			\$345

Seat Cushion Width:	<input type="radio"/> 22"
<input type="radio"/> Libra (E2625) - Recommended for highest risk of or existing skin breakdown	\$595
<input type="radio"/> Vi (E2608) - Recommended for moderate to high risk of skin breakdown	\$475
<input type="radio"/> PS (E2606) - Recommended for low to moderate risk of skin breakdown	\$545

Libra Cushion Accessories:	<input type="radio"/> Libra Fluid Supplement	\$75
	<input type="radio"/> Libra 1" Hemi Wedge	\$65
	<input type="radio"/> Libra 1" Anterior Seat Wedge	\$125
	<input type="radio"/> Libra Tapered Wedge Left	\$50
	<input type="radio"/> Libra Tapered Wedge Right	\$50
	<input type="radio"/> Libra Abductor	\$60

Vi, PSP & PS Cushion Accessories:  Pelvic Obliquity Kit \$50

Cushion Rigidizer: All Cushions	E0992	16" Depth	18" Depth	20" / 22" Depth
		<input type="radio"/> 16" Width <input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
		<input type="radio"/> 18" Width <input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
		<input type="radio"/> 20" Width <input type="radio"/> \$115	<input type="radio"/> \$115	<input type="radio"/> \$115
		22" Width	N/A	<input type="radio"/> \$115 <input type="radio"/> \$115

**Matrix® Backs**

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height:  16"  18" (excludes PB & PB Deep)  20"

Back Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"
<input type="radio"/> E2 (E2615) - Recommended for users seeking optimal immersion and airflow			\$645
<input type="radio"/> E2 Deep (E2620) - Recommended for users seeking optimal immersion w/ contoured support needs			\$775
<input type="radio"/> Elite TR (E2620) - Recommended for tilt in space & recline applications			\$650
<input type="radio"/> Elite (E2620) - Recommended for more active users			\$625
<input type="radio"/> Elite Deep (E2620) - Recommended for active users with contoured support needs			\$760
<input type="radio"/> PB (E2615) - Recommended for less active users seeking comfort			\$625
<input type="radio"/> PB Deep (E2615) - Recommended for less active users seeking comfort & support			\$645
E2 Back Accessories:	<input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair) (E0956)	\$290	
	<input type="radio"/> E2 Chest Harness Interface	\$75	
	<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support	\$55	
	<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32"L	\$245	
	<input type="radio"/> Quad Release Tab for PostureFlex	\$25	
Back Slipcovers:	<input type="radio"/> E2 Infection Control Slipcover	\$65	
	<input type="radio"/> E2 Deep Infection Control Slipcover	\$65	
	<input type="radio"/> Infection Control Slipcover for Matrix Cushions	\$65	

Back Width:	<input type="radio"/> 22"	
<input type="radio"/> Elite TR (E2621) - Recommended for tilt in space & recline applications	\$845	
<input type="radio"/> Elite (E2621) - Recommended for more active users	\$795	
<input type="radio"/> Elite Deep (E2621) - Recommended for active users with contoured support needs	\$930	
Matrix Back	<input type="radio"/> Air Force	N/C
Military Service Decals:	<input type="radio"/> Army	N/C
	<input type="radio"/> Marine	N/C
	<input type="radio"/> Navy	N/C
Privacy Flap	<input type="radio"/> 16" Width Small Privacy Flap	\$40
Accessories:	<input type="radio"/> 18" Width Medium Privacy Flap	\$40
	<input type="radio"/> 20" Width Large Privacy Flap	\$40
	<input type="radio"/> 22" Width Extra Large Privacy Flap	\$40
Lateral Support	<input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L (E0956)	\$110
Accessories:	<input type="radio"/> Elite Back Swing Away Lateral Hardware (E0128)	\$225
	<input type="radio"/> Swing Away Lateral Pad (Each) 4.25"Hx6.5"L (E0956)	\$110

**Matrix® Headrest**

Head Support Mounting:	<input type="radio"/> Elan Head Support Hardware (E0128) - Required with Selection of Head Support	\$220
Choose One	<input type="radio"/> Loxx Head Support Hardware (E0128) - Required with Selection of Head Support	\$215
Support Options:	<input type="radio"/> Head Support Cushion (E0955)	<input type="radio"/> 10" <input type="radio"/> 14" \$235
	<input type="radio"/> Head Support 4 Point Cushion (E0955)	<input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10" \$325
	<input type="radio"/> Headband with 4 Point Cushion (E0955)	<input type="radio"/> Small <input type="radio"/> Large \$450

**Additional Options:**

**Delivery Instructions or Special Requests:**

Total: \$

Accessory Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessory services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
  - Residential Delivery Service:** Delivery to a residential location
  - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
  - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
  - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.*

587V3: SEAT HEIGHT & WHEEL COMBINATIONS									
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB	
14"	20	5	SMALL	1	1	TOP	3	2	2
	22	5	SMALL	1	2	TOP	2	1	1
15.5"	20	6	MEDIUM	2	1	TOP	2	1	1
	22	6	MEDIUM	1	1	TOP	2	1	1
	24	6	SMALL	1	2	TOP	1	1	1
	24	5	MEDIUM	1	2	TOP	1	1	1
16"	24	6	MEDIUM	1	1	TOP	1	1	1
	20	5	SMALL	3	3	TOP	1	1	1
	20	5	MEDIUM	3	3	TOP	1	1	1
	20	6	MEDIUM	3	2	TOP	1	1	1
16.5"	22	5	SMALL	2	3	TOP	1	1	1
	22	6	SMALL	2	2	TOP	1	1	1
	22	5	MEDIUM	2	3	TOP	1	1	1
	22	6	MEDIUM	2	2	TOP	1	1	1
17"	24	5	SMALL	1	3	TOP	1	1	1
	20	6	SMALL	3	3	TOP	1	1	1
	24	6	SMALL	2	3	BOTTOM	3	1	1
	24	5	MEDIUM	2	4	BOTTOM	3	1	1
17.5"	24	6	MEDIUM	2	3	BOTTOM	3	1	1
	20	5	MEDIUM	4	5	BOTTOM	3	1	1
	20	6	MEDIUM	4	4	BOTTOM	3	1	1
	22	5	MEDIUM	3	5	BOTTOM	3	1	1
18"	22	6	MEDIUM	3	4	BOTTOM	3	1	1
	20	8	SMALL	4	3	BOTTOM	3	1	1
	24	8	SMALL	3	3	BOTTOM	2	1	1
	24	6	MEDIUM	3	5	BOTTOM	2	1	1
19"	24	8	MEDIUM	3	3	BOTTOM	2	1	1
	20	8	MEDIUM	5	4	BOTTOM	2	1	1
	22	8	MEDIUM	4	4	BOTTOM	2	1	1

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.  
Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG- 404 587V3 Rev. 1

Knee-to-Heel Range (Comfort Rehab ONLY)			
Front Rigging Model	Description	Minimum	Maximum
Range Adjusts in 3/4" increments			
F-SW 75-587V3	75° Drop-in, black	14.25"	17.25"
F-SW 75ADJ-587V3	75° Drop-in w/ angle adjustable footplates, black	14.25"	17.25"
F-SW 60ADJ-587V3	60° Drop-in w/ angle adjustable footplates, black	15"	18"
F-SW-ELV-587V3	Drop-in, elevating, black	15"	18"
F-SW-ELV-ANGLE-587V3	Drop-in, elevating w/ angle adjustable footplates, black	15"	18"