

<b>PO:</b> _____	<b>Date:</b> _____	<b>Customer #:</b> _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



(shown with optional upgrades)

Standard Features	
• 250lb Maximum Occupant Weight (rocking)	• Rocking Mechanism (for up to 33° of tilt)
• 350lb Maximum Occupant Weight (non-rocking)	• Seat Mount Lower Leg Support Assembly
• Tubular Steel, Powder Coated Frame	• 2 x 20" Front Mag wheels with Front Anti-Tippers
• Adjustable Seat Surface Height	• 2 x 5" Total Lock Rear Caster Wheels
• Independent Adjustable Arm Support Height	
• Up to 25° Infinitely Adjustable Seat Tilt (non-rocking)	
• Up to 33° Infinitely Adjustable Seat Tilt (rocking)	

### Encore Pedal Wheelchair® - Model: 48V4/48RV4 MSRP \$3,874.95

**WC19 Transit Tie Down Vehicle Package\*** - Factory Installed Only, 300lb Maximum Occupant Weight \$375.13

\*Available up to 24" seat width and front mag wheels only.

Seat and Back Support Options			
<b>Rocking Mechanism:</b>	<input type="radio"/> Eliminate Rocking Function	<b>-\$387.50</b>	
<b>Seat Width:</b>	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	N/C	
		\$512.73	
<b>Non-Rocking Only</b>	<input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28"	\$1,161.62	
<b>Seat Surface Height Rocking</b>	<input type="radio"/> 14" <input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19"	N/C	
		N/C	
<b>Seat Surface Height No Rock</b>	<input type="radio"/> 12.5" <input type="radio"/> 13.5" <input type="radio"/> 14.5" <input type="radio"/> 15.5" <input type="radio"/> 16.5" <input type="radio"/> 17.5"	N/C	
		N/C	
<b>Seat Depth:</b>	<input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"	Standard	
		N/C	
		\$122.66	
<b>Standard Back Support Height:</b>	<input type="radio"/> 32" <input type="radio"/> 33.5" <input type="radio"/> 35"	Standard	
		\$122.66	
		\$122.66	
<b>Short Back Support Height</b>	<input type="radio"/> 22" (does not include Adjustable Lateral Shoulder Supports and Assembly, for Neck Support Option - see page 2)	-\$490.67	

Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair	Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair
		Standard
<input type="radio"/> 22" Mag Wheels - Pair	N/C	<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels*
		\$233.60
		*(available with front mags only)
<input type="radio"/> 5" Total Lock Caster Wheels - Pair		<input type="radio"/> 14" Mag Wheels - Pair
		N/C
<input type="radio"/> 5" Total Lock Caster Wheels - Pair	<b>(\$680.01)</b>	<input type="radio"/> 5" Fixed Non-Locking Caster Wheels - Pair

Arm Support Options			
<b>Arm Support Cushion Options:</b>	<input type="radio"/> Arm Support - Pair	Standard	
	<input type="radio"/> Arm Support Trough - Each	\$166.82	
	<input type="radio"/> Cushioned Arm Support* - Pair	\$98.13	
	*Not compatible with trays		
<b>Arm Support Style:</b>	<input type="radio"/> Height Adjustable - Pair	Standard	
	<input type="radio"/> Height and Depth Adjustable - Pair	\$169.27	
<b>Arm Mount</b>	<input checked="" type="radio"/> Seat Mounted Arm Supports - Pair	Standard	

Seating and Frame Color Options			
<b>Comfort Tension Seating Color:</b>	<input checked="" type="radio"/> Black	Standard	N/C
<b>Frame Colors:</b>	<input type="radio"/> Black	<input type="radio"/> Deep Orange	N/C
	<input type="radio"/> Tungsten Sand	<input type="radio"/> Bronze	<input type="radio"/> Zinc Yellow
	<input type="radio"/> Carmine Red	<input type="radio"/> Ocean Blue	<input type="radio"/> Satin Nickel

Upholstery Style Options <i>(Select one)</i>				
<b>Accent Color Vinyl:</b>	<input type="radio"/> Interlace Indigo***	<input type="radio"/> Interlace Espresso***	<input type="radio"/> Interlace Graphite***	
	<input type="radio"/> Savoy Black			
<b>Vinyl:</b>	<input type="radio"/> Beeline Portobello	<input type="radio"/> Beeline Geyser	<input type="radio"/> Beeline Forge	<input type="radio"/> Kaleidoscope Smokehouse
<b>Nylon:</b>	<input type="radio"/> Digital Camo	<input type="radio"/> Retro Camo Jungle		
				N/C
				\$159.47

\*\*Upholstery Colors come with seat and back cushions in black vinyl. \*\*\*Interlace Indigo, Espresso, and Graphite available with seat and back cushions in matching vinyl color

Cushions - Choose Option 1, 2, OR 3				
<input type="radio"/> <b>Option 1 - Broda Standard Cushions</b>				N/C
Includes: Neck Support, Adjustable Lateral Shoulder Supports, Lower Side Covers, Seat and Back Support Cushions with Covers				
<b>Seat &amp; Back Pad Cover</b>	<input type="radio"/> Black Vinyl			Standard
<b>Options:</b>	<input type="radio"/> Accent Color Vinyl ***Option For Interlace: Indigo, Espresso or Graphite only			N/C
	(Select Black Vinyl OR Accent Color Vinyl Above)			
	<input type="radio"/> Full Vinyl Wrap Seat and Back Support Cushions			\$122.66
	(Select Black Vinyl OR Accent Color Vinyl Above)			
	<input type="radio"/> Grey Terry Cover	<input type="radio"/> White Terry Cover		N/C
<b>Lateral Shoulder Support</b>	<input type="radio"/> Lateral Shoulder Supports - Pair			N/C
<b>Options:</b>	<input type="radio"/> Long Lateral Shoulder Supports - Pair			\$166.81

**Optional Substitutions Available Between APP & SPP Package Cushion Components**

**Option 2 - Broda Additional Positioning Padding (APP), 24"-28" seat widths for non-rocking only**

Broda APP Package Includes: Seat, Back and Head Support Cushions (no Head Support or Lateral Shoulder Supports on short back), plus Adjustable, Removable Lateral Shoulder Supports and Lateral Upper Leg Support Cushions

<b>Seat Width:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	<b>Seat Depth:</b>	<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	
	<input type="radio"/> 22"	<input type="radio"/> 24"			<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$451.41
	<input type="radio"/> 26"	<input type="radio"/> 28"			<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$706.54
					<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$836.57

**Option 3 - Broda Slim Posture Padding (SPP), 24"-26" seat widths for non-rocking only**

Broda SPP Cushion Options: Seat, Back Support, and Back Support Extension Cushions

<b>Seat Width:</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	<b>Seat Depth:</b>	<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	
	<input type="radio"/> 22"	<input type="radio"/> 24"			<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$161.97
	<input type="radio"/> 26"				<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$178.64
					<input type="radio"/> 17"	<input type="radio"/> 18.5"	<input type="radio"/> 20"	\$207.22
<b>Lower Back</b>	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"		<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$142.91
<b>Support</b>	<input type="radio"/> 22"	<input type="radio"/> 24"			<input type="radio"/> 22"	<input type="radio"/> 24"		\$107.19
13" Height	<input type="radio"/> 26"				<input type="radio"/> 22"	<input type="radio"/> 24"		\$171.49
				12" Height	<input type="radio"/> 26"			\$207.22
								\$119.09
								\$131.01

**Matrx® Head Support (Available with 22" Short Back Option Only)**

<b>Head Support Mounting:</b>	<input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support			\$235.40
<b>Choose One</b>	<input type="radio"/> Loxx Head Support Hardware (E0128) - Required with Selection of Head Support			\$215
<b>Support Options:</b>	<input type="radio"/> Head Support Cushion	<input type="radio"/> 10"	<input type="radio"/> 14"	\$251.45
	<input type="radio"/> Head Support 4 Point Cushion	<input type="radio"/> 9" x 7"	<input type="radio"/> 11" x 10"	\$347.75
	<input type="radio"/> Head Support with 4 Point Cushion	<input type="radio"/> Small	<input type="radio"/> Large	\$481.50

<i>Drop In Style</i>	<b>Must Choose One Lower Leg Support Assembly Style</b>	<i>Pin and Saddle Style</i>
<input type="radio"/> Lower Leg Support Assembly (Drop-in)- Pair	N/C	<input type="radio"/> Lower Leg Support Assembly (Pin & Saddle) - Pair*
		N/C

\*Pin & Saddle Lower Leg Support Assembly not available with 16" SW chair due to overlapping foot supports

**Must Choose One Seat Mount Lower Leg Support Assembly Location**

- Seat Frame Mount\* N/C |  Sub-Seat Frame Mount\*\* N/C

\*Not available on rocking version with 17" seat depth

\*\*Improves the ease of the rocking mechanism


\*Not available on non-rocking version when choosing 12.5"

seat surface height with pin & saddle style lower leg support assembly

Seat Frame Mount



Sub-Seat Frame Mount



Drop-In Style	Lower Leg Support Options	Pin & Saddle Style
<input type="radio"/> F-SW-IKN: Swing Away Lower Leg Supports Silver - Pair \$116.54	<input type="radio"/> F-SW: Swing Away Lower Leg Supports Black- Pair \$186.45	
<input type="radio"/> F-SW-ELV-IKN: Elevating Swing Away Lower Leg Supports \$484.53	<input type="radio"/> F-SW-ELV: Elevating Swing Away Lower Leg Supports \$250.	
Silver Aluminum w/ Angle Adjustable Flip-up Foot Supports - Pair	Black w/ Flip-up Foot Supports - Pair	
<input type="radio"/> F-SW-IKN-ADJ: Swing Away Lower Leg Supports Black \$239.37	<input type="radio"/> F-SW-ANGLE Swing Away 70 Degree Lower Leg Support w/ Angle Adjustable Footplates Black \$316.48	
Steel w/ Adjustable Flip-up Foot Supports - Pair	<input type="radio"/> LEG-FOOT-ENCORE*: Lower Leg Support Foot Board \$291.	
	*Available with F-SW & F-SW-ELV only	
	<input type="radio"/> LEGPANEL-ENCORE: Lower Leg Support ABS Cushions \$98.13	
<input type="radio"/> 2" Extended Footrest Tubes - <i>Only compatible on the F-SW-IKN and F-SW-IKN-ADJ</i>		\$71.45
Encore Pedal with Seat Mounted Hangers	Compatible with 14.5" Seat Height and above	
Encore Pedal with Sub-Seat Mounted Hangers	Compatible with 15.5" Seat Height and above	
Encore Pedal with Rocker with Seat Mounted Hangers	Compatible with 15" Seat Height and above	
Encore Pedal with Rocker with Sub-Seat Mounted Hangers	Compatible with 16" Seat Height and above	

**Common Accessories**

<b>Head</b>	<input type="radio"/> NR-CV	Head Support Visco Foam Contoured Cushion	<input type="radio"/> Triptex	Omit Standard Head Support	\$130.03
<b>Supports:</b>	<input type="radio"/> NS	Neck Support Soft Visco Cushion			\$130.03
<b>Positioning</b>	<input type="radio"/> SR	Lateral Support Reduction Cushions - Pair			\$122.66
<b>Supports:</b>	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports - Pair			\$578.98
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)			\$711.45
<b>Pelvic Belts:</b>	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$63.79
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> CenterBuckle	<input type="radio"/> Side Buckle	\$63.79
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned			\$93.23
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned and Auto Buckle			\$93.23
<b>Thigh Belts:</b>	<input type="radio"/> T-BELT	Thigh Belt Cushioned			\$149.65
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned and Auto Buckle			\$154.82
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene			\$149.65
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene and Auto Buckle			\$154.82
	<input type="radio"/> T-BELT-NEO-APP	Thigh Belt Cushioned Neoprene for APP & SPP Packages			\$149.65
	<input type="radio"/> T-BELT-NEO-A-API	Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages			\$154.82
<b>Trays:</b>	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$382.71
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)		<input type="radio"/> Black Only	\$549.55
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$394.98
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes lower side cover)		<input type="radio"/> Grey <input type="radio"/> Black	\$525.10



# Encore Pedal Wheelchair Model: 48V4/48RV4

US MSRP Order Form  
July 2022

<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover	\$193.82
<b>Other:</b> <input type="radio"/> IV-POLE	IV Pole & Assembly <input type="radio"/> Front or <input type="radio"/> Rear Mount	\$122.66
<input type="radio"/> O2-HOLDER	Oxygen Holder & Assembly <input type="radio"/> Front or <input type="radio"/> Rear Mount	\$206.03
<input type="radio"/> STABIL-FR	Rear Stabilizer Frame	\$110.40
<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each	\$49.70
<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font <input type="radio"/> Block Font	\$30.68
<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag	\$97.52
<input type="radio"/> CUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)	\$55.20

**Additional Options:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

**Delivery Instructions or Special Requests:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Total: \$
-----------

**Accessorial Services:** Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
  - Residential Delivery Service:** Delivery to a residential location
  - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
  - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
  - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.**