

PO:	Date:	Customer #:
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
<ul style="list-style-type: none"> • Solid Seat Pan • Tubular Steel, Powder Coated Frame • Seat cushion, back & head support not included • Independent Height Adjustable Arm Supports • Rocking Mechanism (<i>allows up to 33° of tilt</i>) • Fold Down and Angle Adjustable Back Canes • 33° Infinitely Adjustable Seat Tilt 	<ul style="list-style-type: none"> • 250lb Maximum Occupant Weight • Adjustable Seat Surface Height • 20" Front Mount Mag Wheels • 2 x 5" Rear Total Lock Caster Wheels • Lower Leg Support Assembly • Removable Arm Supports

Encore Rehab Wheelchair® - Model: 500RW MSRP \$3,896.66

WC19 Transit Tie Down Vehicle Package* - Factory Installed Only \$416.82
 *Compatible only with F-SW-ANGLE footrest

Seat and Back Support Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	N/C N/C \$470.41	
Seat Depth:	<input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20"	N/C N/C \$113.13	
Seat Surface	<input type="radio"/> 14" <input type="radio"/> 15" <input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19"	N/C Standard N/C	
Height:			
	<input checked="" type="radio"/> Solid Seat Pan		N/C
	Back Cane	<input type="radio"/> Ultra-Low Double Bend (16") Fixed	Standard
	Height:	<input type="radio"/> Low (18") Fold Down <input type="radio"/> Medium (21") Fold Down <input type="radio"/> Tall (24") Fold Down	N/C N/C N/C



Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair	Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair
<input type="radio"/> 22" Mag Wheels - Pair	N/C	<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels - Pair
<input checked="" type="radio"/> Front Anti-Tip Devices Required Safety Features	N/C	<input type="radio"/> Parking Brake Extension Handles - Pair
<input type="radio"/> Quick Release Axles - Pair	\$142.91	

Arm Support Options			
Arm Support	<input type="radio"/> Arm Supports - Pair	Standard	
Options:	<input type="radio"/> Arm Support Trough - Each <input type="radio"/> Cushioned Arm Supports* - Pair	\$154.82 \$89.32	
	*Not compatible with a tray		
	Arm Support	<input type="radio"/> Height Adjustable Arm Supports	Standard
	Style:	<input type="radio"/> *Height and Depth Adjustable Arm Support	\$154.82
		*Not compatible with 22" mag wheels in the 14", 15", & 16" seat heights	
	Arm Support Mount	<input checked="" type="radio"/> Seat Mounted Arm Supports	Standard

Frame Color Options			
Frame Color:	<input type="radio"/> Black <input type="radio"/> Orange <input type="radio"/> Carmine Red	<input type="radio"/> Tungsten Sand	<input type="radio"/> Ocean Blue
	<input type="radio"/> Bronze <input type="radio"/> Satin Nickel <input type="radio"/> Yellow Zinc		

Drop In Style	Must Choose Lower Leg Support Assembly Style	Pin and Saddle Style
<input type="radio"/> SM-HANGER-48RV4-RKR-IKN (pair)	N/C	<input type="radio"/> SM-HANGER-48RV4-RKR-PIN (pair) <i>Less \$54</i>

Drop In Style	Lower Leg Support Options	Pin and Saddle Style
<input type="radio"/> F-SW-IKN Swing Away Silver (pair)	\$125.05	<input type="radio"/> F-SW Swing Away Lower Leg Support Black - Pair
<input type="radio"/> F-SW-IKN ADJ Swing Away Lower Leg Support Black Pair w/ Angle Adjustable Flip-up Foot Supports	\$239.37	<input type="radio"/> F-SW-ELV Elevating Swing Lower Leg Support Black - Pair
<input type="radio"/> F-SW-ELV-IKN Elevating Swing Away Lower Leg Support Silver w/ Angle Adjustable Flip-up Foot Supports	\$518.04	<input type="radio"/> F-SW-ANGLE Swing Away, 70 Degree Lower Leg Support w/ Angle Adjustable Footplates Black - Pair
<input checked="" type="radio"/> Heel Loops - Pair (for 3 above drop-in style options)	\$37.10	<input type="radio"/> LEG-FOOT-ENCORE Lower Leg Support Foot Board
<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support Panel	\$89.32	<input type="radio"/> LEGPANEL-ENCORE Lower Leg Support Panel

*Slight ground interference may occur in the lowest seat height. Slightly tilt chair or elevate legrest as needed
 2" Extended Footrest Tubes - *Only compatible on the F-SW-IKN and F-SW-IKN-ADJ* \$71.45

Knee-to-Heel Range			
Front Rigging Model	Description	Minimum	Maximum
Range Adjusts in 1" increments			
F-SW-IKN	70° Drop-in, silver	13"	17" (19" with optional extender)
F-SW-IKN-ADJ	70° Drop-in w/ angle adjustable footplates, black	13"	17" (19" with optional extender)
Range Adjusts in 3/4" increments			
F-SW	70° Pin & Saddle, black	11"	14"
F-SW-ANGLE	70° Pin & Saddle w angle adjustable footplates, black	11"	14.75"
F-SW-ELV	Pin & Saddle, elevating, black	15.25"	18.25"
F-SW-ELV-ANG	Pin & Saddle, elevating w/ angle adjustable footplates, black	15.25"	18.25"
Infinite Range Adjustments			
F-SW-ELV-IKN	Drop-in, elevating w/ angle adjustable footplates, silver	13"	20"

"IKN" footrests are compatible with IKON 40, Encore Pedal, Encore Rehab, and Comfort Rehab Wheelchairs
All other footrests are compatible with all Pedal and Rehab Wheelchairs

Common Accessories			
Pelvic Belts:	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle \$61.34
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle \$61.34
	<input type="radio"/> S-BELT-587V3	Pelvic Belt Grommet Mount	<input type="radio"/> Ctr Buckle <input type="radio"/> Side Buckle \$61.34
Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Strap <input type="radio"/> Grey <input type="radio"/> Black \$361.86
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)	<input type="radio"/> Black Only \$521.32
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap \$374.13
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(with upper side pad)	<input type="radio"/> Grey <input type="radio"/> Black \$496.79
	<input type="radio"/> TRAY-COVER	Tray ABS Cushion Cover	\$184.01
Other:	<input type="radio"/> IV-POLE	IV Pole & Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount	\$116.54
	<input type="radio"/> O2-HOLDER	O2 Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount	\$196.27
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag	\$104.94
	<input type="radio"/> STABIL-FR	Rear Stabilizer Frame	\$104.26

Common Cushion Upgrades	
Positioning: <input type="radio"/> SR	Side Reduction Cushions - Pair \$116.54

Matrx® Seating Series

Matrx® Seat Cushions

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions



Seat Cushion Depth:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" (Only available w/ PS Cushion)
----------------------------	--

Seat Cushion Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"
----------------------------	---

<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown \$745	<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown \$1,025
<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$624	<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$681
<input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown \$645	
<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$569	<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$975

Libra Cushion	<input type="radio"/> Libra Fluid Supplement \$120	Accessories:	<input type="radio"/> Pelvic Obliquity Kit \$85
Accessories:	<input type="radio"/> Libra 1" Hemi Wedge \$101		
	<input type="radio"/> Libra 1" Anterior Seat Wedge \$202		
	<input type="radio"/> Libra Tapered Wedge Left \$85		
	<input type="radio"/> Libra Tapered Wedge Right \$85		
	<input type="radio"/> Libra Abductor \$85		
		Cushion Rigidizer:	
		All Cushions	
			16" Depth <input type="radio"/> \$115 18" Depth <input type="radio"/> \$115 20" / 22" Depth <input type="radio"/> \$115
			16" Width <input type="radio"/> \$115 18" Width <input type="radio"/> \$115 20" Width <input type="radio"/> \$115
			16" Width <input type="radio"/> \$115 18" Width <input type="radio"/> \$115 20" Width <input type="radio"/> \$115

		22" Width	N/A	<input type="radio"/> \$115	<input type="radio"/> \$115
--	--	-----------	-----	-----------------------------	-----------------------------

Matrx® Backs

For Correct Fit: Back Size Should Match Selected Chair Seat Dimensions

Back Height: <input type="radio"/> 16" <input type="radio"/> 18" (excludes PB & PB Deep) <input type="radio"/> 20"	
Back Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Back Width: <input type="radio"/> 22"
<input type="radio"/> E2 - Recommended for users seeking optimal immersion and airflow \$965	<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$1,151
<input type="radio"/> E2 Deep - Recommended for users seeking optimal immersion w/ contoured support needs \$1,125	<input type="radio"/> Elite - Recommended for more active users \$1,051
<input type="radio"/> Elite TR - Recommended for tilt in space & recline applications \$866	<input type="radio"/> Elite Deep - Recommended for active users with contoured support \$1,095
<input type="radio"/> Elite - Recommended for more active users \$766	Matrx Back <input type="radio"/> Air Force N/C
<input type="radio"/> Elite Deep - Recommended for active users with contoured support needs \$918	Military <input type="radio"/> Army N/C
<input type="radio"/> PB - Recommended for less active users seeking comfort \$734	Service <input type="radio"/> Marine N/C
<input type="radio"/> PB Deep - Recommended for less active users seeking comfort & support \$892	Decals: <input type="radio"/> Navy N/C
E2 Back <input type="radio"/> E2 Fixed Lateral Supports 6.5"Hx5"L (Pair) \$385	Privacy Flap <input type="radio"/> 16" Width Small Privacy Flap \$53
Accessories: <input type="radio"/> E2 Chest Harness Interface \$85	Accessories: <input type="radio"/> 18" Width Medium Privacy Flap \$53
<input type="radio"/> E2 Contoured Backbone Lumbar/PSIS Support \$75	<input type="radio"/> 20" Width Large Privacy Flap \$53
<input type="radio"/> E2 PostureFlex Abdominal Support 7"Hx20"-32"L \$295	<input type="radio"/> 22" Width Extra Large Privacy Flap \$53
<input type="radio"/> Quad Release Tab for PostureFlex \$35	Lateral Support <input type="radio"/> Elite Back Fixed Laterals 5"Hx7.5"L \$116
Back <input type="radio"/> E2 Infection Control Slipcover \$95	Support <input type="radio"/> Elite Back Swing Away Lateral Hardware \$243
Slipcovers: <input type="radio"/> E2 Deep Infection Control Slipcover \$95	Accessories: <input type="radio"/> Swing Away Lateral Pad (Each) 4.25"Hx6.5 \$122
<input type="radio"/> Infection Control Slipcover for Matrx Cushions \$95	

Matrx® Headrest

Head Support Mounting: <input type="radio"/> Elan Head Support Hardware - Required with Selection of Head Support \$285
Choose One <input type="radio"/> Loxx Head Support Hardware - Required with Selection of Head Support \$275
Support Options: <input type="radio"/> Head Support Cushion <input type="radio"/> 10" <input type="radio"/> 14" \$310
<input type="radio"/> Head Support 4 Point Cushion <input type="radio"/> 9" x 7" <input type="radio"/> 11" x 10" \$448
<input type="radio"/> Headband with 4 Point Cushion <input type="radio"/> Small <input type="radio"/> Large \$610

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorial services needed to deliver are requested when an order is submitted to Broda.

Lift Gate Service: Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level

Residential Delivery Service: Delivery to a residential location

Notification Before Delivery: Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery

Inside Delivery Service: Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs

Limited Access Service: Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.