

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
• Comfort Tension Seating®	• 250lb Maximum Occupant Weight
• Tubular Steel, Powder Coated Frame	• Adjustable Seat Surface Height
• Up to 18° Infinitely Adjustable Seat Tilt	• 2 x 20" Front Mag Wheels
• Independent Height Adjustable Arm Supports	• 2 x 5" Total Lock Rear Caster Wheels
• Rocking Mechanism (allow for up to 33° of tilt)	

Latitude Pedal Wheelchair - Model: 48R MSRP \$3,282.48

Seat and Back Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C	Seat Surface	<input type="radio"/> 14" N/C
	<input type="radio"/> 22" \$512.73	Height:	<input type="radio"/> 15" Standard
Seat Depth:	<input type="radio"/> 15.5" N/C		<input type="radio"/> 16" <input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" N/C
	<input type="radio"/> 17" Standard	Back Support	<input type="radio"/> 32" Standard
	<input type="radio"/> 18.5" N/C	Height:	<input type="radio"/> 33.5" \$122.66
	<input type="radio"/> 20" \$122.66		<input type="radio"/> 35" \$122.66

Front Options	Wheel Options	Rear Options
<input type="radio"/> 20" Mag Wheels - Pair Standard	<input type="radio"/> 5" Total Lock Caster Wheels - Pair Standard	
<input type="radio"/> 22" Mag Wheels - Pair N/C		
OR		
<input type="radio"/> 5" Total Lock Caster Wheels - Pair* N/C	<input type="radio"/> 20" Mag Wheels - Pair N/C	
<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels - Pair \$233.06	<input type="radio"/> 22" Mag Wheels - Pair N/C	
OR		
<input type="radio"/> 5" Total Lock Caster Wheels - Pair N/C	<input type="radio"/> 5" Fixed Caster Wheels - Pair (\$680.01)	
	<input type="radio"/> Parking Brake Extension Handles - Pair \$78.51	

Arm Support Options			
Arm Support	<input type="radio"/> Arm Supports - Pair Standard	Arm Support	<input type="radio"/> Height Adjustable Arm Support Standard
Options:	<input type="radio"/> Arm Support Trough - Each \$166.82	Style:	<input type="radio"/> Height and Depth Adjustable \$169.27
	<input type="radio"/> Cushioned Arm Supports* - Pair \$98.13	Armtop	<input checked="" type="radio"/> Base Mounted Arm Supports Standard

*Not compatible with a tray

Seating and Frame Color Options			
Strap Color:	<input type="radio"/> Grey <input type="radio"/> Black <input type="radio"/> Blue N/C	Frame Color:	<input type="radio"/> Grey <input type="radio"/> Black N/C

Upholstery Color Options (Select one)			
Traditional Upholstery Colors:	<input type="radio"/> Burgundy <input type="radio"/> Charcoal <input type="radio"/> Sapphire Blue N/C		

*Traditional Upholstery Colors come with seat and back pad in matching vinyl color.

New Upholstery Colors:

Vinyl:	<input type="radio"/> Interlace Indigo**	<input type="radio"/> Interlace Graphite**	<input type="radio"/> Savoy Black	<input type="radio"/> Interlace Espresso**	
	<input type="radio"/> Kaleidoscope Smokehouse	<input type="radio"/> Kaleidoscope Chocolate			N/C
Vinyl:	<input type="radio"/> Beeline Portobello	<input type="radio"/> Beeline Geyser	<input type="radio"/> Beeline Fossil	<input type="radio"/> Beeline Forge	N/C
Nylon:	<input type="radio"/> Digital Camo	<input type="radio"/> Retro Camo			\$159.47

Upholstery Colors come with seat and back cushions in savoy black vinyl. *Interlace: Indigo, Espresso, and Graphite available with seat and back cushions in matching vinyl color



*Traditional Upholstery



**New Upholstery Color (Digital Camo Shown)



***New Upholstery Color (Espresso Shown) with option of seat and back pad in matching vinyl color.

Cushions - Choose Option 1, 2, OR 3

<input type="radio"/> Option 1 - Broda Cushions		N/C
Standard Broda Cushions:		
Includes: Neck Support, Lateral Shoulder Supports, Lower Side Covers, Seat & Back Cushions with Covers		
Seat & Back Cushion	<input type="radio"/> Black Vinyl (with New Upholstery Colors Only)	Standard
Cover Options:	<input type="radio"/> Accent Color	N/C
	*Standard with Traditional Upholstery Colors, or ***Option For Interlace: Indigo, Espresso or Graphite only	
	<input type="radio"/> Full Vinyl Wrap Seat & Back Cushions	\$122.66
	(Select Black Vinyl OR Accent Color Vinyl Above)	
	<input type="radio"/> Grey Terry Cover	N/C
	<input type="radio"/> White Terry Cover	N/C
Lateral Shoulder	<input type="radio"/> Standard Short	N/C
	<input type="radio"/> Rounded Short	N/C
	<input type="radio"/> ACC Short	N/C
Support Options:	<input type="radio"/> Standard Long	\$166.82
	<input type="radio"/> Rounded Long	\$166.82
	<input type="radio"/> ACC Long	\$166.82

Optional Substitutions Available Between APP & SPP Components

<input type="radio"/> Option 2 - Broda Additional Positioning Padding (APP) (not compatible with 35" back support height)		
Broda APP Package Includes: Seat, Back and Head Support Cushions (no Head Support or Lateral Shoulder Supports)		
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	Seat Depth: <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21" <input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"
		\$451.41
		\$706.54
<input type="radio"/> Option 3 - Broda Slim Posture Padding (SPP)		
Broda SPP Cushion Options: Seat, Back Support, and Back Support Extension Cushions		
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	Seat Depth: <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20" <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"
		\$161.97
		\$178.64
Back:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22"	Back Extension: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"
	\$142.91	\$107.19
13" Height	\$171.49	12" Height <input type="radio"/> 22" \$119.09

Common Accessories

<i>Broda recommends ordering chair with Lower Leg Supports for safe patient handling</i>			
Lower Leg	<input type="radio"/> F-SW-ELV	Elevating Swing Away Lower Leg Supports Black - Pair	\$316.48
Supports:	<input type="radio"/> F-SW	Swing Away Footrests Lower Leg Supports Black - Pair	\$252.69
	<input type="radio"/> F-SW-ANGLE	Pin and Saddle, Swing Away, 70 Degree Lower Leg Support with Angle Adjustable Footplates Black - Pair	\$316.48
	<input type="radio"/> LEG-FOOT-48/48R	Lower Leg Support Foot Board	\$291.95
	<input type="radio"/> LEGPANEL-48	Lower Leg Support Panel Pad Black	\$98.13

Lower Leg Support Assembly:	<i>(allows for lower foot support to tilt with chair)(configurations specific - please contact Customer Support for full details)</i>					
<input type="radio"/> SM-HANGER	Lower Leg Support Assembly Flared Seat Mount - Short			\$142.29		
<input type="radio"/> SM-HANGER	Lower Leg Support Assembly Flared Seat Mount - Long			\$142.29		
<input type="radio"/> SM-HANGER	Lower Leg Support Assembly Straight Seat Mount - Short			\$142.29		
<input type="radio"/> SM-HANGER	Lower Leg Support Assembly Straight Seat Mount - Long			\$142.29		
Neck Supports:	<input type="radio"/> NR-CV	Neck Support Visco Foam Contoured Cushion	<input type="radio"/> Triptex	<input type="radio"/> Omit Standard	\$130.03	
	<input type="radio"/> NS	Neck Support Soft Visco Cushion		Neck Support	\$130.03	
Positioning:	<input type="radio"/> SR	Lateral Support Reduction Cushions - Pair			\$122.66	
	<input type="radio"/> SAS	Side Arm Supports - Pair			\$161.92	
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports - Pair			\$578.98	
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED) - Pair			\$711.45	
Pelvic Belts:	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$63.79	
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$63.79	
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned			\$93.23	
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned & Auto Buckle			\$93.23	
Thigh Belts:	<input type="radio"/> T-BELT	Thigh Belt Cushioned			\$149.65	
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned & Auto Buckle			\$154.82	
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene			\$149.65	
	<input type="radio"/> T-BELT-NEO-A	Thigh Belt Cushioned Neoprene Auto Buckle			\$154.82	
	<input type="radio"/> T-BELT-NEO-APP	Thigh Belt Cushioned Neoprene for APP & SPP Packages			\$149.65	
	<input type="radio"/> T-BELT-NEO-A-APF	Thigh Belt Cushioned Neoprene Auto Buckle for APP & SPP Packages			\$154.82	
Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey	<input type="radio"/> Black	\$382.71
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap			\$394.98
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes upper side pad)				\$525.10
	<input type="radio"/> TRAY-COVER	Cushioned ABS Tray Cover				\$193.82
Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one)(pair)			\$110.40	
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers			\$193.82	
Other:	<input type="radio"/> IV-POLE*	IV Pole and Assembly	<input type="radio"/> Front	<input type="radio"/> Rear Mount	\$122.66	
	<input type="radio"/> O2-HOLDER*	Oxygen Holder and Assembly	<input type="radio"/> Front	<input type="radio"/> Rear Mount	\$206.03	
	*IV-POLE and O2-HOLDER not compatible with rear mounted mag wheels					
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each			\$49.70	
	<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font OR <input type="radio"/> Block Font			\$30.68	
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag			\$97.52	
	<input type="radio"/> CUP-HOLDER	Cup Holder (not compatible with front mounted mag wheels)			\$55.20	

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order,**

any accessorial services needed to deliver are requested when an order is submitted to Broda.

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.