

PO: _____

Date: _____

Customer #: _____

Customer Name: _____

Ship To Name: _____

Address: _____

Address: _____

Phone: _____

Phone: _____

Contact: _____

Contact: _____

Email: _____

Tag: _____



Standard Features

- Up to 12° Infinitely Adjustable Seat Tilt (Posterior)
- Up to 10° Infinitely Adjustable Seat Tilt (Anterior)
- Up to 67° Infinitely Adjustable Back Support Recline
- Stainless Steel, Powder Coated Frame
- 3 x 5" Total Lock Stainless Steel Caster Wheels
- 1 x 5" Rear Directional Lock Stainless Steel Caster Wheels
- Angle Adjustable, Articulating Single Lower Leg Support
- Swing Away Removable Arm Supports
- Comfort Tension Seating®
- 400lb Maximum Occupant Weight
- Head Support

Revive Shower Commode - Model: CS385 MSRP \$3,878.62

Seat and Back Support Options

Seat Width: 18" 22" N/C

Seat Depth: 18" Standard

Seat Height: 24.5" Standard

Back Height: 26" Standard

No Anterior Tilt*: N/C

*Max Upright Position - 5 degree posterior tilt

Wheel Options

5" Rear Directional Lock Stainless Steel Caster Wheels & 5" Total Lock Stainless Steel Caster Wheels Standard

Seat Options (Select one)

Seat Commode Cushioned Standard

Seat Transport Cushioned N/C

Seat Small Hole Cushioned N/C

Seat Large Hole Cushioned 22" Seat width Only N/C

Bedpan & Basket Options

Bedpan with Basket Standard

Bedpan - Aqua N/C

Bedpan - Short Handle N/C

Basket Frame - Aqua N/C

Arm Support Options

Arm Support Options: Arm Supports Standard

Arm Support Style: Arm Support Top-Offset-CS385 Seat Width Extending (extends width to 25") \$323.84

Head Support

Head Support: Molded 2" Offset Standard

Strap Color

Strap Color: Grey Blue Green N/C

Common Accessories

Lateral Supports: Molded Laterals - Pair *Not recommended with 18" seat width \$98.13

Pelvic Belt: Pelvic Belt (side buckle) (18"= small, 22"= medium) \$63.79

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. **To ensure a seamless delivery of each Broda order, any accessorial services needed to deliver are requested when an order is submitted to Broda.**

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
- Residential Delivery Service:** Delivery to a residential location
- Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
- Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
- Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p

Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.