

PO: _____

Date: _____

Customer #: _____

Customer Name: _____
Address: _____
Phone: _____
Contact: _____
Email: _____

Ship To Name: _____
Address: _____
Phone: _____
Contact: _____
Tag: _____



Standard Features	
• Comfort Tension Seating®	• 15° Adjustable Back Recline
• Tubular Steel, Powder Coated Frame	• 2 x 5" Total Lock Front Caster Wheels
• Infinitely Adjustable Seat Tilt to 40°	• 2 x 20" Mag Rear Wheels
• 350lb Maximum Occupant Weight	

Elite Rehab Wheelchair - Model: 550 SR MSRP \$4,428.16

- WC-19 Transport Package** WC-19 Compliant Transit Tie Down Vehicle Transport Package - Factory Install Only \$546.63
- Acute Care Configuration:** Includes: Full Vinyl Covers, ACC Lateral Shoulder Supports, ABS Lower Leg Support Foot Cushions (Factory Installed)(Skip lower leg support and cushion cover/wrap sections below when selecting this option, 16"-24" widths only) \$552.00

Seating Frame Color Options

- Strap Color:** Grey Black Blue N/C **Cushion** Burgundy Charcoal N/C
Frame Color: Grey Black N/C **Color:** Meteor Sapphire Blue N/C

Seat Options

- Seat Width:** 16" or 18" N/C **Seat Depth:** 17" Standard
 20" N/C
 22" \$484.53
 24" 26" \$1,097.84
 28" \$1,097.84
Seat Surface 19.5" N/C
Height: 20" (ONLY option w/ WC-19 package) Standard
 20" N/C
 22" N/C

- Seat Option:** Comfort Tension Seating® (CTS) Standard Solid Seat Pan \$166.95
See CTS Seat Options below See Matrx Seat Cushion Options below

Comfort Tension Seating® (CTS) ONLY: Seat Cushion Options

- Broda Seat Cushions Includes:** Vinyl Cover Terry Cover (Grey or White) N/C
 Full Vinyl Wrap Seat \$58.28
- Broda Additional Positioning Padding (APP) Seat Cushion** 16" 18" 20" **Seat** 17" 19" 21" \$239.20
Seat Width 22" 24" 26" 28" **Depth:** 17" 19" 21" \$252.68
- Broda Slim Posture Padding (SPP) Seat Cushion**
- Seat Width:** 16" 18" 20" **Seat Depth:** 17" 18.5" 20" \$211.98
 22" 24" 17" 18.5" 20" \$233.42
 26" 17" 18.5" 20" \$262.00

Solid Seat Pan: Matrx® ONLY: Seat Cushion Options



Please note, when ordering the Seat cushion only please contact Broda for additional shipping costs

Matrx® Seat Cushions

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth: 16" 18" 20" 22" (Only available w/ PS Cushion)

Seat Cushion Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Cushion Width: <input type="radio"/> 22"
<input type="radio"/> Libra - Recommended for highest risk of or existing skin breakdown \$797.15	<input type="radio"/> Libra -Recommended for highest risk or existing skin breakdown \$1,096.75
<input type="radio"/> Vi - Recommended for moderate to high risk of skin breakdown \$667.68	<input type="radio"/> Vi -Recommended for moderate to high risk of skin breakdown \$728.67
<input type="radio"/> PSP - Recommended for moderate to high risk of skin breakdown \$690.15	
<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$608.83	<input type="radio"/> PS - Recommended for low to moderate risk of skin breakdown \$1,043.25

Libra Cushion Accessories:	<input type="radio"/> Libra Fluid Supplement \$128.40	Vi, PSP & PS Cushion	<input type="radio"/> Pelvic Obliquity Kit \$90.95
	<input type="radio"/> Libra 1" Hemi Wedge \$108.07		
	<input type="radio"/> Libra 1" Anterior Seat Wedge \$216.14	Cushion Rigidizer:	
	<input type="radio"/> Libra Tapered Wedge Left \$90.95	All Cushions	
	<input type="radio"/> Libra Tapered Wedge Right \$90.95		
	<input type="radio"/> Libra Abductor \$90.95		

For Correct Fit: Cushion Size Should Match Selected Chair Seat Dimensions

Seat Cushion Depth: 16" 18" 20" 22" (Only available w/ PS Cushion)

Seat Cushion Width: <input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Cushion Width: <input type="radio"/> 22" <input type="radio"/> 24"(PS &Libra) <input type="radio"/> 26" (PS only)
<input type="radio"/> Libra-Recommend for highest risk of or existing skin breakdown \$745.00	<input type="radio"/> Libra-Recommend for highest risk of or existing skin breakdown \$1,025
<input type="radio"/> Vi -Recommend for moderate to high risk of skin breakdown \$624.00	<input type="radio"/> Vi -Recommend for moderate to high risk of skin breakdown \$681
<input type="radio"/> PSP -Recommend for moderate to high risk of skin breakdown \$645.00	
<input type="radio"/> PS - Recommend for low to moderate risk of skin breakdown \$569.00	<input type="radio"/> PS -Recommend for low to moderate risk of skin breakdown \$975

Libra Cushion Accessories:	<input type="radio"/> Libra Fluid Supplement \$120.00	Vi, PSP & PS Cushion Accessories:	<input type="radio"/> Pelvic Obliquity Kit \$85.00
	<input type="radio"/> Libra 1" Hemi Wedge \$101.00		
	<input type="radio"/> Libra 1" Anterior Seat Wedge \$202.00	Cushion Rigidizer:	
	<input type="radio"/> Libra Tapered Wedge Left \$85.00	For All Cushions	
	<input type="radio"/> Libra Tapered Wedge Right \$85.00		
	<input type="radio"/> Libra Abductor \$85.00		

Infection Control Slipcover for Matrx Cushion: \$95

Back Options

Back Height: 32" Standard **Back Seating:** Comfort Tension Seating® (CTS) Standard

Back Recline: Back Recline Semi (0° - 15°) Standard Back Recline Full(>15° - ≤ 80°) N/C

Back Cushion Options

Broda Back Padding N/C

Include: Upper & Lower Side Covers, Back Cushion with Cover, Foot Cushions

Vinyl Back Terry Back (Grey or White) N/C

Full Vinyl Wrap Back \$58.28

Optional Substitutions Available Between APP & SPP Components

<input type="radio"/> Additional Positioning Padding (APP) Back Cushion	Back Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$179.10
			<input type="radio"/> 22"	<input type="radio"/> 24"	\$230.60
			<input type="radio"/> 26"	<input type="radio"/> 28"	\$257.59

Broda SPP Back Cushion

Back:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$169.11	Back Extension:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$125.05
13" Height	<input type="radio"/> 22"	<input type="radio"/> 24"		\$191.73	12" Height	<input type="radio"/> 22"	<input type="radio"/> 24"		\$135.76
	<input type="radio"/> 26"			\$202.46		<input type="radio"/> 26"			\$147.67

Lateral Shoulder	<input type="radio"/> Standard Short	<input type="radio"/> Rounded Short	<input type="radio"/> ACC Short	N/C
Support Options:	<input type="radio"/> Standard Long	<input type="radio"/> Rounded Long	<input type="radio"/> ACC Long	\$159.46
	<input type="radio"/> APP Shoulder Supports			\$116.54

Common Padding Upgrades

Huntington's Padding Package Includes: Long Lateral Shoulder Supports, Thigh Belt & Long Lower Foot Support Footbox \$613.84

Positioning:	<input type="radio"/> NR-CV	Head Support Visco Foam Contoured Pad	<input type="radio"/> Triptex	<input type="radio"/> Omit Standard Headrest	\$122.66		
	<input type="radio"/> NS	Neck Support Soft Visco Cushion			\$122.66		
	<input type="radio"/> Additional Positioning Padding (APP) Head Support		Seat Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$73.61
				<input type="radio"/> 22"	<input type="radio"/> 24"		\$85.87
				<input type="radio"/> 26"	<input type="radio"/> 28"		\$98.13
	<input type="radio"/> SR	Side Reduction Pads (pair)					\$116.54
	<input type="radio"/> SAS	Side Arm Supports (pair)					\$143.58
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports					\$547.09
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)					\$674.66

Front Mount Configuration

Wheel Options

Rear Mount Configuration

<input type="radio"/> 5" Total Lock Caster Wheels - Pair	Standard	<input type="radio"/> 20" Mag Wheels - Pair	Standard
<input type="radio"/> 5" Locking Twin Caster Wheels - Pair	\$233.06	<input type="radio"/> 14" Mag Wheels - Pair	N/C
<input type="radio"/> 5" Aluminum Non-Locking Caster Wheels* Pair *(available with rear mags only)	\$233.06	<input type="radio"/> 22" Mag Wheels - Pair*	N/C
		*(22" mag may slightly limit tilt depending on arm support height)	
<input checked="" type="radio"/> Anti-Tip Devices(factory required safety feature)	N/C	<input type="radio"/> Parking Brake Extension Handles (E0961)- Pair	\$78.50

Lower Leg Support Options (Select One)

<input type="radio"/> Angle Adjustable, Articulating Split Lower Leg Support (not available in 16" seat width)	Standard
Choose Number of Straps (foot support length range)	<input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Angle Adjustable, Articulating Single Lower Leg Support	N/C
Choose Number of Straps (foot support length range)	<input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Angle Adjustable, Articulating ABS Lower Leg Support (includes foot cushions 16"-24" widths)	\$398.66

Foot Cushions:	<input type="radio"/> FB-LONG	Lower Leg Support Footbox - Full Extended Cushion	<input type="radio"/> Omit Foot Support Cushions	\$398.66			
	<input type="radio"/> APP CALF/FOOT - APP Calf and Foot Support		Seat Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	\$177.86
				<input type="radio"/> 22"	<input type="radio"/> 24"		\$220.80
				<input type="radio"/> 26"	<input type="radio"/> 28"		\$252.68

Arm Support Options (Select One)

Armrest Options:	<input type="radio"/> Arm Supports - Pair	Standard
	<input type="radio"/> Long Arm Supports - Pair	\$49.06
	<input type="radio"/> Arm Support Trough - Each	\$159.46
Arm Support Style:	<input type="radio"/> Seat Mounted Adjustable Removable Arm Supports	Standard
	<input type="radio"/> Height and Depth Adjustable - Pair	\$159.46

Common Accessories

Pelvic Belts:	<input type="radio"/> S-BELT	Pelvic Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-A	Pelvic Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$61.34
	<input type="radio"/> S-BELT-PD	Pelvic Belt Cushioned			\$88.32
	<input type="radio"/> S-BELT-PD-A	Pelvic Belt Cushioned & Auto Buckle			\$88.32
Thigh Belts:	<input type="radio"/> T-BELT	Thigh Belt Cushioned			\$141.06
	<input type="radio"/> T-BELT-A	Thigh Belt Cushioned & Auto Buckle			\$147.19
	<input type="radio"/> T-BELT-NEO	Thigh Belt Cushioned Neoprene			\$141.06

Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (<i>select one</i>)(pair)	\$104.26
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers	\$184.01
Other:	<input type="radio"/> IV-POLE	IV Pole and Assembly - Front Mount	\$116.54
	<input type="radio"/> O2-HOLDER	Oxygen Holder and Assembly - Front Mount	\$196.27
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each	\$49.06
	<input type="radio"/> BAG PERSONALIZING	(2 Initials) <input type="radio"/> Brush Font	\$30.67
	<input type="radio"/> O2-HOLDER-BAG	Oxygen Holder, Soft Carry Bag	\$104.94
	<input type="radio"/> CUP-HOLDER	Cup Holder	\$55.19

Additional Options:

Delivery Instructions or Special Requests:

Total: \$

Accessorial Services: Each with an additional charge, are available to accommodate the delivery needs of Broda customers and/or the needs of a shipment's final delivery location. To ensure a seamless delivery of each Broda order,

any accessorial services needed to deliver are requested when an order is submitted to Broda.

- Lift Gate Service:** Delivery is made with a lift gate trailer that moves the shipment from trailer level to ground level
 - Residential Delivery Service:** Delivery to a residential location
 - Notification Before Delivery:** Contact is requested to the Consignee by the Carrier PRIOR to shipment delivery
 - Inside Delivery Service:** Movement of the shipment beyond the unloading position, i.e. move closer to a door, inside a building, upstairs, down a hall, etc. Price is specific per order, based on the customer's inside delivery needs
 - Limited Access Service:** Delivery locations that are not open to the walk-in public, during normal business hours. Mon-Fri. 8a-5p
- Expedited delivery is available by customer request at the time of order, which is priced per order and is based on the delivery location and level of service.