

BRODA[®]

Hygiene Seating Solutions



BRODA[®]

brodaseating.com
1-844-552-7632

When it comes to hygiene care, the right shower commode chair can make the difference between a negative bathing experience and a positive one for both patient and caregiver. That's why Broda shower commodes have unique features to make hygiene comfortable and safe for everyone involved.

TABLE OF CONTENTS

<i>Broda Benefits</i>	3
<i>Why Broda?</i>	4
<i>Sentinel Shower Commode</i>	5
<i>Revive Shower Commode</i>	6
<i>Revive Bariatric Shower Commode</i>	7
<i>Configure Your Shower Commode</i>	8
<i>Accessories</i>	8



BRODA BENEFITS

At Broda, we promise to promote and preserve the wellbeing of the people who use our products. You can be sure your Broda chair or wheelchair will provide you with the following benefits.



**SEATING
COMFORT**



POSITIONING



MOBILITY



**PATIENT
SAFETY**



**CAREGIVER
SAFETY**



**INFECTION
CONTROL**



DURABILITY



PERSONALIZATION



**EASY
ORDERING**

THE IMPORTANCE OF HYGIENE

According to studies on caregiver burden, bathing assistance has been shown to be the most taxing aspect of care in nursing careers.¹


The most stressful aspects of bathing include...

- Up to five manual transfers
- Challenges dressing & undressing
- Awkward postures and movements
- Struggle to maintain upright patient posture



Hygiene Care Challenges

- Transferring a patient in order to bathe them may be difficult and uncomfortable
- Basic stock equipment may be uncomfortable, overused, worn out, and easily broken
- Patients & caregivers face fall risks in tight, slippery spaces
- Patients are vulnerable and need to maintain their dignity



1: Yamamoto, A., Azuma, T., Barrameda, H. C., Kida, N., Goto, A., & Ota, T. (2015). Caregiver's Eye Gaze and Field of View Presumption Method During Bathing Care in Elderly Facility. In V. G. Duffy (Ed.), Digital Human Modeling. Applications in Health, Safety, Ergonomics and Risk Management: Ergonomics and Health (Vol. 9185, pp. 524-532). Springer International Publishing. <https://doi.org/10.1007/978-3-319-21070-4>

THE BRODA HYGIENE DIFFERENCE

1. Safer Lifts and Transfers

Broda shower commodes enable safe transition options. Plus, the ease of maneuvering the patient reduces the need for adjustments and transfers.

2. Cleanliness Leads to Happiness

Comfort and dignity during the bathing process increases satisfaction with personal hygiene. This leads to better overall health and quality of life.

3. Cleanability

Easy-to-remove, minimally textured surfaces make Broda products and their components, such as bed pans and lateral supports, easy to clean.

4. Ergonomic Design

Broda shower commode chairs and accessories can adapt for a patient in multiple environments or stages of mobility, providing safety throughout the care continuum.

5. Durable Surfaces

Broda shower commodes feature durable, rust-resistant components that are easy to clean and disinfect between bathing sessions. Unlike basic shower commodes, Broda products provide safe, comfortable bathing solutions that last for years.

6. 10-Year Frame Warranty

Quality products come with quality warranties. Broda guarantees our shower commode frames with an industry-leading 10-year warranty. While other shower commodes break frequently and need replacement, Broda shower commodes stand the test of time.

How Broda Meets Hygiene Needs

- Dynamic tilt & recline provides safe, thorough, and efficient bathing
 - Maximum durability with stainless steel frames and rust-resistant parts
 - Removable arm supports create more space for safe transfers
 - Unparalleled seating comfort reduces agitation and resistance to bathing
-



*Broda's **Revive Shower Commode** meets our patient's needs. He indicates he feels safe and he is able to participate in his own care each morning. He is **safe at all times** while in the chair, and the caregivers are safe at all times while providing his care. All staff report being pleased with the commode and find it much easier to assist him.*

*- Pam
Director of Nursing*



Sentinel Shower Commode

*SAFE, COMFORTABLE,
INDEPENDENT BATHING*



Shown with standard mag wheel configuration and optional commode pan and basket

PRODUCT FEATURES

BACK ANGLE	Adjustable 5°, 10°, 15°
SEAT ANGLE	Adjustable 0° - 6°
SEAT WIDTH	18", 20", 22"
SEAT DEPTH	Adjustable 16 - 19"
WEIGHT CAPACITY	Up to 300 lbs.



Features a range of simple adjustments for seat depth, wheel position, seat height, and more.



2-in-1 solid seat and commode is highly-versatile and easily mountable in all four directions



Configurable with either user-propelled mag wheels or cost-effective casters



Flip-up arm supports and swing-away leg supports allow for easy bathing and transfers

Revive Shower Commode

COMPREHENSIVE BATHING
SYSTEM WITH TILT & RECLINE
POSITIONING FEATURES



*Anterior and posterior tilt aid in
fall prevention and safe
patient handling*



*Multi-positional functionality
with tool-free adjustability for
ease of use*



*Multiple configurations available to
accommodate unique needs*



*Rust-resistant frame and
components, with a 10-year
warranty*

PRODUCT FEATURES

DEGREES OF TILT	Up to 12° Posterior, Up to 10° Anterior
DEGREES OF RECLINE	Up to 67°
SEAT WIDTH	18", 22"
SEAT DEPTH	18"
WEIGHT CAPACITY	Up to 400 lbs.

Revive *Bariatric Shower Commode*

*BARIATRIC HYGIENE CARE IN
A CLASS ALL ITS OWN*



Anterior tilt aids in patient transfers and ease of patient care



Mid-mounted Mag Wheel provides increased maneuverability



Optional concave back accommodates and supports gluteal shelf



Rust-resistant frame and components are durable and long-lasting

PRODUCT FEATURES


DEGREES OF TILT	Up to 14° Posterior, Up to 8° Anterior
DEGREES OF RECLINE	Up to 67°
SEAT WIDTH	27", 33" with width adjustable armrests
SEAT DEPTH	21", 23"
WEIGHT CAPACITY	Up to 800 lbs.

Configure Your Hygiene Wheelchair

Features, Functionality, and Customization Level


Bariatric

Revive Bariatric Shower Commode




More Independent

More Dependent



Sentinel Shower Commode



Revive Shower Commode

300 lbs. weight capacity 5°, 10°, 15° Tilt	400 lbs. weight capacity Up to 12° Posterior Tilt Up to 67° Recline
---	---

OPTIONS & ACCESSORIES



Molded Lateral Supports

COMPANION WHEELCHAIR FOR DAILY USE



Encore Rehab Wheelchair

The Encore Rehab Wheelchair offers a universal seating system, customizable for each user's seating and positioning needs. With 33 degrees of infinitely-adjustable tilt and an optional dynamic rocking feature, the Encore Rehab is perfect for a wide variety of care environments.



Commode pan with Basket



Pelvic Positioning Strap



Broda USA

15455 Conway Rd, Suite 200
Chesterfield, MO 63017, USA

Contract # 36F79721D0039

BRODA[®]

Essential Comfort

brodaseating.com
1-844-552-7632

Broda Canada

560 Bingemans Centre Dr
Kitchener, ON N2B 3X9, Canada

 **Ontario**
Assistive Devices Program