

# BRODA<sup>®</sup> WC19 Transport Package

## *Elite Positioning Wheelchair*



*Rear mount mag wheels facilitate hand propulsion allowing patients to engage in their surroundings*



*Adjustable tilt and recline helps decrease postural deviations such as slumping, lateral lean and head drop*



*Posterior tilt effectively opens the diaphragm improving digestion, oxygenation, and organ function*



*Comfort Tension Seating<sup>®</sup> improves postural support and sitting tolerance*

### ***Patient, Facility, & Caregiver Benefits***

- Facilitates positioning for comfort and respiratory function.
- Tilt and recline provides ease of use for frequent re-positioning.
- Proper seating alignment decreases postural deviations.

### ENHANCED POSITIONING AND FUNCTIONALITY

#### UPHOLSTERY OPTIONS



#### FRAME OPTIONS



#### SPECIFICATIONS

|                                   |  |
|-----------------------------------|--|
| Degrees of Tilt                   | up to 40°  |
| Degrees of Recline                | up to 90°  |
| Weight Capacity                   | 300 lbs.<br>(158 kg)   |
| Seat Widths                       | 16", 18", 20", 22", 24", 26", 28"<br>(40cm, 45 cm, 50cm, 55cm, 60cm, 65cm) |
| Seat Depths                       | 15.5", 17", 18.5", 20"<br>(37.5cm, 43cm, 46cm, 50cm)                       |
| Floor-to-Seat Heights             | 19", 20.5", 22", 23.5"<br>(48cm, 51cm, 55cm, 59cm)                         |
| Wheels                            | 5" (12cm) Total Lock Front Casters<br>20" or 22" (50cm or 55cm) Rear Mags  |
| Lifetime Warranty Frame           | 10 years on 16 Gauge tubular steel frame                                   |
| Pressure-Reducing Seat & Backrest | Integrated with Comfort Tension Seating®                                   |
| Standard Padding*                 | Lower side pads, upper side pads, seat and back pad with cover             |

Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps

[brodaseating.com](http://brodaseating.com) | 1-844-552-7632