

# BRODA<sup>®</sup>

## Elite

### *Positioning Wheelchair*



*Tilt-in-Space* functionality allows easy repositioning of users




*Removable armrests* make transfers safer for caregivers and patients



*Posterior Tilt* effectively opens the diaphragm improving digestion, oxygenation, and organ function

### ***Patient, Facility, & Caregiver Benefits***

- Facilitates positioning for comfort and respiratory function
- Tilt and recline provides ease of use for frequent re-positioning
- Proper seating alignment decreases postural deviations



*Comfort Tension Seating<sup>®</sup>* increases support and sitting tolerance

**brodaseating.com | 1-844-552-7632**

### ENHANCED POSITIONING AND FUNCTIONALITY


#### UPHOLSTERY OPTIONS



#### SPECIFICATIONS

Degrees of Tilt	Up to 40° Posterior
Weight Capacity	Up to 350 lbs. (159 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45 cm, 50cm, 55cm, 60cm, 65cm)
Seat Depths	15.5", 17", 18.5", 20", 22" (37.5cm, 43cm, 46cm, 50cm, 56cm)
Seat Heights	19", 20.5", 22", 23.5" (48cm, 51cm, 55cm, 59cm)
Wheels	5" Total Lock Front Casters (13cm) 20" or 22" Rear Mag Wheels (50cm or 55cm) 6" Multi-directional Locking Casters (15cm)
Lifetime Frame Warranty	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Lower side pads, upper side pads, seat and back pad with cover neckrest, upper laterals and, sole/calf pads

#### FRAME OPTIONS



#### OPTIONAL ACCESSORIES



O<sup>2</sup> Holder



Lateral Support



Half Tray



Clear Tray



Padded Thigh Strap

Additional accessories available: foot supports, trays, and positioning straps