

BRODA[®] WC19 Transport Package

Encore Pedal Wheelchair



Shown with short back and Matrix® Elan headrest support

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning



New frame design, including adjustable wings and tall back option



Front pivot seat tilt maintains a proper position for self-propulsion



Optional dynamic rocking feature calms and creates contentment




Comfort Tension Seating® increases support as well as sitting tolerance

brodaseating.com | 1-844-552-7632

MOBILITY WITH COMFORT,
FUNCTIONALITY, & STYLE

UPHOLSTERY OPTIONS



SPECIFICATIONS

Tilt	Up to 40°
Weight Capacity	Up to 350 lbs. (114 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm) 16", 18", 20", 22" with Dynamic Rocking (40cm, 45cm, 50cm, 55cm)
Seat Depths	17", 18.5", 20" (43cm, 46cm, 50cm) 15.5", 17", 18.5", 20" with Dynamic Rocking (39cm, 43cm, 46cm, 50cm)
Back Lengths	22", 32", 33.5", 35" (55 cm, 81cm, 85cm, 89cm)
Seat Surface Heights	Adjustable 12.5" - 17.5" (31cm-44cm) Adjustable 14" - 19" with Dynamic Rocking (35-48cm)
Wheels	5" Front Casters: Total Lock (12.5cm) 20" or 22" Front Mag (50cm or 55cm) 14" Rear Mag (36cm) 5" Rear Casters (12.5cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Tall back: Neckrest, shoulder bolsters, lower side pads, seat and back pad with cover Short back: lower side pads, seat and back pad with cover

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps

brodaseating.com | 1-844-552-7632