



Front pivot seat tilt maintains user comfort with smooth, effortless tilt




Adjustable arm-rests aid in user comfort and positioning



Adjustable seat surface height allows for proper foot propulsion


Patient, Facility, & Caregiver Benefits

- *Tilt enables ease of use for frequent repositioning*
- *Design for safety and stability helping to reduce risk of injury*
- *Provides stimulation and contentment - excellent for memory care*



Comfort Tension Seating[®] improves postural support and sitting tolerance

UPHOLSTERY OPTIONS




Enhanced Mobility and Freedom

SPECIFICATIONS

Degrees of Tilt	up to 18°
Weight Capacity	350 lbs. (158 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 63cm, 66cm,)
Seat Depths	15.5", 17", 18.5", 20" (39cm, 43cm, 46cm, 50cm)
Seat Surface Heights	Adjustable: 12.5", 13.5", 14.5", 15.5", 16.5", 17.5" (31cm, 34cm, 36cm, 39cm, 41cm, 44cm)
Wheels	5" (12cm) Aluminum Non-locking Front Caster Wheels 5" (12cm) Fixed Front Caster Wheels 20" or 22" (50cm or 55cm) Front or Rear Mag Wheels 5" (12cm) Total Lock Rear Caster Wheels
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Back Support	Integrated with Comfort Tension Seating®
Standard Padding	Head Support, Adjustable Lateral Shoulder Supports, Lower Side Cushions, Seat & Back Cushions with Covers

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning belts