

## Customer Details

### Billing Address

Contact Name \_\_\_\_\_ Facility: \_\_\_\_\_

Contact Phone: \_\_\_\_\_ Contact Email: \_\_\_\_\_

Street Address: \_\_\_\_\_

### Shipping Address

Contact Name \_\_\_\_\_ Facility: \_\_\_\_\_

Contact Phone: \_\_\_\_\_ Contact Email: \_\_\_\_\_

Street Address: \_\_\_\_\_

### Order Information

Date: \_\_\_\_\_

PO Number: \_\_\_\_\_

Customer Number: \_\_\_\_\_

Tag/Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_



### Standard Features

- Comfort Tension Seating®
- 250lb Maximum Occupant Weight
- Tubular Steel, Powder Coated Frame
- 2 x 20" Rear Mag Wheels
- Adjustable Seat Surface Height
- 2 x 5" Front Total Lock Caster Wheels
- Height Adjustable Arm Supports
- Up to 25° Infinitely Adjustable Seat Tilt
- Angle Adjustable, Articulating Lower Leg Support
- Up to 45° Infinitely Adjustable Back Recline

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

## Product Options

### Chair Model

- |   |                   |   |                   |
|---|-------------------|---|-------------------|
| <input type="checkbox"/> Centric Positioning Wheelchair | <b>\$3,475.00</b> | <input type="checkbox"/> Centric Positioning Wheelchair with Mag Wheels | <b>\$4,138.00</b> |
|---|-------------------|---|-------------------|

### Seat Width

- |   |                   |   |                   |
|---|-------------------|---|-------------------|
| <input type="checkbox"/> 16" Seat Width | <b>\$0.00</b>     | <input type="checkbox"/> 18" Seat Width | <b>\$0.00</b>     |
| <input type="checkbox"/> 20" Seat Width | <b>\$0.00</b>     | <input type="checkbox"/> 22" Seat Width | <b>\$500.00</b>   |
| <input type="checkbox"/> 24" Seat Width | <b>\$1,131.00</b> | <input type="checkbox"/> 26" Seat Width | <b>\$1,131.00</b> |
| <input type="checkbox"/> 28" Seat Width | <b>\$1,131.00</b> |   |                   |

### Seat Depth

- |   |                 |   |                 |
|---|-----------------|---|-----------------|
| <input type="checkbox"/> 15.5" Seat Depth (Lower Leg Support Assembly Option N/A)   | <b>\$121.00</b> | <input type="checkbox"/> 17" Seat Depth (Check Lower Leg Support Assembly Option) | <b>\$0.00</b>   |
| <input type="checkbox"/> 18.5" Seat Depth (Check Lower Leg Support Assembly Option) | <b>\$0.00</b>   | <input type="checkbox"/> 20" Seat Depth (Check Lower Leg Support Assembly Option) | <b>\$121.00</b> |
| <input type="checkbox"/> Custom Seat Depth (Specify)                                | <b>\$0.00</b>   |   |                 |

### Seat Height

- |  |               |   |               |
|--|---------------|---|---------------|
| <input type="checkbox"/> 19.5" Seat Height | <b>\$0.00</b> | <input type="checkbox"/> 20.5" Seat Height            | <b>\$0.00</b> |
| <input type="checkbox"/> 21.5" Seat Height | <b>\$0.00</b> | <input type="checkbox"/> Custom Seat Height (Specify) | <b>\$0.00</b> |

### Back Height

- |  |                 |   |                 |
|--|-----------------|---|-----------------|
| <input type="checkbox"/> 32" Back Height | <b>\$0.00</b>   | <input type="checkbox"/> 33.5" Back Height            | <b>\$121.00</b> |
| <input type="checkbox"/> 35" Back Height | <b>\$121.00</b> | <input type="checkbox"/> Custom Back Height (Specify) | <b>\$0.00</b>   |

## Upholstery Color

<input type="checkbox"/> Renaissance Burgundy Vinyl	\$0.00	<input type="checkbox"/> Renaissance Charcoal Vinyl	\$0.00
<input type="checkbox"/> Silvertex Sapphire Vinyl	\$0.00		

## Strap Color

<input type="checkbox"/> Grey Strap	\$0.00	<input type="checkbox"/> Blue Strap	\$0.00
<input type="checkbox"/> Green Strap	\$0.00		

## Wheel Options

<input type="checkbox"/> 20" Front Mag Wheels	\$0.00	<input type="checkbox"/> 20" Rear Mag Wheels	\$0.00
<input type="checkbox"/> 22" Front Mag Wheels	\$0.00	<input type="checkbox"/> 22" Rear Mag Wheels	\$0.00
<input type="checkbox"/> 24" Front Mag Wheels	\$0.00	<input type="checkbox"/> 24" Rear Mag Wheels	\$0.00
<input type="checkbox"/> 1 Rear Directional Lock Caster	\$0.00		

## Wheel Lock Extensions

<input type="checkbox"/> Parking Brake Extension Handles (Pair) (E0961)	\$82.00
---	---------

## Comfort Co. Lateral Support

<input type="checkbox"/> Comfort Co. Lateral Support (E0956) Compatible with seat widths up to 24"	\$564.00
---	----------

## Arm Pad

<input type="checkbox"/> Narrow Arm Supports	\$0.00	<input type="checkbox"/> Cushioned Arm Support (Vinyl) Not compatible with a tray	\$95.00
<input type="checkbox"/> Soft Arm Supports	\$0.00	<input type="checkbox"/> Arm Trough (Pair)	\$329.00
<input type="checkbox"/> Arm Trough (L) & Soft Arm Support (R)	\$165.00	<input type="checkbox"/> Arm Trough (R) & Soft Arm Support (L)	\$165.00

## Legrest Type

<input type="checkbox"/> Strapped Lower Leg Support	\$0.00	<input type="checkbox"/> ABS Lower Leg Support	\$411.00
<input type="checkbox"/> Contracture Sling Frame In place of lower leg support	\$316.00	<input type="checkbox"/> Flared Lower Leg Support Assembly Long (Pin & Saddle)	\$139.00
<input type="checkbox"/> Flared Lower Leg Support Assembly Short (Pin & Saddle)	\$139.00	<input type="checkbox"/> Straight Lower Leg Support Assembly Long (Pin & Saddle)	\$139.00
<input type="checkbox"/> Straight Lower Leg Support Assembly Short (Pin & Saddle)	\$139.00		

## Strapped Footrest

<input type="checkbox"/> Single Lower Leg Support 2 Strap (12"-19")	\$0.00	<input type="checkbox"/> Single Lower Leg Support 3 Strap (14"-21")	\$0.00
<input type="checkbox"/> Single Lower Leg Support 4 Strap (16"-23")	\$0.00	<input type="checkbox"/> Single Lower Leg Support 5 Strap (18"-25")	\$0.00
<input type="checkbox"/> Split Lower Leg Support 2 Strap (12"-19") Not available in 16" seat width	\$0.00	<input type="checkbox"/> Split Lower Leg Support 3 Strap (14"-21") Not available in 16" seat width	\$0.00
<input type="checkbox"/> Split Lower Leg Support 4 Strap (16"-23") Not available in 16" seat width	\$0.00	<input type="checkbox"/> Split Lower Leg Support 5 Strap (18"-25") Not available in 16" seat width	\$0.00

## Swing Away Footrest

<input type="checkbox"/> Angle-Adjustable Swing-Away Lower Leg Support [Style 07/08]	\$310.00	<input type="checkbox"/> Swing-Away Lower Leg Support [Style 03/04]	\$247.00
--	----------	---	----------

## Neckrest Pads

<input type="checkbox"/> Neck Support	\$0.00	<input type="checkbox"/> Contoured Visco Neck Support (Vinyl)	\$51.00
<input type="checkbox"/> Contoured Visco Neck Support (Triptex)	\$51.00	<input type="checkbox"/> No Neck Support	\$0.00
<input type="checkbox"/> Soft Visco Neck Support Cushion	\$51.00		

## Shoulder Bolsters

<input type="checkbox"/> Shoulder Support Short	\$0.00	<input type="checkbox"/> Shoulder Support Long	\$165.00
<input type="checkbox"/> Shoulder Support ACC Short	\$0.00	<input type="checkbox"/> Shoulder Support ACC Long	\$165.00
<input type="checkbox"/> Shoulder Support Round Short	\$0.00	<input type="checkbox"/> Shoulder Support Round Long	\$165.00

## Lateral Covers

<input type="checkbox"/> Short Adjustable Shoulder Support Covers (HSP)	\$209.00	<input type="checkbox"/> Long Adjustable Shoulder Support Covers (HSP)	\$272.00
---	----------	--	----------

## Seat &amp; Back Pads

<input type="checkbox"/> Grey Terry	\$0.00	<input type="checkbox"/> Removable Vinyl Seat & Back Covers	\$0.00
<input type="checkbox"/> Full Cover Seat & Back Pads	\$121.00	<input type="checkbox"/> Full Cover Seat & Back Pads With Slot	\$121.00
<input type="checkbox"/> Removable Startex Seat & Back Covers	\$0.00		

## Footrest Pads

<input type="checkbox"/> Sole Cushion	\$0.00	<input type="checkbox"/> Neoprene Foot Sole Pad (Black)	\$0.00
<input type="checkbox"/> Neoprene Foot Sole Pad with Vinyl	\$0.00	<input type="checkbox"/> Sole Cover With Neoprene and Heel Cover for ABS Foorest Retainment	\$0.00

## Foot Box

<input type="checkbox"/> Long Footbox (HSP)	\$411.00	<input type="checkbox"/> Long Footbox With pocket	\$443.00
---	----------	---	----------

**Pelvic Positioning Belt**

<input type="checkbox"/> Pelvic Positioning Belt Side Buckle	<b>\$64.00</b>	<input type="checkbox"/> Pelvic Positioning Belt Center Buckle	<b>\$64.00</b>
<input type="checkbox"/> Pelvic Positioning Belt Side Auto-Style Buckle	<b>\$64.00</b>	<input type="checkbox"/> Pelvic Positioning Belt Center Auto-Style Buckle	<b>\$64.00</b>
<input type="checkbox"/> Padded Pelvic Positioning Belt (Black)	<b>\$91.00</b>	<input type="checkbox"/> Padded Pelvic Positioning Belt (Matches Selected Accent Color)	<b>\$91.00</b>
<input type="checkbox"/> Padded Pelvic Positioning Belt with Auto-Style Buckle (Black)	<b>\$91.00</b>	<input type="checkbox"/> Padded Pelvic Positioning Belt with Auto-Style Buckle (Matches Selected Accent Color)	<b>\$91.00</b>

**Thigh Positioning Belt**

<input type="checkbox"/> Cushioned Thigh Positioning Belt (HSP)	<b>\$146.00</b>	<input type="checkbox"/> Cushioned Thigh Positioning Belt with Auto-Style Buckle	<b>\$152.00</b>
<input type="checkbox"/> Neoprene Thigh Positioning Belt	<b>\$146.00</b>	<input type="checkbox"/> Neoprene Thigh Positioning Belt with Auto-Style Buckle	<b>\$152.00</b>
<input type="checkbox"/> APP Extended Neoprene Thigh Positioning Belt	<b>\$146.00</b>	<input type="checkbox"/> APP Extended Neoprene Thigh Positioning Belt with Auto-Style Buckle	<b>\$152.00</b>

**Side Reduction Pad**

<input type="checkbox"/> Side Reduction Pad (Black)	<b>\$121.00</b>	<input type="checkbox"/> Side Reduction Pad (Matches Selected Accent Color)	<b>\$121.00</b>
---	-----------------	---	-----------------

**Side Arm Support**

<input type="checkbox"/> Side Arm Support (Black)	<b>\$158.00</b>	<input type="checkbox"/> Side Arm Support (Matches Selected Accent Color)	<b>\$158.00</b>
---	-----------------	---	-----------------

**Front Stabilizers**

<input type="checkbox"/> Front Stabilizers	<b>\$108.00</b>
--	-----------------

**Rear Stabilizers**

<input type="checkbox"/> Rear Stabilizers	<b>\$108.00</b>
---	-----------------

**Lateral Stabilizers**

<input type="checkbox"/> Lateral Stabilizers	<b>\$190.00</b>
--	-----------------

**Lateral Torso Supports**

Lateral Torso Supports **\$695.00**

**Accessories****Accessories - Add-On**

Oxygen Carry Bag **\$109.00**  Cup Holder (Grey Swivel) **\$57.00**

**Acute Care Configuration (ACC)****Acute Care Configuration (ACC) Package - Chair Model**

Includes: Full Vinyl Covers, ABS Lower Leg Support with Foot Support Cushions (Factory installed, 16"-24" widths)

ACC Package - For Centric/30VT **\$552.00**

Includes: Full Vinyl Covers, Acute Care Short Shoulder Supports, ABS Lower Leg Support with Foot Support Cushions. Factory installed. Skip lower leg support cushion package section when selecting this option. 16"-24" widths.

**Additional Positioning Padding (APP)-Slim Posture Padding (SPP)****Additional Positioning Padding (APP) - Package Price**

Broda APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Covers (Not compatible with 35" back support height)

APP Package - Positioning - 16" 18" 20" Widths **\$641.00**  APP Package - Positioning - 22" 24" Widths **\$771.00**

APP Package - Positioning - 26" 28" Widths **\$901.00**

## Slim Posture Padding (SPP) - Seat Cushion

<input type="checkbox"/> SPP Seat Cushion 16" Seat Width X 17" Seat Depth (E2601) <b>\$212.00</b>	<input type="checkbox"/> SPP Seat Cushion 16" Seat Width X 18.5" Seat Depth (E2601) <b>\$212.00</b>
<input type="checkbox"/> SPP Seat Cushion 16" Seat Width X 20" Seat Depth (E2601) <b>\$212.00</b>	<input type="checkbox"/> SPP Seat Cushion 18" Seat Width X 17" Seat Depth (E2601) <b>\$212.00</b>
<input type="checkbox"/> SPP Seat Cushion 18" Seat Width X 18.5" Seat Depth (E2601) <b>\$212.00</b>	<input type="checkbox"/> SPP Seat Cushion 18" Seat Width X 20" Seat Depth (E2601) <b>\$212.00</b>
<input type="checkbox"/> SPP Seat Cushion 20" Seat Width X 17" Seat Depth (E2601) <b>\$212.00</b>	<input type="checkbox"/> SPP Seat Cushion 20" Seat Width X 18.5" Seat Depth (E2601) <b>\$212.00</b>
<input type="checkbox"/> SPP Seat Cushion 20" Seat Width X 20" Seat Depth (E2601) <b>\$212.00</b>	<input type="checkbox"/> SPP Seat Cushion 22" Seat Width X 17" Seat Depth (E2602) <b>\$234.00</b>
<input type="checkbox"/> SPP Seat Cushion 22" Seat Width X 18.5" Seat Depth (E2602) <b>\$234.00</b>	<input type="checkbox"/> SPP Seat Cushion 22" Seat Width X 20" Seat Depth (E2602) <b>\$234.00</b>
<input type="checkbox"/> SPP Seat Cushion 24" Seat Width X 17" Seat Depth (E2602) <b>\$234.00</b>	<input type="checkbox"/> SPP Seat Cushion 24" Seat Width X 18.5" Seat Depth (E2602) <b>\$234.00</b>
<input type="checkbox"/> SPP Seat Cushion 24" Seat Width X 20" Seat Depth (E2602) <b>\$234.00</b>	<input type="checkbox"/> SPP Seat Cushion 26" Seat Width X 17" Seat Depth (E2602) <b>\$262.00</b>
<input type="checkbox"/> SPP Seat Cushion 26" Seat Width X 18.5" Seat Depth (E2602) <b>\$262.00</b>	<input type="checkbox"/> SPP Seat Cushion 26" Seat Width X 20" Seat Depth (E2602) <b>\$262.00</b>

## Accessory Bag

## Accessory Bag

<input type="checkbox"/> Accessory Bag (Single) <b>\$51.00</b>	<input type="checkbox"/> Accessory Bag (Pair) <b>\$102.00</b>
<input type="checkbox"/> Embroidered Accessory Bag (Single) <b>\$83.00</b>	<input type="checkbox"/> Embroidered Accessory Bag (Pair) <b>\$165.00</b>

## Accessory Bag - Embroidery

<input type="checkbox"/> Accessory Bag Embroidery - Block Font <b>\$0.00</b>	<input type="checkbox"/> Accessory Bag Embroidery - Brush Font <b>\$0.00</b>
--	--

## Accessory Bag - Characters

- |   |               |   |               |
|---|---------------|---|---------------|
| <input type="checkbox"/> Accessory Bag - 2 Characters | <b>\$0.00</b> | <input type="checkbox"/> Accessory Bag - 1 Character Only | <b>\$0.00</b> |
|---|---------------|---|---------------|

## Accessory Bag - 1st Character

## Accessory Bag - 2nd Character

## Huntington's Special Padding (HSP)

## Huntington's Special Padding (HSP) Package - Chair Model

- |  |                 |
|--|-----------------|
| <input type="checkbox"/> HSP Package - For 30  | <b>\$651.00</b> |
| Includes Long Lateral Shoulder Supports, Thigh Belt & Long Lower Leg Support Footbox |                 |

## Huntington's Special Padding (HSP) Package - Shoulder Bolster

- |  |               |  |               |
|--|---------------|--|---------------|
| <input type="checkbox"/> HSP Package - Shoulder Bolster Long (Matches Selected Upholstery Color) (E0956)<br>Available in all upholstery colors | <b>\$0.00</b> | <input type="checkbox"/> HSP Package - Shoulder Bolster ACC Long (Matches Selected Upholstery Color) (E0956)<br>Available in all upholstery colors | <b>\$0.00</b> |
| <input type="checkbox"/> HSP Package - Shoulder Bolster Standard + Wing Cover Long (Matches Selected Upholstery Color) (E0956)                 | <b>\$0.00</b> | <input type="checkbox"/> HSP Package - Wing Cover Long (Matches Selected Upholstery Color) (E0956)   | <b>\$0.00</b> |

## Huntington's Special Padding (HSP) Package - Thigh Positioning Belt

- |  |               |  |               |
|--|---------------|--|---------------|
| <input type="checkbox"/> HSP Package - Padded Thigh Positioning Belt (Matches Selected Upholstery Color) (Standard) (E0957)<br>Available in all upholstery colors except Beeline Espresso and Beeline Portobello | <b>\$0.00</b> | <input type="checkbox"/> HSP Package - Padded Thigh Positioning Belt with Auto-Style Buckle (Matches Selected Upholstery Color) (E0957)<br>Available in all upholstery colors except Beeline Espresso and Beeline Portobello | <b>\$0.00</b> |
| <input type="checkbox"/> HSP Package - Neoprene Thigh Positioning Belt (E0957)   | <b>\$0.00</b> | <input type="checkbox"/> HSP Package - Neoprene Thigh Positioning Belt with Auto-Style Buckle (E0957)  | <b>\$0.00</b> |
| <input type="checkbox"/> HSP Package - Neoprene Extended Thigh Positioning Belt (Seat Pan or FVC Without Slot) APP (E0957)   | <b>\$0.00</b> | <input type="checkbox"/> HSP Package - Neoprene Extended Thigh Positioning Belt with Auto-Style Buckle (Seat Pan or FVC Without Slot) (E0957)  | <b>\$0.00</b> |

**IV Pole****Add IV-Pole**

IV-Pole **\$121.00**

**IV Pole - Factory Installed on Chair**

IV Pole - For 30 Without Mag Wheels **\$0.00**  
If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

IV Pole - For 30 With Front Mag Wheels **\$0.00**  
If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

IV Pole - For 30 With Rear Mag Wheels **\$0.00**

**IV Pole - Frame Color**

IV Pole - Grey Standard Frame **\$0.00**

**IV Pole - Mount**

IV Pole - Standard IV-Pole Bracket Only **\$0.00**

**IV Pole - Location**

IV Pole - Front Mount **\$0.00**

IV Pole - Rear Mount **\$0.00**

**O2 Holder****Add Oxygen Holder**

Oxygen Holder **\$203.00**

**O2 Holder - Factory Installed on Chair**

Oxygen Holder - For 30 Without Mag Wheels **\$0.00**  
If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

Oxygen Holder - For 30 With Front Mag Wheels **\$0.00**  
If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

Oxygen Holder - For 30 With Rear Mag Wheels **\$0.00**

**O2 Holder - Frame Color**

Oxygen Holder - Grey Frame **\$0.00**

**O2 Holder - Mounting Type**

Oxygen Holder - Standard Oxygen Holder Clamps Only **\$0.00**

If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

**O2 Holder - Mounting Location**

Oxygen Holder - Front Mount **\$0.00**

If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

Oxygen Holder - Rear Mount **\$0.00**

**Lap Tray****Add Lap Tray**

Lap Tray **\$0.00**

**Lap Tray - Width**

Full Tray - 16" Seat Width (E0950) **\$373.00**

Full Tray - 18" Seat Width (E0950) **\$373.00**

Full Tray - 20" Seat Width (E0950) **\$373.00**

Full Tray - 22" Seat Width (E0950) **\$373.00**

Full Tray - 24" Seat Width (E0950) **\$537.00**

Full Tray - 26" Seat Width (E0950) **\$537.00**

Full Tray - 28" Seat Width (E0950) **\$537.00**

**Lap Tray - Top Color**

Full Tray (Black) **\$0.00**

Full Tray (Gray) **\$0.00**

Full Tray (Clear) **\$13.00**

**Lap Tray - Upholstery Color**

Full Tray COVER - Modena Black Vinyl **\$190.00**

Full Tray COVER - Renaissance Burgundy Vinyl **\$190.00**

Full Tray COVER - Renaissance Charcoal Vinyl **\$190.00**

Full Tray COVER - Silvertex Sapphire Vinyl **\$190.00**

**Lap Tray - Retaining Strap**

- |   |        |  |        |
|---|--------|--|--------|
| <input type="checkbox"/> Full Tray - Standard Strap - Center Buckle     | \$0.00 | <input type="checkbox"/> Full Tray - Standard Strap - Left-Hand Buckle | \$0.00 |
| <input type="checkbox"/> Full Tray - Standard Strap - Right-Hand Buckle | \$0.00 | <input type="checkbox"/> Full Tray - Sew-On Strap - Left-Hand Buckle   | \$0.00 |
| <input type="checkbox"/> Full Tray - Sew-On Strap - Right-Hand Buckle   | \$0.00 |  |        |

**Veteran****Military Veteran Order**

- |  |        |
|--|--------|
| <input type="checkbox"/> VA Order "Honor Our Veterans" | \$0.00 |
|--|--------|

**Branch of Service Crest**

- |   |        |  |        |
|---|--------|--|--------|
| <input type="checkbox"/> Crest U.S. Air Force   | \$0.00 | <input type="checkbox"/> Crest U.S. Army         | \$0.00 |
| <input type="checkbox"/> Crest U.S. Coast Guard | \$0.00 | <input type="checkbox"/> Crest U.S. Marine Corps | \$0.00 |
| <input type="checkbox"/> Crest U.S. Navy        | \$0.00 |  |        |

## Shipping &amp; Freight

## Initial Shipping Costs

- |  |                 |  |                 |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Ship to Authorized Broda Dealer | <b>\$165.00</b> | <input type="checkbox"/> Ship Direct to Customer | <b>\$200.00</b> |
|--|-----------------|--|-----------------|

## Additional Delivery Services

- |  |                 |   |                         |
|--|-----------------|---|-------------------------|
| <input type="checkbox"/> Lift Gate Services<br>Delivery is made with a lift gate trailer that moves from trailer level to ground level.              | <b>\$200.00</b> | <input type="checkbox"/> Residential Delivery Service<br>Delivery is to a residential location.   | <b>\$200.00</b>         |
| <input type="checkbox"/> Notification Before Delivery<br>Contact is requested to the cosignee by the carrier PRIOR to shipment delivery.             | <b>\$85.00</b>  | <input type="checkbox"/> Inside Delivery Service<br>Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | <b>Min<br/>\$250.00</b> |
| <input type="checkbox"/> Limited Access Delivery Service<br>Delivery is made with a lift gate trailer that moves from trailer level to ground level. | <b>\$165.00</b> |   |                         |