

Customer Details

Billing Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Shipping Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Order Information

Date: _____

PO Number: _____

Customer Number: _____

Tag/Comments: _____



Standard Features

- Tubular Steel, Powder Coated Frame
- Drop In, Swing Away Lower Leg Supports
- 45 Degrees, Infinitely Adjustable Seat Tilt
- Composite Foot Support Plates
- Fixed Height Arm Supports
- Cushioned Side Panels
- Adjustable Seat Depth 16" - 20"
- Adjustable Backrest Angle
- Adjustable Seat Surface Height 14" - 19"
- 250 lbs. Maximum Occupant Weight
- Choice of Wheel, Caster and Fork Sizes
- Seat, back and head support cushions not included
- Push Lock Parking Brakes

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

Product Options

Chair Model

Comfort Rehab Wheelchair (E1161) **\$3,045.00**

Seat Type

Solid Seat Pan (Required) (E2231) **\$167.00**

Seat Width

16" Seat Width **\$0.00** 18" Seat Width **\$0.00**

20" Seat Width (E2201) **\$217.00** 22" Seat Width (E2201) **\$217.00**

Seat Depth

16" Seat Depth **\$0.00** 17" Seat Depth **\$0.00**

18" Seat Depth **\$0.00** 19" Seat Depth **\$0.00**

20" Seat Depth (E2203) **\$0.00**

Seat Height

14" Seat Height **\$0.00** 15" Seat Height **\$0.00**

15.5" Seat Height **\$0.00** 16" Seat Height **\$0.00**

16.5" Seat Height **\$0.00** 17" Seat Height **\$0.00**

17.5" Seat Height **\$0.00** 18" Seat Height **\$0.00**

19" Seat Height **\$0.00**

Back Rail Height

<input type="checkbox"/> 16" Folding Back Rail Height	\$224.00	<input type="checkbox"/> 16" Non-Folding Back Rail Height	\$0.00
<input type="checkbox"/> 18" Folding Back Rail Height	\$224.00	<input type="checkbox"/> 21" Folding Back Rail Height	\$224.00
<input type="checkbox"/> 24" Folding Back Rail Height	\$224.00		

Back Angle

<input type="checkbox"/> 10° Recline Angle	\$0.00	<input type="checkbox"/> 15° Recline Angle	\$0.00
<input type="checkbox"/> 20° Recline Angle	\$0.00	<input type="checkbox"/> 5° Recline Angle	\$0.00

Base Frame Color

<input type="checkbox"/> Black Base Frame	\$0.00	<input type="checkbox"/> Bronze Base Frame	\$0.00
<input type="checkbox"/> Carmine Red Base Frame	\$0.00	<input type="checkbox"/> Deep Orange Base Frame	\$0.00
<input type="checkbox"/> Forest Green Base Frame	\$0.00	<input type="checkbox"/> Ocean Blue Base Frame	\$0.00
<input type="checkbox"/> Passion Pink Base Frame	\$0.00	<input type="checkbox"/> Royal Purple Base Frame	\$0.00
<input type="checkbox"/> Satin Nickel Base Frame	\$0.00	<input type="checkbox"/> Tungsten Sand Base Frame	\$0.00
<input type="checkbox"/> Zinc Yellow Base Frame	\$0.00		

Wheel Options

- | | |
|---|---|
| <input type="checkbox"/> 5" Medium Fork 20" Mag Wheels MP #2 CF #1 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 5" Medium Fork 20" Mag Wheels MP #3 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 5" Medium Fork 20" Mag Wheels MP #4 CF #5 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 5" Medium Fork 22" Mag Wheels MP #1 CF #1 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 5" Medium Fork 22" Mag Wheels MP #2 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 5" Medium Fork 22" Mag Wheels MP #3 CF #5 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 5" Medium Fork 24" Mag Wheels MP #1 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 5" Medium Fork 24" Mag Wheels MP #2 CF #4 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 5" Small Fork 20" Mag Wheels MP #1 CF #1 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights | <input type="checkbox"/> 5" Small Fork 20" Mag Wheels MP #2 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |
| <input type="checkbox"/> 5" Small Fork 20" Mag Wheels MP #3 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights | <input type="checkbox"/> 5" Small Fork 22" Mag Wheels MP #1 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |
| <input type="checkbox"/> 5" Small Fork 22" Mag Wheels MP #2 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights | <input type="checkbox"/> 5" Small Fork 24" Mag Wheels MP #1 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |
| <input type="checkbox"/> 6" Medium Fork 20" Mag MP #2 CF #1 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 6" Medium Fork 20" Mag MP #3 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 6" Medium Fork 20" Mag MP #4 CF #4 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 6" Medium Fork 22" Mag MP #1 CF #1 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |

- | | |
|--|--|
| <input type="checkbox"/> 6" Medium Fork 22" Mag MP #2 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 6" Medium Fork 22" Mag MP #3 CF #4 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 6" Medium Fork 24" Mag MP #1 CF #1 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 6" Medium Fork 24" Mag MP #2 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 6" Medium Fork 24" Mag MP #3 CF #5 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 6" Small Fork 20" Mag MP #3 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |
| <input type="checkbox"/> 6" Small Fork 22" Mag MP #2 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights | <input type="checkbox"/> 6" Small Fork 24" Mag MP #1 CF #2 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |
| <input type="checkbox"/> 6" Small Fork 24" Mag MP #2 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights | <input type="checkbox"/> 8" Medium Fork 20" Mag MP #5 CF #4 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 8" Medium Fork 22" Mag MP #4 CF #4 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 8" Medium Fork 24" Mag MP #3 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights |
| <input type="checkbox"/> 8" Medium Fork 24" Mag MP #4 CF #5 With Adjustable Anti-Tip Devices \$0.00
Recommended for large mag wheel (24") and higher seat surface heights | <input type="checkbox"/> 8" Small Fork 20" Mag MP #4 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |
| <input type="checkbox"/> 8" Small Fork 22" Mag MP #3 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights | <input type="checkbox"/> 8" Small Fork 24" Mag MP #3 CF #3 With Adjustable Anti-Tip Devices \$0.00
Recommended for small mag wheel (20" or 22") and lower seat surface heights |

Mag Wheel Location

- Mag Wheel Position "A" (18.5" from Caster Swivel Axis) **\$0.00**
Most forward- smaller turning radius
- Mag Wheel Position "B" (19.5" from Caster Swivel Axis) **\$0.00**
- Mag Wheel Position "C" (20.5" from Caster Swivel Axis) **\$0.00**
- Mag Wheel Position "D" (21.5" from Caster Swivel Axis) **\$0.00**
Most rearward- larger turning radius. Not recommended with 20" mag wheels

Wheel Lock Extensions

- With Parking Brake Extensions (Pair) (E0961) **\$73.00**

Anti-Tippers

- Rear Anti-Tip Devices (Required) (Pair) (E0971) **\$150.00**

Arm Bottom

- Adjustable Height Arm Supports (Pair) (E0973) **\$278.00**
- Fixed Height Arm Supports (Pair) **\$0.00**

Swing Away Footrest

- Elevating Lower Leg Support with Angle-Adjustable Foot Plate (Drop-In) (Style 18/19) (Pair) (E0990 + K0040) **\$667.00**
- Lower Leg Support 60° angle adj foot plate (12/13) (Pair) (K0040) **\$0.00**
- Lower Leg Support 75° angle adj foot plate (10/11) (Pair) (K0040) **\$0.00**
- Lower Leg Support 75° fixed foot plate (14/15) (Pair) **\$0.00**
- No Lower Leg Support **\$0.00**
- Swing-Away Lower Leg Support Small Angle-Adjustable Foot Plate With Heel Loops (Drop-In) [Style 36/ (E0951 + K0040)] **\$239.00**

Footrest Extension

- Ikon Extended Footrest Adjustment Tube (Pair) **\$69.00**
Only compatible on the F-SW-IKN and F-SW-IKN-ADJ

Footrest Heel Loop

- With Heel Loop (E0951) **\$22.00**

Footrest Accessories

- Cushion Lower Leg Support Panel Black (E0954) **\$73.00**
For calf support

Footrest Hanger / Front Rigging

- Offset Lower Leg Support Assembly (Drop-In) (Pair) **\$0.00**
- Offset Lower Leg Support Assembly (Pin & Saddle) (Pair) **\$0.00**
- Straight Lower Leg Support Assembly (Drop-In) (Pair) **\$0.00**
- Straight Lower Leg Support Assembly (IKN Drop-In) (Pair) **\$0.00**
- Straight Lower Leg Support Assembly (Pin & Saddle) (Pair) **\$0.00**

Pelvic Positioning Belt

- 587 Pelvic Positioning Belt Center Buckle (E0978) **\$64.00**

Side Reduction Pad

- Side Reduction Pad (Black) (Pair) (E0956) **\$123.00**

Accessories**Accessories - Add-On**

- Cup Holder (Grey Swivel) **\$57.00**
- Oxygen Carry Bag (E2208) **\$98.00**

IV Pole**Add IV-Pole**

- IV-Pole **\$115.00**

IV Pole - Factory Installed on Chair

- IV Pole - For 587/Comfort Rehab Wheel Chair **\$0.00**
If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

IV Pole - Frame Color

IV Pole - Black Standard Frame **\$0.00**

IV Pole - Mount

IV Pole - Standard IV-Pole Bracket Only **\$0.00**

IV Pole - Location

IV Pole - Front Mount **\$0.00**

Matrx**Motion Concepts (MC) - Seat Model**

Matrx Libra Seat Cushion 16" 18" 20" Seat Width (E2624) **\$485.00**
Recommended for highest risk of or existing skin breakdown

Matrx Libra Seat Cushion 22" 24" Seat Width (E2625) **\$595.00**
Recommended for highest risk of or existing skin breakdown

Matrx Multi Seat Cushion 16" 17" 18" 20" Seat Width (E2607) **\$445.00**
Recommended for moderate to high risk of skin breakdown

Matrx Multi Seat Cushion 22" Seat Width (E2608) **\$595.00**
Recommended for moderate to high risk of skin breakdown

Matrx PS Positioning Seat Cushion 16" 18" 20" Seat Width (E2605) **\$345.00**
Recommended for low to moderate risk of skin breakdown

Matrx PS Positioning Seat Cushion 22" 24" 26" Seat Width (E2606) **\$545.00**
Recommended for low to moderate risk of skin breakdown

Matrx PSP Positioning Seat Cushion 16" 18" 20" Seat Width (E2607) **\$395.00**
Recommended for moderate to high risk of skin breakdown

Matrx Vi Seat Cushion 16" 18" 20" Seat Width (E2607) **\$425.00**
Recommended for moderate to high risk of skin breakdown

Matrx Vi Seat Cushion 22" Seat Width (E2608) **\$475.00**
Recommended for moderate to high risk of skin breakdown

Motion Concepts (MC) - Seat Width

Matrx Seat Cushion 16" Seat Width **\$0.00**

Matrx Seat Cushion 18" Seat Width **\$0.00**

Matrx Seat Cushion 20" Seat Width **\$0.00**

Matrx Seat Cushion 22" Seat Width **\$0.00**

Motion Concepts (MC) - Seat Depth

<input type="checkbox"/> Matrx Seat Cushion 16" Deep	\$0.00	<input type="checkbox"/> Matrx Seat Cushion 18" Deep	\$0.00
<input type="checkbox"/> Matrx Seat Cushion 20" Deep	\$0.00	<input type="checkbox"/> Matrx Seat Cushion 22" Deep Only available w/ Multi & PS Cushions	\$0.00

Motion Concepts (MC) - Seat Rigidizer

<input type="checkbox"/> Matrx Seat Rigidizer 16 Large (20" 22" Seat Depth) (E0992)	\$115.00	<input type="checkbox"/> Matrx Seat Rigidizer 16 Medium (18" Seat Depth) (E0992)	\$115.00
<input type="checkbox"/> Matrx Seat Rigidizer 16 Small (16" Seat Depth) (E0992)	\$115.00	<input type="checkbox"/> Matrx Seat Rigidizer 18 Large (20" 22" Seat Depth) (E0992)	\$115.00
<input type="checkbox"/> Matrx Seat Rigidizer 18 Medium (18" Seat Depth) (E0992)	\$115.00	<input type="checkbox"/> Matrx Seat Rigidizer 18 Small (16" Seat Depth) (E0992)	\$115.00
<input type="checkbox"/> Matrx Seat Rigidizer 20 Large (20" 22" Seat Depth) (E0992)	\$115.00	<input type="checkbox"/> Matrx Seat Rigidizer 20 Medium (18" Seat Depth) (E0992)	\$115.00
<input type="checkbox"/> Matrx Seat Rigidizer 20 Small (16" Seat Depth) (E0992)	\$115.00	<input type="checkbox"/> Matrx Seat Rigidizer 22 Large (20" 22" Seat Depth) (E0992)	\$115.00
<input type="checkbox"/> Matrx Seat Rigidizer 22 Medium (18" Seat Depth) (E0992)	\$115.00		

Motion Concepts (MC) - Seat Accessories

<input type="checkbox"/> Matrx 1" Anterior Wedge	\$129.00	<input type="checkbox"/> Matrx 1" Hemi Wedge	\$67.00
<input type="checkbox"/> Matrx Libra Abductor	\$62.00	<input type="checkbox"/> Matrx Libra Fluid Supplement	\$78.00
<input type="checkbox"/> Matrx Pelvic Obliquity Kit	\$52.00	<input type="checkbox"/> Matrx Seat Slip Cover	\$67.00
<input type="checkbox"/> Matrx Tapered Wedge (Left)	\$52.00	<input type="checkbox"/> Matrx Tapered Wedge (Right)	\$52.00

Motion Concepts (MC) - Back Model

<input type="checkbox"/> Matrx E2 Back - 16" 17" 18" 19" 20" Seat Width (E2615) \$645.00 Recommended for users seeking optimal immersion and airflow	<input type="checkbox"/> Matrx E2 Deep Back - 16" 17" 18" 19" 20" Seat Width \$775.00 (E2620) Recommended for users seeking optimal immersion w/ contoured support needs
<input type="checkbox"/> Matrx Elite Back - 16" 18" 20" Seat Width (E2620) \$625.00 Recommended for more active users	<input type="checkbox"/> Matrx Elite Back - 22" 24" Seat Width (E2621) \$795.00 Recommended for more active users
<input type="checkbox"/> Matrx Elite Deep Back - 16" 18" 20" Seat Width (E2620) \$760.00 Recommended for more active users with contoured support needs	<input type="checkbox"/> Matrx Elite Deep Back - 22" 24" Seat Width (E2621) \$930.00 Recommended for more active users with contoured support needs
<input type="checkbox"/> Matrx Elite TR Back - 16" 18" 20" Seat Width (E2620) \$650.00 Recommended for tilt in space & recline applications	<input type="checkbox"/> Matrx Elite TR Back - 22" 24" Seat Width (E2621) \$845.00 Recommended for tilt in space & recline applications
<input type="checkbox"/> Matrx MAC Foam Back - 16" 18" 20"SW (E2620) \$955.00	<input type="checkbox"/> Matrx MAC Hug Back - 16" 18" 20"SW (E2620) \$1,055.00
<input type="checkbox"/> Matrx Posture Back (PB) - 16" 18" 20" Seat Width (E2615) \$625.00 Recommended for less active users seeking comfort	<input type="checkbox"/> Matrx Posture Back (PB) Deep - 16" 18" 20" Seat Width (E2615) \$645.00 Recommended for less active users seeking comfort & support

Motion Concepts (MC) - Back Width

<input type="checkbox"/> Matrx Back Width - 16" \$0.00	<input type="checkbox"/> Matrx Back Width - 18" \$0.00
<input type="checkbox"/> Matrx Back Width - 20" \$0.00	<input type="checkbox"/> Matrx Back Width - 22" \$0.00

Motion Concepts (MC) - Back Height

<input type="checkbox"/> Matrx Back Height 16" \$0.00	<input type="checkbox"/> Matrx Back Height 18" \$0.00 Excludes PB & PB Deep
<input type="checkbox"/> Matrx Back Height 20" \$0.00	<input type="checkbox"/> Matrx Back Height 22" \$0.00 MAC Back only

Motion Concepts (MC) - Back MAC Lateral LH

- | | | | |
|---|---------------|---|---------------|
| <input type="checkbox"/> Left-Side Lateral Support for Matrx MAC Back - 3.5" Contour Depth (Each) | \$0.00 | <input type="checkbox"/> Left-Side Lateral Support for Matrx MAC Back - 6.0" Contour Depth (Each) | \$0.00 |
| <input type="checkbox"/> Left-Side Lateral Support for Matrx MAC Back - 9.0" Contour Depth (Each) | \$0.00 | | |

Motion Concepts (MC) - Back MAC Lateral RH

- | | | | |
|--|---------------|--|---------------|
| <input type="checkbox"/> Right-Side Lateral Support for Matrx MAC Back - 3.5" Contour Depth (Each) | \$0.00 | <input type="checkbox"/> Right-Side Lateral Support for Matrx MAC Back - 6.0" Contour Depth (Each) | \$0.00 |
| <input type="checkbox"/> Right-Side Lateral Support for Matrx MAC Back - 9.0" Contour Depth (Each) | \$0.00 | | |

Motion Concepts (MC) - Back Privacy Flap

- | | | | |
|--|----------------|--|----------------|
| <input type="checkbox"/> Matrx Back Flap Extra Large (22"SW) | \$42.00 | <input type="checkbox"/> Matrx Back Flap Large (20"SW) | \$42.00 |
| <input type="checkbox"/> Matrx Back Flap Medium (18"SW) | \$42.00 | <input type="checkbox"/> Matrx Back Flap Small (16"SW) | \$42.00 |

Motion Concepts (MC) - Back Accessories

<input type="checkbox"/> Matrx Back Meshtex Cover for MAC Hug	\$176.00	<input type="checkbox"/> Matrx Back Shape Insert Customizable interface for asymmetries. Works best with HUG model	\$98.00
<input type="checkbox"/> Matrx Backbone Lumbar/PSIS Support	\$57.00	<input type="checkbox"/> Matrx E2 Chest Harness	\$78.00
<input type="checkbox"/> Matrx E2 PostureFlex Abdominal Support (7"H x 20"-32"L)	\$253.00	<input type="checkbox"/> Matrx EasySet Fixed Hardware Upgrade	\$52.00
<input type="checkbox"/> Matrx Extended Range Hardware Upgrade Recommended with SHAPE insert	\$78.00	<input type="checkbox"/> Matrx Military Service Decal - US Air Force	\$56.00
<input type="checkbox"/> Matrx Military Service Decal - US Army	\$56.00	<input type="checkbox"/> Matrx Military Service Decal - US Marines	\$56.00
<input type="checkbox"/> Matrx Military Service Decal - US Navy	\$56.00	<input type="checkbox"/> Matrx MiniSet Fixed Hardware Upgrade (16"BH)	\$52.00
<input type="checkbox"/> Matrx Quad Release Tab for PostureFlex	\$26.00	<input type="checkbox"/> Meshtex Cover for Matrx MAC Back	\$150.00
<input type="checkbox"/> Slip Cover for Matrx E2 Back Support	\$67.00	<input type="checkbox"/> Slip Cover for Matrx E2 Deep Back Support	\$67.00

Motion Concepts (MC) - Headrest Hardware

<input type="checkbox"/> Matrx Elan Flipp Head Support Hardware (E0128) Required with selection of head support	\$525.00	<input type="checkbox"/> Matrx Elan Head Support Hardware (E1028) Required with selection of head support	\$220.00
<input type="checkbox"/> Matrx Loxx Flipp Head Support Hardware (E0128) Required with selection of head support	\$515.00	<input type="checkbox"/> Matrx Loxx Head Support Hardware (E1028) Required with selection of head support	\$215.00

Motion Concepts (MC) - Headrest Pad

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Matrx Elan Head Support 10" Standard Cushion (E0955) | \$235.00 | <input type="checkbox"/> Matrx Elan Head Support 14" Standard Cushion (E0955) | \$235.00 |
| <input type="checkbox"/> Matrx Elan Head Support Large 4-Point Cushion (E0955) | \$325.00 | <input type="checkbox"/> Matrx Elan Head Support Small 4-Point Cushion (E0955) | \$325.00 |
| <input type="checkbox"/> Matrx Elan Headband Head Support with Large 4-Point Cushion (E0955) | \$450.00 | <input type="checkbox"/> Matrx Elan Headband Head Support with Small 4-Point Cushion (E0955) | \$450.00 |
| <input type="checkbox"/> Matrx Elan Hug Head Support - 10" Cushion (E0955) | \$235.00 | <input type="checkbox"/> Matrx Elan HugHead Support - 14" Cushion (E0955) | \$235.00 |
| <input type="checkbox"/> Matrx Head Support Hardware Only - No Cushion | \$0.00 | | |

Motion Concepts (MC) - Lateral Torso Supports

- | | | | |
|---|-----------------|--|-----------------|
| <input type="checkbox"/> Matrx E2 PR Lateral EasyFix Fixed (5"x6.5") (Pair) (E0956) | \$290.00 | <input type="checkbox"/> Matrx PR Lateral Adjustable Swingaway Shoulder Support With Cushion (4.25"x6.5") (Pair) (E0956/E0128) | \$335.00 |
| <input type="checkbox"/> Matrx PR Lateral Fixed (5"x7.5") (Pair) (E0956/E0128) | \$110.00 | | |

O2 Holder**Add Oxygen Holder**

- | | |
|--|-----------------|
| <input type="checkbox"/> Oxygen Holder | \$193.00 |
|--|-----------------|

O2 Holder - Factory Installed on Chair

- | | |
|--|---------------|
| <input type="checkbox"/> Oxygen Holder - For 587 | \$0.00 |
|--|---------------|
- If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

O2 Holder - Frame Color

- | | |
|--|---------------|
| <input type="checkbox"/> Oxygen Holder - Black Frame | \$0.00 |
|--|---------------|

O2 Holder - Mounting Type

- | | |
|---|---------------|
| <input type="checkbox"/> Oxygen Holder - Single Clamp and Adapter Mount | \$0.00 |
|---|---------------|

O2 Holder - Mounting Location

Oxygen Holder - Front Mount **\$0.00**

Lap Tray**Add Lap Tray**

Lap Tray **\$0.00**

Lap Tray - Top Color

Full Tray (Black) **\$0.00**

Lap Tray - Upholstery Color

Full Tray COVER - Modena Black Vinyl **\$200.00**

Lap Tray - Retaining Strap

Full Tray - Sew-On Strap - Left-Hand Buckle **\$0.00**

Full Tray - Sew-On Strap - Right-Hand Buckle **\$0.00**

Full Tray - Standard Strap - Center Buckle **\$0.00**

Full Tray - Standard Strap - Left-Hand Buckle **\$0.00**

Full Tray - Standard Strap - Right-Hand Buckle **\$0.00**

Veteran**Military Veteran Order**

VA Order "Honor Our Veterans" **\$0.00**

Branch of Service Crest

Crest U.S. Air Force **\$56.00**

Crest U.S. Army **\$56.00**

Crest U.S. Coast Guard **\$56.00**

Crest U.S. Marine Corps **\$56.00**

Crest U.S. Navy **\$56.00**

Shipping & Freight

Initial Shipping Costs

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Ship to Authorized Broda Dealer | \$165.00 | <input type="checkbox"/> Ship Direct to Customer | \$200.00 |
|--|-----------------|--|-----------------|

Additional Delivery Services

- | | | | |
|--|-----------------|---|-------------------------|
| <input type="checkbox"/> Lift Gate Services
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$200.00 | <input type="checkbox"/> Residential Delivery Service
Delivery is to a residential location. | \$200.00 |
| <input type="checkbox"/> Notification Before Delivery
Contact is requested to the cosignee by the carrier PRIOR to shipment delivery. | \$85.00 | <input type="checkbox"/> Inside Delivery Service
Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | Min
\$250.00 |
| <input type="checkbox"/> Limited Access Delivery Service
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$165.00 | | |