

Customer Details

Billing Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Shipping Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Order Information

Date: _____

PO Number: _____

Customer Number: _____

Tag/Comments: _____



Standard Features

- 250lb Maximum Occupant Weight (rocking)
- Up to 40° Infinitely Adjustable Seat Tilt (non-rocking)
- 350lb Maximum Occupant Weight (non-rocking)
- Rocking Mechanism (for up to 33° of tilt)
- 500 lb Maximum Occupant Weight with the Heavy Duty Package (no rocking only)
- Up to 33° Infinitely Adjustable Seat Tilt (rocking)
- Tubular Steel, Powder Coated Frame
- Seat Mount Lower Leg Support Assembly
- Adjustable Seat Surface Height
- 2 x 20" Front Mag Wheels with Front Anti-Tippers
- Independent Adjustable Arm Support Height
- 2 x 5" Total Lock Rear Caster Wheels

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

Product Options

Chair Model

- | | | | |
|--|-------------------|--|-------------------|
| <input type="checkbox"/> Encore Pedal Wheelchair (E1161) | \$3,067.00 | <input type="checkbox"/> Encore Pedal Wheelchair with Mag Wheels (E1161) | \$3,731.00 |
| <input type="checkbox"/> Encore Pedal Wheelchair with Mag Wheels (w/ WC19) (E1161) | \$4,160.00 | | |
- Available up to 24" seat width, unless selected with the Heavy Duty package. Available with mag wheels only.

Weight Capacity

- | | | | |
|---|---------------|--|-----------------|
| <input type="checkbox"/> Standard Weight Capacity (350 LBS) | \$0.00 | <input type="checkbox"/> Heavy Duty (HD) Weight Capacity (500 LBS) | \$567.00 |
|---|---------------|--|-----------------|
- Available in 22"-28" seat widths only. Non-rocking only. Aluminum front casters not available. Sub-mount seat hanger style only with pin & saddle, angle-adjustable footplates (F-SW-ANGLE) only. If combined with WC19 package, only mag wheels are available (20", 22" front or 14" rear with twin front casters).

Seat Type

- | | | | |
|--|-----------------|---|---------------|
| <input type="checkbox"/> Seat With Strap | \$0.00 | <input type="checkbox"/> Seat Without Strap | \$0.00 |
| <input type="checkbox"/> Seat Pan (Seat Without Strap) (E2231) | \$184.00 | | |

Back Type

- | | | | |
|--|---------------|---|---------------|
| <input type="checkbox"/> Back With Strap | \$0.00 | <input type="checkbox"/> Back Without Strap | \$0.00 |
|--|---------------|---|---------------|

Seat Width

- | | | | |
|---|-------------------|---|-------------------|
| <input type="checkbox"/> 16" Seat Width | \$0.00 | <input type="checkbox"/> 18" Seat Width | \$0.00 |
| <input type="checkbox"/> 20" Seat Width (E2201) | \$0.00 | <input type="checkbox"/> 22" Seat Width (E2201) | \$500.00 |
| <input type="checkbox"/> 24" Seat Width (E2202) | \$1,131.00 | <input type="checkbox"/> 26" Seat Width (E2202) | \$1,131.00 |
| <input type="checkbox"/> 28" Seat Width | \$1,131.00 | | |

Seat Depth

- | | |
|--|--|
| <input type="checkbox"/> 17" Seat Depth (Check Lower Leg Support Assembly Option) \$0.00 | <input type="checkbox"/> 18.5" Seat Depth (Check Lower Leg Support Assembly Option) \$0.00 |
| <input type="checkbox"/> 20" Seat Depth (Check Lower Leg Support Assembly Option) (E2203) \$121.00 | |

Seat Height

- | | |
|--|---|
| <input type="checkbox"/> 12.5" Seat Height \$0.00 | <input type="checkbox"/> 13.5" Seat Height \$0.00 |
| <input type="checkbox"/> 14.5" Seat Height \$0.00 | <input type="checkbox"/> 15.5" Seat Height \$0.00 |
| <input type="checkbox"/> 16.5" Seat Height \$0.00 | <input type="checkbox"/> 17.5" Seat Height \$0.00 |
| <input type="checkbox"/> 17.5" Seat Height-40 Degree Tilt \$0.00 | |

Back Height

- | | |
|---|---|
| <input type="checkbox"/> 22" Back Height -\$505.00 | <input type="checkbox"/> 32" Back Height \$0.00 |
| <input type="checkbox"/> 33.5" Back Height \$121.00 | <input type="checkbox"/> 35" Back Height \$121.00 |

Base Frame Color

- | | |
|--|---|
| <input type="checkbox"/> Black Base Frame \$0.00 | <input type="checkbox"/> Ocean Blue Base Frame \$0.00 |
| <input type="checkbox"/> Bronze Base Frame \$0.00 | <input type="checkbox"/> Satin Nickel Base Frame \$0.00 |
| <input type="checkbox"/> Deep Orange Frame \$0.00 | <input type="checkbox"/> Carmine Red Base Frame \$0.00 |
| <input type="checkbox"/> Tungsten Sand Base Frame \$0.00 | <input type="checkbox"/> Zinc Yellow Base Frame \$0.00 |
| <input type="checkbox"/> Passion Pink Base Frame \$0.00 | <input type="checkbox"/> Royal Purple Base Frame \$0.00 |
| <input type="checkbox"/> Forest Green Base Frame \$0.00 | |

Upholstery Color

<input type="checkbox"/> Modena Black Cushions	\$0.00	<input type="checkbox"/> Espresso Cushions	\$0.00
<input type="checkbox"/> Interlace Graphite Cushions	\$0.00	<input type="checkbox"/> Interlace Indigo Cushions	\$0.00
<input type="checkbox"/> Beeline Portobello Cushions	\$0.00	<input type="checkbox"/> Beeline Geyser Cushions	\$0.00
<input type="checkbox"/> Kaleidoscope Inclination Smokehouse Cushions	\$0.00	<input type="checkbox"/> Milano Stitch Bark Cushions	\$0.00
<input type="checkbox"/> Tweed Mirage Cushions	\$0.00		

Wheel Options

<input type="checkbox"/> 20" Front Mag Wheels (Pair)	\$0.00	<input type="checkbox"/> 22" Front Mag Wheels (Pair)	\$0.00
<input type="checkbox"/> 22"- Heavy Duty Front Mag Wheels (Pair)	\$0.00	<input type="checkbox"/> 1 Rear Directional Lock Caster	\$0.00
<input type="checkbox"/> 2 Rear Fixed Caster Wheels (Pair)	\$0.00	<input type="checkbox"/> 14" Rear Mag Wheels (Pair)	\$0.00

Mag Wheel Location

<input type="checkbox"/> Frontmost Fore-aft Mag Wheel Mounting Port	\$0.00	<input type="checkbox"/> Middle Fore-aft Mag Wheel Mounting Port	\$0.00
<input type="checkbox"/> Rearmost Fore-aft Mag Wheel Mounting Port	\$0.00		

Wheel Lock Position

<input type="checkbox"/> Rear Mount Wheel Locks (Pair)	\$0.00	<input type="checkbox"/> Front Mount Wheel Locks (Pair)	\$0.00
--	--------	---	--------

Wheel Lock Extensions

<input type="checkbox"/> With Parking Brake Extensions (Pair) (E0961)	\$81.00
---	---------

Anti-Tippers

<input type="checkbox"/> Anti-Tip Devices (Required) (Pair) (E0971)	\$0.00	<input type="checkbox"/> Extended Front Anti-Tip Devices (Required) (Pair) (E0971)	\$0.00
---	--------	--	--------

Casters

- | | | | |
|---|---------------|--|-----------------|
| <input type="checkbox"/> 2 X 5" Total Lock Caster Wheels (Pair) | \$0.00 | <input type="checkbox"/> 2 X 5" Aluminum Hub Forked Caster Wheels (Pair) | \$241.00 |
| <input type="checkbox"/> 2 X 5" Twin Wheel Caster Wheels (Pair) | \$0.00 | | |

Lateral Supports

- | | | | |
|--|---------------|--------------------------------------|---------------|
| <input type="checkbox"/> Lateral Mount Wide (Pair) (E1028) | \$0.00 | <input type="checkbox"/> No Laterals | \$0.00 |
|--|---------------|--------------------------------------|---------------|

Comfort Co. Lateral Support

- | | | | |
|---|-----------------|--|--|
| <input type="checkbox"/> Comfort Co. Lateral Support (E0956)
Compatible with seat widths up to 24" | \$564.00 | | |
|---|-----------------|--|--|

Arm Top

- | | | | |
|--|-----------------|---|---------------|
| <input type="checkbox"/> Depth Adjustable Arm Top (Pair) (E0973) | \$165.00 | <input type="checkbox"/> Height Adjustable Arm Top Short (Pair) (E0973) | \$0.00 |
|--|-----------------|---|---------------|

Arm Pad

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Narrow Arm Supports (Pair) | \$0.00 | <input type="checkbox"/> Cushioned Vinyl Arm Support (Black) (Pair) | \$95.00 |
| <input type="checkbox"/> Soft Arm Supports (Pair) | \$0.00 | <input type="checkbox"/> Arm Trough (Pair) (E2209) | \$329.00 |
| <input type="checkbox"/> Arm Trough (L) & Soft Arm Support (R) (E2209) | \$165.00 | <input type="checkbox"/> Arm Trough (R) & Soft Arm Support (L) (E2209) | \$165.00 |

Flip-Down Armrest Pad

- | | | | |
|---|-----------------|--|--|
| <input type="checkbox"/> Padded Foam Overlay - Armrest-Regular (Pair) | \$147.00 | | |
|---|-----------------|--|--|

Swing Away Footrest

- | | | | |
|---|-----------------|---|-----------------|
| <input type="checkbox"/> No Lower Leg Support | \$0.00 | <input type="checkbox"/> Swing-Away Lower Leg Support (Drop-In) [Style 24/25] (Pair) | \$133.00 |
| <input type="checkbox"/> Swing-Away Lower Leg Support Angle-Adjustable Foot Plate (Drop-In) [Style 34/35] (Pair) | \$247.00 | <input type="checkbox"/> Swing-Away Lower Leg Support Small Angle-Adjustable Foot Plate (Drop-In) [Style 34/35+430609004] (Pa | \$247.00 |
| <input type="checkbox"/> Angle-Adjustable Lower Leg Support Small (Drop-In) [Style 30/31] (Pair) (E0990) | \$550.00 | <input type="checkbox"/> Angle-Adjustable Lower Leg Support - Steel Button (Drop-In) [Style 32/33] (Pair) (E0990) | \$550.00 |
| <input type="checkbox"/> Swing-Away Lower Leg Support (Pin & Saddle) [Style 03/04] (Pair) | \$203.00 | <input type="checkbox"/> Swing-Away Lower Leg Support with Angle-Adjustable Foot Plate (Pin & Saddle) [Style 05/06] (Pair) | \$310.00 |
| <input type="checkbox"/> HD Modified Swing-Away Lower Leg Support Angle-Adjustable Foot Plate (Pin & Saddle) [Style 05/06] (P | \$310.00 | <input type="checkbox"/> Angle-Adjustable Lower Leg Support (Pin & Saddle) [Style 07/08] (Pair) (E0990) | \$266.00 |
| <input type="checkbox"/> Angle-Adjustable Lower Leg Support Angle-Adjustable Foot Plate (Pin & Saddle) [Style 22/23] (Pair) (E0990) | \$651.00 | | |

Footrest Extension

- | | |
|--|----------------|
| <input type="checkbox"/> Ikon Extended Footrest Adjustment Tube (Pair) | \$74.00 |
|--|----------------|

Foot Plate

- | | |
|--|---------------|
| <input type="checkbox"/> With Angle-Adjustable Foot Plate (Pair) (K0040) | \$0.00 |
|--|---------------|

Footrest Heel Loop

- | | |
|--|---------------|
| <input type="checkbox"/> With Heel Loop (Pair) (E0951) | \$0.00 |
|--|---------------|

Footrest Hanger / Front Rigging

- | | |
|--|--|
| <input type="checkbox"/> Lower Leg Support Assembly With Receiver (Sub-Seat) (Pin & Saddle) (18.5+20) (Pair) \$0.00 | <input type="checkbox"/> Lower Leg Support Assembly Short Seat Mount (Pin & Saddle) (Pair) \$0.00 |
| <input type="checkbox"/> Lower Leg Support Assembly With Receiver (Sub-Seat) (Drop-In) (17+18.5+20) (Pair) \$0.00 | <input type="checkbox"/> Lower Leg Support Assembly Short Seat Mount (Drop-In) (Pair) \$0.00 |
| <input type="checkbox"/> Lower Leg Support Assembly Receiver Only (Pair) \$0.00 | <input type="checkbox"/> No Lower Leg Support Assembly + No Receiver \$0.00 |

Padding Package

- | | |
|---|---|
| <input type="checkbox"/> Standard Cushions \$0.00
Includes: Neck Support, Adjustable Lateral Shoulder Supports, Lower Side Covers, Seat and Back Support Cushions with Covers | <input type="checkbox"/> Additional Padding Package (APP) or Slim Posture Padding (SPP) Only (Pricing Varies by Size) \$0.00
Includes: Seat, Back, and Head Support Cushions (no Head Support or Lateral Shoulder Supports on short back) plus adjustable, removable lateral shoulder supports and lateral upper leg support cushions |
| <input type="checkbox"/> Standard Cushion + Additional Padding Package (APP) or Slim Posture Padding (SPP) (Pricing Varies by) \$0.00 | |

Neckrest Pads

- | | |
|---|--|
| <input type="checkbox"/> Neck Support (Matches Selected Accent Color) (E0955) \$0.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello | <input type="checkbox"/> Contoured Visco Neck Support (Triptex) (E0955) \$51.00 |
| <input type="checkbox"/> Contoured Visco Neck Support (Vinyl) (Matches Selected Accent Color) (E0955) \$51.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello | <input type="checkbox"/> Soft Visco Neck Support Cushion (E0955) \$51.00 |

Shoulder Bolsters

- | | |
|---|---|
| <input type="checkbox"/> Standard Shoulder Supports (Matches Selected Accent Color) (Pair) (E0956) \$0.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello | <input type="checkbox"/> Long Shoulder Supports (Matches Selected Accent Color) (Pair) (E0956) \$134.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello |
|---|---|

Seat Pad

<input type="checkbox"/> Full Seat Cover (Black)	\$60.00	<input type="checkbox"/> Full Seat Cover (Matches Selected Accent Color)	\$60.00
		Available in all upholstery colors except Beeline Espresso and Beeline Portobello	
<input type="checkbox"/> Full Seat Cover With Slot (Black)	\$60.00	<input type="checkbox"/> Full Seat Cover With Slot (Matches Selected Accent Color)	\$60.00
		Available in all upholstery colors except Beeline Espresso and Beeline Portobello	
<input type="checkbox"/> Terry Cloth Seat Cover with Cushion (Grey)	\$0.00	<input type="checkbox"/> Removable Seat Covers With Cushion (Black)	\$0.00
<input type="checkbox"/> Removable Seat Covers With Cushion (Matches Selected Accent Color)	\$0.00	<input type="checkbox"/> Removable Seat Covers (Startex) With Cushion	\$0.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello			

Back Pad

<input type="checkbox"/> Full Back Cover (Black)	\$61.00	<input type="checkbox"/> Full Back Cover (Matches Selected Accent Color)	\$61.00
		Available in all upholstery colors except Beeline Espresso and Beeline Portobello	
<input type="checkbox"/> Back Terry Cloth Cover with Cushion (Grey)	\$0.00	<input type="checkbox"/> Removable Back Covers (Black) With Cushion	\$0.00
<input type="checkbox"/> Removable Back Covers With Cushion (Matches Selected Accent Color)	\$0.00	<input type="checkbox"/> Removable Back Covers (Startex) With Cushion	\$0.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello			

Pelvic Positioning Belt

- | | |
|---|--|
| <input type="checkbox"/> Pelvic Positioning Belt Side Buckle (E0978) \$64.00 | <input type="checkbox"/> Pelvic Positioning Belt Center Buckle (E0978) \$64.00 |
| <input type="checkbox"/> Pelvic Positioning Belt Side Auto-Style Buckle (E0978) \$64.00 | <input type="checkbox"/> Pelvic Positioning Belt Center Auto-Style Buckle (E0978) \$64.00 |
| <input type="checkbox"/> Padded Pelvic Positioning Belt (Black) (E0978) \$91.00 | <input type="checkbox"/> Padded Pelvic Positioning Belt (Matches Selected Accent Color) (E0978) \$91.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello |
| <input type="checkbox"/> Padded Pelvic Positioning Belt with Auto-Style Buckle (Black) (E0978) \$91.00 | <input type="checkbox"/> Padded Pelvic Positioning Belt with Auto-Style Buckle (Matches Selected Accent Color) (E0978) \$91.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello |

Thigh Positioning Belt

- | | |
|---|--|
| <input type="checkbox"/> Thigh Positioning Belt (Black) (E0957) \$146.00 | <input type="checkbox"/> Thigh Positioning Belt (Matches Selected Accent Color) (E0957) \$146.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello |
| <input type="checkbox"/> Thigh Positioning Belt with Auto-Style Buckle (Black) (E0957) \$152.00 | <input type="checkbox"/> Thigh Positioning Belt (Neoprene) (E0957) \$146.00 |
| <input type="checkbox"/> Thigh Positioning Belt with Auto-Style Buckle (Matches Selected Accent Color) (E0957) \$152.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello | <input type="checkbox"/> Thigh Positioning Belt with Auto-Style Buckle (Neoprene) (E0957) \$152.00 |
| <input type="checkbox"/> APP Extended Neoprene Thigh Positioning Belt (E0957) \$146.00 | <input type="checkbox"/> APP Extended Neoprene Thigh Positioning Belt with Auto-Style Buckle (E0957) \$152.00 |

Side Reduction Pad

- | | |
|---|--|
| <input type="checkbox"/> Side Reduction Pad (Black) (E0956) \$121.00 | <input type="checkbox"/> Side Reduction Pad (Matches Selected Accent Color) (E0956) \$121.00
Available in all upholstery colors except Beeline Espresso and Beeline Portobello |
|---|--|

Rear Stabilizers

- | |
|--|
| <input type="checkbox"/> Rear Stabilizers (Pair) \$108.00 |
|--|

Lateral Torso Supports

- Lateral Torso Supports (Pair) (E0956) **\$695.00**
Factory installed recommended

Accessories

Accessories - Add-On

- Oxygen Carry Bag **\$109.00** Cup Holder (Grey Swivel) **\$57.00**

Additional Positioning Padding (APP)-Slim Posture Padding (SPP)

Additional Positioning Padding (APP) - Package Price

Broda APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Covers (Not compatible with 35" back support height)

- APP Package - Mobility - 16" 18" 20" Widths **\$464.00**
Includes Seat, Back and Head Support Cushions and Lateral Shoulder Supports. Not compatible with 35" back support height
- APP Package - Mobility - 22" 24" Widths **\$708.00**
Includes Seat, Back and Head Support Cushions and Lateral Shoulder Supports. Not compatible with 35" back support height
- APP Package - Mobility - 26" 28" Widths **\$839.00**
Includes Seat, Back and Head Support Cushions and Lateral Shoulder Supports. Not compatible with 35" back support height

Slim Posture Padding (SPP) - Seat Cushion

<input type="checkbox"/> SPP Seat Cushion 16" Seat Width X 17" Seat Depth (E2601) \$212.00	<input type="checkbox"/> SPP Seat Cushion 16" Seat Width X 18.5" Seat Depth (E2601) \$212.00
<input type="checkbox"/> SPP Seat Cushion 16" Seat Width X 20" Seat Depth (E2601) \$212.00	<input type="checkbox"/> SPP Seat Cushion 18" Seat Width X 17" Seat Depth (E2601) \$212.00
<input type="checkbox"/> SPP Seat Cushion 18" Seat Width X 18.5" Seat Depth (E2601) \$212.00	<input type="checkbox"/> SPP Seat Cushion 18" Seat Width X 20" Seat Depth (E2601) \$212.00
<input type="checkbox"/> SPP Seat Cushion 20" Seat Width X 17" Seat Depth (E2601) \$212.00	<input type="checkbox"/> SPP Seat Cushion 20" Seat Width X 18.5" Seat Depth (E2601) \$212.00
<input type="checkbox"/> SPP Seat Cushion 20" Seat Width X 20" Seat Depth (E2601) \$212.00	<input type="checkbox"/> SPP Seat Cushion 22" Seat Width X 17" Seat Depth (E2602) \$234.00
<input type="checkbox"/> SPP Seat Cushion 22" Seat Width X 18.5" Seat Depth (E2602) \$234.00	<input type="checkbox"/> SPP Seat Cushion 22" Seat Width X 20" Seat Depth (E2602) \$234.00
<input type="checkbox"/> SPP Seat Cushion 24" Seat Width X 17" Seat Depth (E2602) \$234.00	<input type="checkbox"/> SPP Seat Cushion 24" Seat Width X 18.5" Seat Depth (E2602) \$234.00
<input type="checkbox"/> SPP Seat Cushion 24" Seat Width X 20" Seat Depth (E2602) \$234.00	<input type="checkbox"/> SPP Seat Cushion 26" Seat Width X 17" Seat Depth (E2602) \$262.00
<input type="checkbox"/> SPP Seat Cushion 26" Seat Width X 18.5" Seat Depth (E2602) \$262.00	<input type="checkbox"/> SPP Seat Cushion 26" Seat Width X 20" Seat Depth (E2602) \$262.00

Accessory Bag**Accessory Bag**

<input type="checkbox"/> Accessory Bag (Single) \$51.00	<input type="checkbox"/> Accessory Bag (Pair) \$102.00
<input type="checkbox"/> Embroidered Accessory Bag (Single) \$83.00	<input type="checkbox"/> Embroidered Accessory Bag (Pair) \$165.00

Accessory Bag - Embroidery

<input type="checkbox"/> Accessory Bag Embroidery - Block Font \$0.00	<input type="checkbox"/> Accessory Bag Embroidery - Brush Font \$0.00
--	--

Accessory Bag - Characters

- | | | | |
|---|---------------|---|---------------|
| <input type="checkbox"/> Accessory Bag - 2 Characters | \$0.00 | <input type="checkbox"/> Accessory Bag - 1 Character Only | \$0.00 |
|---|---------------|---|---------------|

Accessory Bag - 1st Character

Accessory Bag - 2nd Character

IV Pole**Add IV-Pole**

- | | |
|----------------------------------|-----------------|
| <input type="checkbox"/> IV-Pole | \$121.00 |
|----------------------------------|-----------------|

IV Pole - Factory Installed on Chair

- | | | | |
|--|---------------|---|---------------|
| <input type="checkbox"/> IV Pole - For Encore Pedal Without Mag Wheels | \$0.00 | <input type="checkbox"/> IV Pole - For Encore Pedal With Mag Wheels | \$0.00 |
| | | Front mount required if selecting rear mag wheels | |

IV Pole - Frame Color

- | | |
|---|---------------|
| <input type="checkbox"/> IV Pole - Black Standard Frame | \$0.00 |
|---|---------------|

IV Pole - Mount

- | | |
|--|---------------|
| <input type="checkbox"/> IV Pole - Standard IV-Pole Bracket Only | \$0.00 |
|--|---------------|

IV Pole - Location

- | | | | |
|---|---------------|---|---------------|
| <input type="checkbox"/> IV Pole - Front Mount | \$0.00 | <input type="checkbox"/> IV Pole - Rear Mount | \$0.00 |
| Front mount required if selecting rear mag wheels | | Front mount required if selecting rear mag wheels | |

Matrx

Motion Concepts (MC) - Seat Model

<input type="checkbox"/> Matrx Libra Seat Cushion 16" 18" 20" Seat Width (E2624) \$798.00 Recommended for highest risk of or existing skin breakdown	<input type="checkbox"/> Matrx Libra Seat Cushion 22" 24" Seat Width (E2625) \$1,097.00 Recommended for highest risk of or existing skin breakdown
<input type="checkbox"/> Matrx Vi Seat Cushion 16" 18" 20" Seat Width (E2607) \$668.00 Recommended for moderate to high risk of skin breakdown	<input type="checkbox"/> Matrx Vi Seat Cushion 22" Seat Width (E2608) \$729.00 Recommended for moderate to high risk of skin breakdown
<input type="checkbox"/> Matrx PS Positioning Seat Cushion 16" 18" 20" Seat Width (E2605) \$609.00 Recommended for low to moderate risk of skin breakdown	<input type="checkbox"/> Matrx PS Positioning Seat Cushion 22" 24" 26" Seat Width (E2606) \$1,044.00 Recommended for low to moderate risk of skin breakdown
<input type="checkbox"/> Matrx PSP Positioning Seat Cushion 16" 18" 20" Seat Width (E2607) \$691.00 Recommended for moderate to high risk of skin breakdown	<input type="checkbox"/> Matrx Multi Seat Cushion 16" 17" 18" 20" Seat Width (E2607) \$745.00 Recommended for moderate to high risk of skin breakdown
<input type="checkbox"/> Matrx Multi Seat Cushion 22" Seat Width (E2608) \$995.00 Recommended for moderate to high risk of skin breakdown	

Motion Concepts (MC) - Seat Width

<input type="checkbox"/> Matrx Seat Cushion 16" Seat Width \$0.00	<input type="checkbox"/> Matrx Seat Cushion 18" Seat Width \$0.00
<input type="checkbox"/> Matrx Seat Cushion 20" Seat Width \$0.00	<input type="checkbox"/> Matrx Seat Cushion 22" Seat Width \$0.00
<input type="checkbox"/> Matrx Seat Cushion 24" Seat Width \$0.00	<input type="checkbox"/> Matrx Seat Cushion 26" Seat Width \$0.00

Motion Concepts (MC) - Seat Depth

<input type="checkbox"/> Matrx Seat Cushion 16" Deep \$0.00	<input type="checkbox"/> Matrx Seat Cushion 18" Deep \$0.00
<input type="checkbox"/> Matrx Seat Cushion 20" Deep \$0.00	<input type="checkbox"/> Matrx Seat Cushion 22" Deep \$0.00 Only available w/ Multi & PS Cushions

Motion Concepts (MC) - Seat Rigidizer

<input type="checkbox"/> Matrx Seat Rigidizer 16 Small (16" Seat Depth) (E0992)	\$124.00	<input type="checkbox"/> Matrx Seat Rigidizer 16 Medium (18" Seat Depth) (E0992)	\$124.00
<input type="checkbox"/> Matrx Seat Rigidizer 16 Large (20" 22" Seat Depth) (E0992)	\$124.00	<input type="checkbox"/> Matrx Seat Rigidizer 18 Small (16" Seat Depth) (E0992)	\$124.00
<input type="checkbox"/> Matrx Seat Rigidizer 18 Medium (18" Seat Depth) (E0992)	\$124.00	<input type="checkbox"/> Matrx Seat Rigidizer 18 Large (20" 22" Seat Depth) (E0992)	\$124.00
<input type="checkbox"/> Matrx Seat Rigidizer 20 Small (16" Seat Depth) (E0992)	\$124.00	<input type="checkbox"/> Matrx Seat Rigidizer 20 Medium (18" Seat Depth) (E0992)	\$124.00
<input type="checkbox"/> Matrx Seat Rigidizer 20 Large (20" 22" Seat Depth) (E0992)	\$124.00	<input type="checkbox"/> Matrx Seat Rigidizer 22 Medium (18" Seat Depth) (E0992)	\$124.00
<input type="checkbox"/> Matrx Seat Rigidizer 22 Large (20" 22" Seat Depth) (E0992)	\$124.00		

Motion Concepts (MC) - Seat Accessories

<input type="checkbox"/> Matrx 1" Anterior Wedge	\$217.00	<input type="checkbox"/> Matrx 1" Hemi Wedge	\$109.00
<input type="checkbox"/> Matrx Libra Abductor	\$91.00	<input type="checkbox"/> Matrx Libra Fluid Supplement	\$129.00
<input type="checkbox"/> Matrx Pelvic Obliquity Kit	\$91.00	<input type="checkbox"/> Matrx Tapered Wedge (Left)	\$91.00
<input type="checkbox"/> Matrx Tapered Wedge (Right)	\$91.00	<input type="checkbox"/> Matrx Seat Slip Cover	\$95.00
<input type="checkbox"/> Extra Matrx Seat Cover - PS (Match Seat Width/Depth)	\$129.00	<input type="checkbox"/> Extra Matrx Seat Cover - PS (Match Seat Width/Depth - 22" Wide or Greater)	\$268.00
<input type="checkbox"/> Extra Matrx Seat Cover - PSP (Match Seat Width/Depth)	\$129.00	<input type="checkbox"/> Extra Matrx Seat Cover - VI (Match Seat Width/Depth)	\$129.00
<input type="checkbox"/> Extra Matrx Seat Cover - Libra (Match Seat Width/Depth)	\$134.00	<input type="checkbox"/> Extra Matrx Seat Cover - Libra (Match Seat Width/Depth - 22" Wide or Greater)	\$161.00
<input type="checkbox"/> Extra Matrx Seat Cover - Multi (Match Seat Width/Depth) (E2619)	\$120.00	<input type="checkbox"/> Extra Matrx Seat Cover - Multi (Match Seat Width/Depth - 22" Wide or Greater) (E2619)	\$250.00

Motion Concepts (MC) - Back Model

<input type="checkbox"/> Matrx E2 Back - 16" 17" 18" 19" 20" Seat Width (E2615) Recommended for users seeking optimal immersion and airflow	\$965.00	<input type="checkbox"/> Matrx E2 Deep Back - 16" 17" 18" 19" 20" Seat Width (E2620) Recommended for users seeking optimal immersion w/ contoured support needs	\$1,125.00
<input type="checkbox"/> Matrx MAC Foam Back - 16" 18" 20"SW (E2620)	\$1,495.00	<input type="checkbox"/> Matrx MAC Hug Back - 16" 18" 20"SW (E2620)	\$1,645.00

Motion Concepts (MC) - Back Width

<input type="checkbox"/> Matrx Back Width - 16"	\$0.00	<input type="checkbox"/> Matrx Back Width - 18"	\$0.00
<input type="checkbox"/> Matrx Back Width - 20"	\$0.00	<input type="checkbox"/> Matrx Back Width - 22"	\$0.00

Motion Concepts (MC) - Back Height

- | | | | |
|--|--------|---|--------|
| <input type="checkbox"/> Matrx Back Height 16" | \$0.00 | <input type="checkbox"/> Matrx Back Height 18"
Excludes PB & PB Deep | \$0.00 |
| <input type="checkbox"/> Matrx Back Height 20" | \$0.00 | <input type="checkbox"/> Matrx Back Height 22"
MAC Back only | \$0.00 |

Motion Concepts (MC) - Back MAC Lateral LH

- | | | | |
|--|--------|--|--------|
| <input type="checkbox"/> Left-Side Lateral Support for Matrx MAC Back - 3.5"
Contour Depth (Each) | \$0.00 | <input type="checkbox"/> Left-Side Lateral Support for Matrx MAC Back - 6.0"
Contour Depth (Each) | \$0.00 |
| <input type="checkbox"/> Left-Side Lateral Support for Matrx MAC Back - 9.0"
Contour Depth (Each) | \$0.00 | | |

Motion Concepts (MC) - Back MAC Lateral RH

- | | | | |
|---|--------|---|--------|
| <input type="checkbox"/> Right-Side Lateral Support for Matrx MAC Back - 3.5"
Contour Depth (Each) | \$0.00 | <input type="checkbox"/> Right-Side Lateral Support for Matrx MAC Back - 6.0"
Contour Depth (Each) | \$0.00 |
| <input type="checkbox"/> Right-Side Lateral Support for Matrx MAC Back - 9.0"
Contour Depth (Each) | \$0.00 | | |

Motion Concepts (MC) - Back Privacy Flap

- | | | | |
|--|---------|--|---------|
| <input type="checkbox"/> Matrx Back Flap Sall (16"SW) | \$57.00 | <input type="checkbox"/> Matrx Back Flap Medium (18"SW) | \$57.00 |
| <input type="checkbox"/> Matrx Back Flap Large (20"SW) | \$57.00 | <input type="checkbox"/> Matrx Back Flap Extra Large (22"SW) | \$57.00 |
| <input type="checkbox"/> Matrx Back Flap Extra Extra Large (24"SW) | \$57.00 | | |

Motion Concepts (MC) - Back Accessories

<input type="checkbox"/> Matrx E2 Chest Harness	\$85.00	<input type="checkbox"/> Matrx Backbone Lumbar/PSIS Support	\$75.00
<input type="checkbox"/> Matrx E2 PostureFlex Abdominal Support (7"H x 20"-32"L)	\$295.00	<input type="checkbox"/> Matrx Quad Release Tab for PostureFlex	\$35.00
<input type="checkbox"/> Matrx MiniSet Fixed Hardware Upgrade (16" BH)	\$60.00	<input type="checkbox"/> Matrx EasySet Fixed Hardware Upgrade	\$60.00
<input type="checkbox"/> Slip Cover for Matrx E2 Back Support	\$95.00	<input type="checkbox"/> Slip Cover for Matrx E2 Deep Back Support	\$95.00
<input type="checkbox"/> Meshtex Cover for Matrx MAC Back	\$195.00	<input type="checkbox"/> Matrx Back Meshtex Cover for MAC Hug	\$265.00
<input type="checkbox"/> Matrx Back Shape Insert Customizable interface for asymmetries. Works best with HUG model	\$150.00	<input type="checkbox"/> Matrx Extended Range Hardware Upgrade Recommended with SHAPE insert	\$122.00
<input type="checkbox"/> Matrx Military Service Decal - US Army	\$0.00	<input type="checkbox"/> Matrx Military Service Decal - US Navy	\$0.00
<input type="checkbox"/> Matrx Military Service Decal - US Air Force	\$0.00	<input type="checkbox"/> Matrx Military Service Decal - US Marines	\$0.00

Motion Concepts (MC) - Headrest Hardware

<input type="checkbox"/> Matrx Elan Head Support Hardware (E1028) Required with selection of head support	\$305.00	<input type="checkbox"/> Matrx Elan Flipp Head Support Hardware (E0128) Required with selection of head support	\$715.00
<input type="checkbox"/> Matrx Loxx Head Support Hardware (E1028) Required with selection of head support	\$275.00	<input type="checkbox"/> Matrx Loxx Flipp Head Support Hardware (E0128) Required with selection of head support	\$705.00

Motion Concepts (MC) - Headrest Pad

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Matrx Elan Headband Head Support with Small 4-Point Cushion (E0955) | \$653.00 | <input type="checkbox"/> Matrx Elan Headband Head Support with Large 4-Point Cushion (E0955) | \$653.00 |
| <input type="checkbox"/> Matrx Elan Head Support Small 4-Point Cushion (E0955) | \$480.00 | <input type="checkbox"/> Matrx Elan Head Support Large 4-Point Cushion (E0955) | \$480.00 |
| <input type="checkbox"/> Matrx Elan Head Support 10" Standard Cushion (E0955) | \$332.00 | <input type="checkbox"/> Matrx Elan Head Support 14" Standard Cushion (E0955) | \$332.00 |
| <input type="checkbox"/> Matrx Elan Hug Head Support - 10" Cushion (E0955) | \$332.00 | <input type="checkbox"/> Matrx Elan HugHead Support - 14" Cushion (E0955) | \$332.00 |
| <input type="checkbox"/> Matrx Head Support Hardware Only - No Cushion | \$0.00 | | |

Motion Concepts (MC) - Lateral Torso Supports

- | | | | |
|---|-----------------|--|-----------------|
| <input type="checkbox"/> Matrx PR Lateral Fixed (5"x7.5") (Pair) (E0956/E0128) | \$128.00 | <input type="checkbox"/> Matrx PR Lateral Adjustable Swingaway Shoulder Support With Cushion (4.25"x6.5") (Pair) (E0956/E0128) | \$403.00 |
| <input type="checkbox"/> Matrx E2 PR Lateral EasyFix Fixed (5"x6.5") (Pair) (E0956) | \$385.00 | | |

O2 Holder**Add Oxygen Holder**

- | | |
|--|-----------------|
| <input type="checkbox"/> Oxygen Holder | \$203.00 |
|--|-----------------|

O2 Holder - Factory Installed on Chair

- | | | | |
|---|---------------|--|---------------|
| <input type="checkbox"/> Oxygen Holder - For 48Synthesis Without Mag Wheels | \$0.00 | <input type="checkbox"/> Oxygen Holder - For 48Synthesis With Mag Wheels | \$0.00 |
| | | Front mount required if selecting rear mag wheels | |

O2 Holder - Frame Color

- | | |
|--|---------------|
| <input type="checkbox"/> Oxygen Holder - Black Frame | \$0.00 |
|--|---------------|

O2 Holder - Mounting Type

- | | |
|---|---------------|
| <input type="checkbox"/> Oxygen Holder - Single Clamp and Adapter Mount | \$0.00 |
|---|---------------|

O2 Holder - Mounting Location

- | | | | |
|---|---------------|--|---------------|
| <input type="checkbox"/> Oxygen Holder - Front Mount
Front mount required if selecting rear mag wheels | \$0.00 | <input type="checkbox"/> Oxygen Holder - Rear Mount
Front mount required if selecting rear mag wheels | \$0.00 |
|---|---------------|--|---------------|

Half Tray**Add Half Tray**

- | | |
|------------------------------------|-----------------|
| <input type="checkbox"/> Half Tray | \$512.00 |
|------------------------------------|-----------------|

Half Tray - Location

- | | | | |
|---|---------------|--|---------------|
| <input type="checkbox"/> Half Tray - Right Side | \$0.00 | <input type="checkbox"/> Half Tray - Left Side | \$0.00 |
|---|---------------|--|---------------|

Lap Tray**Add Lap Tray**

- | | |
|-----------------------------------|---------------|
| <input type="checkbox"/> Lap Tray | \$0.00 |
|-----------------------------------|---------------|

Lap Tray - Width

- | | | | |
|---|-----------------|---|-----------------|
| <input type="checkbox"/> Full Tray - 16" Seat Width (E0950) | \$373.00 | <input type="checkbox"/> Full Tray - 18" Seat Width (E0950) | \$373.00 |
| <input type="checkbox"/> Full Tray - 20" Seat Width (E0950) | \$373.00 | <input type="checkbox"/> Full Tray - 22" Seat Width (E0950) | \$373.00 |
| <input type="checkbox"/> Full Tray - 24" Seat Width (E0950) | \$537.00 | <input type="checkbox"/> Full Tray - 26" Seat Width (E0950) | \$537.00 |
| <input type="checkbox"/> Full Tray - 28" Seat Width (E0950) | \$537.00 | | |

Lap Tray - Top Color

- | | | | |
|--|----------------|---|---------------|
| <input type="checkbox"/> Full Tray (Black) | \$0.00 | <input type="checkbox"/> Full Tray (Gray) | \$0.00 |
| <input type="checkbox"/> Full Tray (Clear) | \$13.00 | | |

Lap Tray - Upholstery Color

- | | |
|---|-----------------|
| <input type="checkbox"/> Full Tray COVER - Modena Black Vinyl | \$190.00 |
|---|-----------------|

Lap Tray - Retaining Strap

- | | | | |
|---|---------------|--|---------------|
| <input type="checkbox"/> Full Tray - Standard Strap - Center Buckle | \$0.00 | <input type="checkbox"/> Full Tray - Standard Strap - Left-Hand Buckle | \$0.00 |
| <input type="checkbox"/> Full Tray - Standard Strap - Right-Hand Buckle | \$0.00 | <input type="checkbox"/> Full Tray - Sew-On Strap - Left-Hand Buckle | \$0.00 |
| <input type="checkbox"/> Full Tray - Sew-On Strap - Right-Hand Buckle | \$0.00 | | |

Veteran**Military Veteran Order**

- | | |
|--|---------------|
| <input type="checkbox"/> VA Order "Honor Our Veterans" | \$0.00 |
|--|---------------|

Branch of Service Crest

- | | | | |
|---|---------------|--|---------------|
| <input type="checkbox"/> Crest U.S. Air Force | \$0.00 | <input type="checkbox"/> Crest U.S. Army | \$0.00 |
| <input type="checkbox"/> Crest U.S. Coast Guard | \$0.00 | <input type="checkbox"/> Crest U.S. Marine Corps | \$0.00 |
| <input type="checkbox"/> Crest U.S. Navy | \$0.00 | | |

Shipping & Freight

Initial Shipping Costs

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Ship to Authorized Broda Dealer | \$165.00 | <input type="checkbox"/> Ship Direct to Customer | \$200.00 |
|--|-----------------|--|-----------------|

Additional Delivery Services

- | | | | |
|--|-----------------|---|-------------------------|
| <input type="checkbox"/> Lift Gate Services
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$200.00 | <input type="checkbox"/> Residential Delivery Service
Delivery is to a residential location. | \$200.00 |
| <input type="checkbox"/> Notification Before Delivery
Contact is requested to the cosignee by the carrier PRIOR to shipment delivery. | \$85.00 | <input type="checkbox"/> Inside Delivery Service
Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | Min
\$250.00 |
| <input type="checkbox"/> Limited Access Delivery Service
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$165.00 | | |