

## Customer Details

### Billing Address

Contact Name \_\_\_\_\_ Facility: \_\_\_\_\_

Contact Phone: \_\_\_\_\_ Contact Email: \_\_\_\_\_

Street Address: \_\_\_\_\_

### Shipping Address

Contact Name \_\_\_\_\_ Facility: \_\_\_\_\_

Contact Phone: \_\_\_\_\_ Contact Email: \_\_\_\_\_

Street Address: \_\_\_\_\_

### Order Information

Date: \_\_\_\_\_

PO Number: \_\_\_\_\_

Customer Number: \_\_\_\_\_

Tag/Comments:

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_



### Standard Features

- Up to 14° Infinitely Adjustable Seat Tilt (Posterior)
- Up to 800lb Maximum Occupant Weight
- Up to 8° Infinitely Adjustable Seat Tilt (Anterior)
- Contoured Backrest
- Up to 65° Infinitely Adjustable Back Support Recline
- Comfort Tension Seating®
- Stainless Steel, Powder Coated Frame
- Width Adjusting Removable Arm Supports
- 4 x 6" Total Lock Stainless Steel Front & Rear Caster Wheels
- Head Support
- 2 x 8" Locking Middle Wheels with One-Touch Emergency Brake

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

## Product Options

### Revive bariatric tilt/recline shower/commode

Revive Bariatric Shower Commode **\$7,764.00**

### Anterior Tilt

With Anterior Tilt **\$0.00**  Without Anterior Tilt **\$0.00**

### Tilt Lifting Force (Stainless Gas Spring)

Regardless of cylinder selection, weight capacity will always be 800 lbs. Guideline only. Contact Broda for additional weight accommodation options

<input type="checkbox"/> 135 Lbs 2xPAB19 Cylinders 300N (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 135 Lbs. 2xPAB7 Cylinders 300N (With Anterior Tilt) <b>\$0.00</b>
<input type="checkbox"/> 270 Lbs 2xPAB20 Cylinders 600N (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 270 Lbs 2xPAB8 Cylinders 600N (With Anterior Tilt) <b>\$0.00</b>
<input type="checkbox"/> 360 Lbs 2xPAB21 Cylinders 800N (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 360 Lbs 2xPAB9 Cylinders 800N (With Anterior Tilt) <b>\$0.00</b>
<input type="checkbox"/> 405 Lbs 2xPAB19 Cylinders 300N 2xPAB20 Cylinders 600N (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 405 Lbs 2xPAB7 Cylinders 300N 2xPAB8 Cylinders 600N (With Anterior Tilt) <b>\$0.00</b>
<input type="checkbox"/> 540 Lbs 4xPAB20 Cylinders 600N (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 540 Lbs 4xPAB8 Cylinders 600N (With Anterior Tilt) <b>\$0.00</b>
<input type="checkbox"/> 630 Lbs 2xPAB20 Cylinders 600N 2xPAB21 800N (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 630 Lbs 2xPAB8 Cylinders 600N 2xPAB9 800N (With Anterior Tilt) <b>\$0.00</b>
<input type="checkbox"/> 720 Lbs 4xPAB21 800N Cylinders (Without Anterior Tilt) <b>\$0.00</b>	<input type="checkbox"/> 720 Lbs 4xPAB9 Cylinders 800N (With Anterior Tilt) <b>\$0.00</b>

### Seat Type

Commode Seat With pan & Basket **\$0.00**  Transport Seat **\$0.00**

### Seat Depth

21" Seat Depth **\$0.00**  23" Seat Depth **\$0.00**

**Seat Height**

- |  |        |  |        |
|--|--------|--|--------|
| <input type="checkbox"/> 21" Seat Height | \$0.00 | <input type="checkbox"/> 23" Seat Height | \$0.00 |
|--|--------|--|--------|

**Strap Color**

- |                                     |        |                                      |        |
|-------------------------------------|--------|--------------------------------------|--------|
| <input type="checkbox"/> Blue Strap | \$0.00 | <input type="checkbox"/> Green Strap | \$0.00 |
|-------------------------------------|--------|--------------------------------------|--------|

- |                                     |        |  |  |
|-------------------------------------|--------|--|--|
| <input type="checkbox"/> Grey Strap | \$0.00 |  |  |
|-------------------------------------|--------|--|--|

**Arm Pad**

- |  |        |  |        |
|--|--------|--|--------|
| <input type="checkbox"/> Narrow Arm Supports | \$0.00 | <input type="checkbox"/> Soft Arm Supports<br>Extends width to 33" | \$0.00 |
|--|--------|--|--------|

**Bed Pan**

- |   |        |  |        |
|---|--------|--|--------|
| <input type="checkbox"/> Bed Pan - Short Handle<br>For clearance in washing apparatus | \$0.00 | <input type="checkbox"/> Bed Pan - Standard Handle | \$0.00 |
|---|--------|--|--------|

**Accessories**

- |   |         |  |         |
|---|---------|--|---------|
| <input type="checkbox"/> Molded Bolsters (Pair) | \$95.00 | <input type="checkbox"/> Pelvic Positioning Belt (Side Buckle) | \$64.00 |
|---|---------|--|---------|

**Veteran****Military Veteran Order**

- |  |        |  |  |
|--|--------|--|--|
| <input type="checkbox"/> VA Order "Honor Our Veterans" | \$0.00 |  |  |
|--|--------|--|--|

## Shipping &amp; Freight

## Initial Shipping Costs

- |  |                 |  |                 |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Ship to Authorized Broda Dealer | <b>\$165.00</b> | <input type="checkbox"/> Ship Direct to Customer | <b>\$200.00</b> |
|--|-----------------|--|-----------------|

## Additional Delivery Services

- |  |                 |   |                         |
|--|-----------------|---|-------------------------|
| <input type="checkbox"/> Lift Gate Services<br>Delivery is made with a lift gate trailer that moves from trailer level to ground level.              | <b>\$200.00</b> | <input type="checkbox"/> Residential Delivery Service<br>Delivery is to a residential location.   | <b>\$200.00</b>         |
| <input type="checkbox"/> Notification Before Delivery<br>Contact is requested to the cosignee by the carrier PRIOR to shipment delivery.             | <b>\$85.00</b>  | <input type="checkbox"/> Inside Delivery Service<br>Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | <b>Min<br/>\$250.00</b> |
| <input type="checkbox"/> Limited Access Delivery Service<br>Delivery is made with a lift gate trailer that moves from trailer level to ground level. | <b>\$165.00</b> |   |                         |