

Customer Details

Billing Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Shipping Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Order Information

Date: _____

PO Number: _____

Customer Number: _____

Tag/Comments: _____



Standard Features

- Up to 14° Infinitely Adjustable Seat Tilt (Posterior)
- Up to 800lb Maximum Occupant Weight
- Up to 8° Infinitely Adjustable Seat Tilt (Anterior)
- Contoured Backrest
- Up to 65° Infinitely Adjustable Back Support Recline
- Comfort Tension Seating®
- Stainless Steel, Powder Coated Frame
- Width Adjusting Removable Arm Supports
- 4 x 6" Total Lock Stainless Steel Front & Rear Caster Wheels
- Head Support
- 2 x 8" Locking Middle Wheels with One-Touch Emergency Brake

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

Product Options

Revive bariatric tilt/recline shower/commode

Revive Bariatric Shower Commode **\$8,198.00**

Anterior Tilt

With Anterior Tilt **\$0.00** Without Anterior Tilt **\$0.00**

Tilt Lifting Force (Stainless Gas Spring)

Regardless of cylinder selection, weight capacity will always be 800 lbs. Guideline only. Contact Broda for additional weight accommodation options

- | | |
|--|---|
| <input type="checkbox"/> 135 Lbs 2xPAB19 Cylinders 300N (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 135 Lbs. 2xPAB7 Cylinders 300N (With Anterior Tilt) \$0.00 |
| <input type="checkbox"/> 270 Lbs 2xPAB20 Cylinders 600N (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 270 Lbs 2xPAB8 Cylinders 600N (With Anterior Tilt) \$0.00 |
| <input type="checkbox"/> 360 Lbs 2xPAB21 Cylinders 800N (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 360 Lbs 2xPAB9 Cylinders 800N (With Anterior Tilt) \$0.00 |
| <input type="checkbox"/> 405 Lbs 2xPAB19 Cylinders 300N 2xPAB20 Cylinders 600N (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 405 Lbs 2xPAB7 Cylinders 300N 2xPAB8 Cylinders 600N (With Anterior Tilt) \$0.00 |
| <input type="checkbox"/> 540 Lbs 4xPAB20 Cylinders 600N (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 540 Lbs 4xPAB8 Cylinders 600N (With Anterior Tilt) \$0.00 |
| <input type="checkbox"/> 630 Lbs 2xPAB20 Cylinders 600N 2xPAB21 800N (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 630 Lbs 2xPAB8 Cylinders 600N 2xPAB9 800N (With Anterior Tilt) \$0.00 |
| <input type="checkbox"/> 720 Lbs 4xPAB21 800N Cylinders (Without Anterior Tilt) \$0.00 | <input type="checkbox"/> 720 Lbs 4xPAB9 Cylinders 800N (With Anterior Tilt) \$0.00 |

Seat Type

Commode Seat With pan & Basket **\$0.00** Transport Seat **\$0.00**

Seat Depth

21" Seat Depth **\$0.00** 23" Seat Depth **\$0.00**

Seat Height

- | | | | |
|--|--------|--|--------|
| <input type="checkbox"/> 21" Seat Height | \$0.00 | <input type="checkbox"/> 23" Seat Height | \$0.00 |
|--|--------|--|--------|

Strap Color

- | | | | |
|-------------------------------------|--------|--------------------------------------|--------|
| <input type="checkbox"/> Blue Strap | \$0.00 | <input type="checkbox"/> Green Strap | \$0.00 |
|-------------------------------------|--------|--------------------------------------|--------|

- | | | | |
|-------------------------------------|--------|--|--|
| <input type="checkbox"/> Grey Strap | \$0.00 | | |
|-------------------------------------|--------|--|--|

Arm Pad

- | | | | |
|--|--------|--|--------|
| <input type="checkbox"/> Narrow Arm Supports | \$0.00 | <input type="checkbox"/> Soft Arm Supports
Extends width to 33" | \$0.00 |
|--|--------|--|--------|

Bed Pan

- | | | | |
|---|--------|--|--------|
| <input type="checkbox"/> Bed Pan - Short Handle
For clearance in washing apparatus | \$0.00 | <input type="checkbox"/> Bed Pan - Standard Handle | \$0.00 |
|---|--------|--|--------|

Accessories

- | | | | |
|---|----------|--|---------|
| <input type="checkbox"/> Molded Bolsters (Pair) | \$102.00 | <input type="checkbox"/> Pelvic Positioning Belt (Side Buckle) | \$66.00 |
|---|----------|--|---------|

Veteran**Military Veteran Order**

- | | | | |
|--|--------|--|--|
| <input type="checkbox"/> VA Order "Honor Our Veterans" | \$0.00 | | |
|--|--------|--|--|

Shipping & Freight

Initial Shipping Costs

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Ship to Authorized Broda Dealer | \$165.00 | <input type="checkbox"/> Ship Direct to Customer | \$200.00 |
|--|-----------------|--|-----------------|

Additional Delivery Services

- | | | | |
|--|-----------------|---|-------------------------|
| <input type="checkbox"/> Lift Gate Services
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$200.00 | <input type="checkbox"/> Residential Delivery Service
Delivery is to a residential location. | \$200.00 |
| <input type="checkbox"/> Notification Before Delivery
Contact is requested to the cosignee by the carrier PRIOR to shipment delivery. | \$85.00 | <input type="checkbox"/> Inside Delivery Service
Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | Min
\$250.00 |
| <input type="checkbox"/> Limited Access Delivery Service
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$165.00 | | |