

Customer Details

Billing Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Shipping Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Order Information

Date: _____

PO Number: _____

Customer Number: _____

Tag/Comments: _____



Standard Features

- 2 x 24" Rear Mag Wheels
- Swing Away Footrests w/ Heel Loops (14"â€"18" knee-to-heel range)
- 2 x 5" Total Lock Front Stainless Caster Wheels
- Adjustable Seat Surface Height
- 10° Back Angle (adjustable to 5° and 15°)
- Adjustable Seat Depth
- Swing Up Arm Supports
- Adjustable Rear Mag Position (14"â€"21")
- Stainless Steel Adjustable Anti-tippers
- Adjustable Seat Angle (0° to -6°)
- 300lb Maximum Occupant Weight
- Tension Adjustable Back Straps
- Stainless Steel, Powder Coated Frame

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

Product Options

Chair Model

- | | | | |
|--|------------|--|------------|
| <input type="checkbox"/> Sentinel Shower Commode | \$1,854.00 | <input type="checkbox"/> Sentinel Shower Commode with Mag Wheels | \$2,266.00 |
|--|------------|--|------------|

Seat Type

- | | | | |
|---|--------|--------------------------------------|----------|
| <input type="checkbox"/> Molded Seat With Lid | \$0.00 | <input type="checkbox"/> Padded Seat | \$206.00 |
|---|--------|--------------------------------------|----------|

Back Type

- | | | | |
|---|--------|--|--------|
| <input type="checkbox"/> Back Sling With Slots (Grey) | \$0.00 | <input type="checkbox"/> Back Sling Without Slots (Grey) | \$0.00 |
|---|--------|--|--------|

Seat Width

- | | | | |
|---|----------|---|----------|
| <input type="checkbox"/> 18" Seat Width | \$0.00 | <input type="checkbox"/> 20" Seat Width | \$412.00 |
| <input type="checkbox"/> 22" Seat Width | \$515.00 | | |

Seat Depth

- | | | | |
|---|----------|---|----------|
| <input type="checkbox"/> 16" Seat Depth | \$103.00 | <input type="checkbox"/> 17" Seat Depth | \$0.00 |
| <input type="checkbox"/> 18" Seat Depth | \$103.00 | <input type="checkbox"/> 19" Seat Depth | \$103.00 |

Seat Height

- | | | | |
|--|--------|--|--------|
| <input type="checkbox"/> 19.5" Seat Height | \$0.00 | <input type="checkbox"/> 23.5" Seat Height | \$0.00 |
|--|--------|--|--------|

Back Height

- | | | | |
|--|--------|--|--------|
| <input type="checkbox"/> 16" Back Height | \$0.00 | <input type="checkbox"/> 18" Back Height | \$0.00 |
| <input type="checkbox"/> 20" Back Height | \$0.00 | | |

Wheel Options

- | | | | |
|---|--------|--|--------|
| <input type="checkbox"/> 1 Rear Directional Lock Caster | \$0.00 | <input type="checkbox"/> 2 Rear Total Lock Casters | \$0.00 |
| <input type="checkbox"/> Rear 24" Mag Wheels | \$0.00 | | |

Wheel Lock Extensions

- | | |
|---|---------|
| <input type="checkbox"/> Parking Brake Extensions | \$81.00 |
|---|---------|

Arm Pad

- | | | | |
|--|--------|--|--------|
| <input type="checkbox"/> Narrow Arm Supports | \$0.00 | <input type="checkbox"/> Soft Arm Supports | \$0.00 |
|--|--------|--|--------|

Swing Away Footrest

- | | | | |
|---|--------|--|--------|
| <input type="checkbox"/> With Footrests | \$0.00 | <input type="checkbox"/> Without Footrests | \$0.00 |
|---|--------|--|--------|

Bed Pan

- | | | | |
|--|---------|---|--------|
| <input type="checkbox"/> With Bed Pan Lid and Basket | \$64.00 | <input type="checkbox"/> Without Bed Pan Lid and Basket | \$0.00 |
|--|---------|---|--------|

Accessories

- | | | | |
|---|----------|--|---------|
| <input type="checkbox"/> Molded Bolsters (Pair)
Pre-cut slots in back included | \$102.00 | <input type="checkbox"/> Pelvic Positioning Belt (Center Buckle) | \$66.00 |
| <input type="checkbox"/> Pelvic Positioning Belt (Side Buckle) | \$66.00 | | |

Shipping & Freight

Initial Shipping Costs

- | | | | |
|--|-----------------|--|-----------------|
| <input type="checkbox"/> Ship to Authorized Broda Dealer | \$165.00 | <input type="checkbox"/> Ship Direct to Customer | \$200.00 |
|--|-----------------|--|-----------------|

Additional Delivery Services

- | | | | |
|--|-----------------|---|-------------------------|
| <input type="checkbox"/> Lift Gate Services
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$200.00 | <input type="checkbox"/> Residential Delivery Service
Delivery is to a residential location. | \$200.00 |
| <input type="checkbox"/> Notification Before Delivery
Contact is requested to the cosignee by the carrier PRIOR to shipment delivery. | \$85.00 | <input type="checkbox"/> Inside Delivery Service
Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | Min
\$250.00 |
| <input type="checkbox"/> Limited Access Delivery Service
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$165.00 | | |