

Customer Details

Billing Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Shipping Address

Contact Name _____ Facility: _____

Contact Phone: _____ Contact Email: _____

Street Address: _____

Order Information

Date: _____

PO Number: _____

Customer Number: _____

Tag/Comments:



Standard Features

- Crash Tested and WC19 Certified for Vehicle Transport
- Pivot Assistâ„¢ System
- Seat Width: 20â€ (Arm Support Width Adjustable to 24â€)
- Seat Surface Height: 22.5â€â€“32â€ (Foot Pump Adjustable)
- Seat Depth: 18.5â€
- Back Height: 32â€
- Maximum Occupant Weight: 450 lbs
- Chest Harness
- Up to 90° Infinitely Adjustable Back Support Recline
- Up to 20° Infinitely Adjustable Seat Tilt

Use our [Anatomical Assessment Form](#) or [Bariatric Anatomical Assessment Form](#) to properly fit your wheelchair.

Options

Chair Model

- Traversa Wheelchair **\$9,400.00**

Accessories

Neckrest Pads

- Standard Neck Support (Black) **\$0.00** Contoured Visco Neck Support (Vinyl) (Black) **\$95.00**
- Contoured Visco Neck Support (Triptex) **\$95.00**

WC19 Belt Extender

- WC19 Pelvic Positioning Belt Extender **\$62.00**

Side Reduction Pad

- Side Reduction Pads (Black) **\$91.00**

Accessories - Add-On

- Total Transport Package Discount - Includes IV Pole O2 Holder Transfer Board & Sleeve **-\$168.00**
Save 10%

Transfer Board

- Transfer Board With Holder **\$1,113.00** Transfer Board Without Holder **\$1,082.00**
- Transfer Board Holder Only **\$31.00**

Additional Positioning Padding (APP)-Slim Posture Padding (SPP)

Additional Positioning Padding (APP) - Package Price

Broda APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Covers (Not compatible with 35" back support height)

APP Package - Positioning - 16" 18" 20" Widths **\$641.00**

Slim Posture Padding (SPP) - Seat Cushion

SPP Seat Cushion 20" Seat Width X 18.5" Seat Depth (E2601) **\$212.00**

Armrest**Arm Top**

Standard Adjustable Arm (Top) **\$0.00**

Long Adjustable Arm (Top) **\$0.00**
Replaces standard arm supports

Accessory Bag**Accessory Bag**

Accessory Bag (Single) **\$51.00**

Accessory Bag (Pair) **\$102.00**

Embroidered Accessory Bag (Single) **\$83.00**

Embroidered Accessory Bag (Pair) **\$165.00**

Accessory Bag - Embroidery

Accessory Bag Embroidery - Block Font **\$0.00**

Accessory Bag Embroidery - Brush Font **\$0.00**

Accessory Bag - Characters

Accessory Bag - 2 Characters **\$0.00**

Accessory Bag - 1 Character Only **\$0.00**

Accessory Bag - 1st Character

Accessory Bag - 2nd Character

IV Pole**Add IV-Pole - MT**

IV-Pole - MT **\$284.00**

IV Pole - Factory Installed on Chair

IV Pole - For Traversa **\$0.00**

IV Pole - Frame Color

IV Pole - Black MT Side Bracket **\$0.00**

IV Pole - Mount

IV Pole - Traversa IV-Pole and Cam Handle **\$164.00**

IV Pole - Location

IV Pole - Side Mount **\$0.00**

O2 Holder**Add Oxygen Holder - MT**

Oxygen Holder - MT **\$280.00**

O2 Holder - Factory Installed on Chair

Oxygen Holder - For Traversa **\$0.00**

O2 Holder - Frame Color

Oxygen Holder - Black MT With Velcro + Neo Wrap **\$77.00**

O2 Holder - Mounting Type

Oxygen Holder - Standard Oxygen Holder Clamps Only **\$0.00**
If selecting rear mag wheels, the IV-POLE and O2-HOLDER must be front-mounted

O2 Holder - Mounting Location

Oxygen Holder - Rear Mount **\$0.00**

Shipping & Freight**Initial Shipping Costs**

Freight Charges **\$350.00** Freight Charges (Ontario Province Only) **\$150.00**

Additional Delivery Services

- | | | | |
|--|-----------------|---|---------------------|
| <input type="checkbox"/> Lift Gate Services
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$200.00 | <input type="checkbox"/> Residential Delivery Service
Delivery is to a residential location. | \$200.00 |
| <input type="checkbox"/> Notification Before Delivery
Contact is requested to the cosignee by the carrier PRIOR to shipment delivery. | \$85.00 | <input type="checkbox"/> Inside Delivery Service
Movement of the shipment beyond the unloading position, i.e. move closer to door, inside building, upstairs, down a hall, etc. Price is specified per order, based on the customer's inside delivery needs. | Min \$250.00 |
| <input type="checkbox"/> Limited Access Delivery Service
Delivery is made with a lift gate trailer that moves from trailer level to ground level. | \$165.00 | | |